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Sterisart[®] Universal-Pumpe 16419

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Contents

1	General Safety Information	4
2	Intended Use of the Sterisart® Universal Pump	5
3	Important Notes on the Safe and Proper Usage of the Sterisart® Universal Pump	5
4	Unpacking the Pump and Checking the Equipment Supplied.	6
5	Description of the Sterisart® Universal Pump.	7
6	Startup.	9
7	Application.	10
8	Operating the Sterisart® Universal Pump 16419	11
9	Cleaning and Sanitization	12
10	Installation in an Isolator.	13
11	Important Notes Concerning the Warranty	13
12	Troubleshooting.	14
13	Technical Specifications.	15
14	Ordering Information for the Sterisart® Universal Pump	16
15	Ordering Information for Consumables	17
16	Literature.	18
17	Storage and Shipping.	18
	17.1 Storage	18
	17.2 Returning Devices and Parts	18
18	Disposal	19
	18.1 Information on Decontamination	19
	18.2 Disposing of Device and Parts	19
	18.2.1 Information on Disposal	19
	18.2.2 Disposal	19
19	Conformity.	20
20	Legal information about licenses and trademarks	20
	20.1 Millipore™	20

1 General Safety Information

Explanation of Safety Symbols Used

These operating instructions contain important information which is visibly marked. This information must be observed in order to avoid hazards to users as well as to avoid damage and | or malfunctions to the device.

Symbols in the operating instructions:

Important Note! Important information
Danger! General danger with risk of injury

Laser warning

General information and | or help

Symbols on the manufacturer's ID label

AC voltage

Return and disposal

Observe instructions for use

Manufacturer

2 Intended Use of the Sterisart® Universal Pump

The new Sterisart® Universal pump is a peristaltic pump equipped with appropriate devices for attaching a closed sterility test system and a sample container.

The pump draws liquid from the sample container into the specially designed double tubing, distributing the liquid equally to both sample containers of the sterility test system. The sample is rinsed in a similar way. After the nutrient media is added, the sample containers are incubated for a prescribed time and evaluated.

3 Important Notes on the Safe and Proper Usage of the Sterisart® Universal Pump

- The Sterisart® Universal pump must be switched off before any electrical connections are made or broken (e.g. when plugging in or unplugging the power cable or foot switch).
 - The Sterisart® Universal pump must only be used with the original current-carrying parts supplied with it (power cable with ferrite core, D-Sub connecting cable and foot switch). The D-Sub protective cap supplied must cover the port when the D-Sub connecting cable is not connected.
 - The Sterisart® Universal pump is designed for use with Sterisart® NF filtration systems. The suitability of any other tubing must be approved by Sartorius.
 - Only plastic and glass bottles with a filling volume of 500 ml or less are permitted to be placed in the delivered bottle holder.
 - Only authorized Sartorius service technicians are allowed to open the pump housing.
 - Do not expose the Sterisart® Universal pump to extreme temperatures.
 - Improper use or handling will invalidate the manufacturer's guarantee.
 - Sterisart® Universal pump may only be used with a grounded, readily accessible and fuse-protected electrical outlet that complies with national standards.
-

4 Unpacking the Pump and Checking the Equipment Supplied

When unpacking the Sterisart® Universal pump 16419 make sure not to leave any parts in the packaging. The equipment supplied includes the following accessories (one of each):

- Sterisart® Universal pump 16419
- Power cable (depends on country rating)
- Container for draining
- Pump cover for container for draining (Sterisart® units)
- Drain tubing
- Fastening arm
- Bottle holder
- Support rod
- Bow
- Clamping lever with sleeve
- Operating instructions German | English

5 Description of the Sterisart[®] Universal Pump


Fig. 1

The Sterisart[®] Universal pump 16419 (see Fig. 1) consists of pump housing (1) with pump head (2), lever (3), bow (4) and handles (5). On the right handle, there is also a bolt with a tilting lock for the container for draining (11). The selecting knob (6) is at the front on the right. The stand (7) with the holding ring for bottles (8) and fastening arm (9) and the optional ampoule breaker (10) is located on the left side of the pump.

The container for draining (Fig. 1; 11) for the sterility test containers is located on the right side of the pump. The drain tubing can be attached to the appropriate connector (barb) on the container for draining (Fig. 1; 12).

Laser warning

The internal barcode scanner emits a laser beam.
Do not look directly into the beam.

Fig. 2

Located on the rear panel (see Figure 2): a laser warning label (1), On | Off switch (2), AC power socket (3), socket for the foot switch (4), an RS-232 interface port for connecting a PC/MMC card reader (5) and the manufacturer's ID label (6).

6 Startup

Set up the Sterisart® universal pump ergonomically in a suitable area (cleanroom, laminar flow work bench, isolator).

Before initial start-up, the accessories supplied must be assembled on the Sterisart® Universal pump. Please proceed in the following order:

- At the back left, push the support rod (Fig.1, 7) onto the pin of the upper part of the pump until it is firmly in place. The holes in the foot of the support rod allow the positioning of the support rod to be altered.
 - Fasten the bottle holder (Fig. 1, 8) and the fastening arm (Fig. 1, 9) to the support rod. If necessary, attach the ampoule breaker (Fig. 1, 10) beforehand.
 - Insert the bow (Fig. 1, 4) and the clamping lever (Fig. 1, 3) on the corresponding bolts on the top of the pump. The lever is jacketed in a plastic sleeve, which can be removed for cleaning. If signs of wear appear on this sleeve after it has been in use for some time, it can be re-ordered as a spare part (see chapter "Ordering information for the Sterisart® Universal Pump").
 - After connecting the drain tubing to the barb (Fig. 1, 12) of the container for draining (Fig. 1, 11); attach the container for draining to the appropriate bolts on the right handle (Fig. 1, 5). The container for draining and the cover plate are attached by moving the tilting lock.
- Fasten the end of the drain tubing to the appropriate drain.

Very Important Note!

Always make sure that there is a sufficient liquid level differential (> 5 cm) between the container for draining and the drain. This makes safe drainage possible and prevents backflow into the container for draining. Make sure that the drain tubing is kept lower than the pump to prevent the drain pan from overflowing.

Treat the waste as a biohazard. Only dispose of it in accordance with current regulations.

- Plug the foot switch into the designated socket located on the rear panel (see Fig. 2, 4). Fasten the plug to the socket by turning it clockwise.
- Plug the power cable into the male power socket on the back of the Sterisart® Universal pump and into a grounded electrical outlet (Fig. 2, 3).
- Use the On|Off switch (Fig. 2, 2) at the rear side of the Sterisart® Universal pump to turn it on.

7 Application

Caution!
General danger with risk of injury

Do not operate the Sterisart® Universal pump until the tubing has been properly installed. The pump may only be operated with the tubing inserted and only by one person at a time. Never pump any flammable, corrosive or explosive fluids. Users should always wear tight-fitting clothing (no ties, etc.).

To be able to operate the pump, Sterisart® sterility test systems must be installed beforehand. To do this, follow the steps given below:

Remove the sterility test systems from their respective packaging and place the sterility test sample containers in the container for draining on the right side of the pump.

- Now open the lever and the bow.
- Keep the tubing around the rotor of the pump (the arrow on the fixed part surrounding the pump head indicates the subsequent direction of rotation), so that no kinks occur and the tubing is not under tension (Fig. 3). The tubing is attached by simply inserting it (with it lying on the ground) and fastening the lever (Fig. 4). In this way, the tubing is automatically correctly positioned.

Important Note!
When threading the tubing, make sure that it is not kinked or twisted.

Fig. 3

- To start the pumping process, the tubing must be fully inserted and the lever must be engaged. A safety function ensures that the pump does not start if the lever or the bow are open.

- After threading the tubing and closing the lever, insert the removal spike/needle of the Sterisart® NF system as described in the instructions for the Sterisart® NF sterility test system. Then attach the test liquid container (e.g., bag or bottle) by hanging it on the support rod (bag) or by placing it in the holding ring on the support rod (bottle). The holding ring accepts bottles of various sizes. The holder must be turned so that the sterility test system tubing is not under tension after the spike | needle has been inserted and so that the spike | needle cannot slip out of the container.

For testing ampoules, an ampoule breaker is also attached to the support rod. Use the integrated knife to scratch the neck of the ampoule, then snap off the top. The broken-off top of the ampoule falls into the collection container (which can be emptied later).

- For the following pumping process, make sure that the Sterisart® NF sample containers are correctly positioned in the intended holders, and that the tubing has been inserted without tension or kinking.

Important Note!

When the lever is opened, the Sterisart® Universal pump is turned off immediately by a safety function. In this event, the Sterisart® Universal pump can be restarted after adequately closing all components.

8 Operating the Sterisart® Universal Pump 16419

The Sterisart® Universal Pump 16419 can be started up either by pressing the selecting knob or by activating the optional foot switch. The speed can be smoothly adjusted using the selecting knob. The pump is switched off by again pressing the selecting knob or by releasing the foot switch.

Fig. 4

9 Cleaning and Sanitization

Caution!
Ports on the back panel are not tight.

The Sterisart® Universal pump is not suited for sanitization (gaseous or atomized) with peracetic acid, hydrogen peroxide, or other decontamination procedures commonly used in isolators, unless it is built into an isolator countertop. These types of disinfectant can be used when the pump has been professionally installed in the isolator.

- Before cleaning the pump rotor, turn off the main power switch and unplug the power cable.

The following parts of the pump have to be removed for cleaning:

- Bow, lever, support rod and container for draining
- Disconnect the pump head:
 - The rotor must be turned by hand until the marking on the upper side (black dot) is pointing backwards.
 - The marking indicates the position of the fastening screw that is located in the rotor mid-section.
 - Loosen this screw with a 4-mm Allen wrench. Do not turn the screw out completely.
 - Pull the rotor upwards.
- Pull the speed-selecting knob upwards.

Comply with the IP protection rating! No liquid may enter the Sterisart® Universal pump. Disconnect the pump from line voltage before any cleaning, maintenance, or repair work.

Cleaning stainless steel surfaces:

The stainless steel housing can be cleaned at regular intervals. Clean stainless steel parts with a wet cloth. Customary cleaning agents suitable for stainless steel can be used without risk. For surface sanitization, use commercially available disinfectants that are suitable for stainless steel.

Do not use any cleaning agents on the stainless steel housing that contain sodium hydroxide or acetic, sulfuric or citric acid. The use of steel wool cleaning sponges is prohibited.

All stainless steel surfaces can be cleaned by wiping with a moist cloth and with the following products:

- Peracetic acid, H₂O₂ and acetic acid (2.5%)
- Quaternary ammonium
- Formaldehyde and glutaraldehyde
- Isopropyl alcohol
- Ethanol 70%

Autoclavable at 121°C for 20 min. are bow, lever, speed selecting knob, support rod, bottle holder, counter-bottle holder and container for draining (complete).

Never autoclave the pump rotor and never submerge it into any liquid. Only use a slightly moistened cloth | tissue for cleaning | disinfecting.

After cleaning, all parts should be secured as described in Chapter "6 Startup," page 9. The rotor is assembled as follows:

- The rotor shaft has a groove on the side turned towards the back.
- Push the rotor upwards onto the shaft.
- Turn the marking on the rotor towards the back, so that the groove on the rotor shaft aligns with the marking on the rotor
- Tighten the screw in the rotor to guarantee that the rotor is firmly seated.

10 Installation in an Isolator

The installation must be carried out by the isolator manufacturer.

To install the Sterisart® Universal pump in an isolator, it must be opened up first. Make sure to unplug it from AC current (mains supply) before opening it!

During assembly, please follow the installation instructions precisely (part of the installation kit for isolators, order number 1ZE---0050). The instructions are available separately from the manufacturer.

11 Important Notes Concerning the Warranty

The correct functioning and safe use of the Sterisart® Universal pump is only guaranteed when it is operated according to the information given in these instructions.

The equipment manufacturer is responsible only for effects on safety, dependability and performance of the instrument if assembly, extensions, new settings, changes or repairs are carried out by persons authorized by him to do so, if the electrical installation in the room is in conformance with IEC requirements and if the equipment is used in accordance with these operating instructions.

Any changes made to the Sterisart® Universal pump will invalidate the manufacturer's warranty unless they are previously agreed to by authorized Sartorius Service staff, or carried out by them.

12 Troubleshooting

This section lists a few problems that could arise in connection with the Sterisart[®] pump. Several solutions are given for each problem. The first solution given is usually the most likely one. If this does not solve the problem, go down the list of suggestions until the problem has been solved.

Important Note!

The pump must be switched off before electrical connections are unplugged or plugged in!

Should it not be possible to locate and solve the problem in this way, contact the appropriate service center (see the back cover of this manual).

Problem	Solution
No function; no display	<p>The pump is not switched on, or is not correctly connected to the electric supply.</p> <p>Plug the power cable into the electric outlet correctly and to the Sterisart[®] Universal pump. Switch on the on-off switch.</p> <p>Check the electric supply; switch on multiple outlet sockets; check the wiring.</p> <p>If these possibilities are not feasible, the primary fuse in the main line filter may be defective.</p> <p>Changing this, with the power cable unplugged, must only be performed by authorized personnel!</p>
The pump rotor does not rotate.	<p>The locking lever is open or the foot switch is not properly connected.</p> <p>The bow has been removed from the pump.</p> <p>Check the connection of the foot switch to the pump.</p>

13 Technical Specifications

Dimensions pump in mm	(W × D × H) approx. 336 × 260 × 210 (with lever)
Dimensions pump with holding ring for bottles, container for draining in mm	(W × D × H) 440 × 365 × 485
Weight: Basic model 16419	Approx. 13.5 kg
Maximum operating temperature	40°C
Further data	IP41 in accordance to IEC 60529/EN 60529

Electrical connection

Voltage	100 – 240 VAC ± 10%
Frequency	50 – 60 Hz
Power consumption	100 VA

Pump

Pump flow rate	70 – 650 ml/min (depending on the tubing)
----------------	---

Tubing

Correct functioning of the Sterisart® Universal pump is only guaranteed when used with tubing that has been approved by Sartorius.

Materials

Base plate, lever, bow, rotor	Stainless steel 316L (1.4435)
Housing, connection port, support rod, bottle holder	Stainless steel 316 (1.4401)
Connecting cable	PVC
Container for draining, plastic sleeve for lever, hydrodynamic bearing for bow	PET
Hydrodynamic bearing for lever	Iglidur X
Seals, gaskets, drain tubing	Silicone

14 Ordering Information for the Sterisart® Universal Pump

Order No.	Description
16419	Sterisart® Universal pump, basic model
16973	Ampoule tongs
1ZG---0009	Bottle holder
1ZF---0004	Fastening arm
1ZA---0002	Drain tubing
1ZW---0002	Ampoule breaker
1ZE---0033	Foot switch
1ZF---0013	Stand
1ZG---0028	Container for draining
1ZB---0001	Plastic sleeve for lever
1ZF---0007	Wing nut
1ZS---0005	Bow
1ZS---0002	Lever
1ZL---0003	Rotor
1ZG---0023	Cover for container for draining (Sterisart® sterility test units)
1ZG---0024	Cover for container for draining (Millipore™ sterility test units)
1ZG---0014	Adaptor for Sterisart® container for draining for use with a Millipore™ Equinox pump
1ZG---0039	Transport trolley
1ZE---0041	Installation kit isolators

Additional accessories are available on request.

15 Ordering Information for Consumables

Recommended disposable sterility test units for use with pump.

Order No.	Description
16466-ACD	Sterisart [®] NF alpha, dual-needle metal spike for closed containers (box of 10, individually sterilized with ETO; single-packed).
16467-ACD	Sterisart [®] NF alpha, 6-cm metal needle for open containers (box of 10, individually sterilized with ETO; single-packed)
16468-ACD	Sterisart [®] NF alpha, system for medical devices with Luer lock connectors (box of 10, individually sterilized with ETO; single-packed)
16466-GBD	Sterisart [®] NF gamma, dual-needle metal spike for closed containers (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16467-GBD	Sterisart [®] NF gamma, 6-cm metal needle for open containers (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16468-GBD	Sterisart [®] NF gamma, system for medical devices with Luer or Luer lock connectors (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16469-GBD	Sterisart [®] NF gamma, system with adapter for prefilled syringes (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16470-GBD	Sterisart [®] NF gamma, system for slightly soluble powders in unvented vials (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16475-GBD	Sterisart [®] NF gamma, system for lyophilized or freely soluble powders in unvented vials (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16476-GBD	Sterisart [®] NF gamma, system with short dual-needle metal spike for closed containers (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16477-GBD	Sterisart [®] NF gamma, system for liquids in plastic containers (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16478-GBD	Sterisart [®] NF gamma, system for medical containers with Luer lock male connectors (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.

Additional Sterisart[®] units are available on request.

16 Literature

Data sheet Sterisart® Universal Pump Publ. No.: SLD2010-e	English Order No.: 85030-531-17
--	------------------------------------

Brochure Sterisart® Universal Pump Publ. No.: SLD1003-e	English Order No.: 85030-531-14
--	------------------------------------

17 Storage and Shipping

17.1 Storage

Procedure

- ▶ If the device is in operation:
 - ▶ Decommission the device.
- ▶ Clean the device.
- ▶ Store the device according to the ambient conditions (see Chapter "13 Technical Specifications," page 15).

17.2 Returning Devices and Parts

Defective devices or parts can be sent back to Sartorius. Returned devices must be clean, decontaminated, and properly packed.

Transport damage as well as measures for subsequent cleaning and disinfection of the device or parts by Sartorius shall be charged to sender.

WARNING

Risk of injury due to contaminated devices!

Devices contaminated with hazardous materials [nuclear, biological or chemical (NBC)] will **not** be accepted for repair or disposal.

- ▶ Observe the information on decontamination (see Chapter "13 Technical Specifications," page 15).

Procedure

- ▶ Decommission the device.
- ▶ Contact Sartorius Service for instructions on how to return equipment or parts (please refer to our website at www.sartorius.com for return instructions).
- ▶ Pack the device and the parts in an un-used original packaging.
The original packaging can be ordered (contact Sartorius Service).

18 Disposal

18.1 Information on Decontamination

The device does **not** contain any hazardous materials that would necessitate special disposal actions. The media (e.g. acids, bases) used during the process are potentially hazardous materials that could cause biological or chemical hazards.

According to the EU Directive on Hazardous Substances, the owners of devices that come into contact with hazardous substances are responsible for properly disposing of these devices and for declaring such devices when transporting them.

WARNING

Risk of injury due to contaminated devices!

Devices contaminated with hazardous materials [nuclear, biological or chemical (NBC)] will **not** be accepted for repair or disposal.

18.2 Disposing of Device and Parts

18.2.1 Information on Disposal

The device and the device accessories must be disposed of properly by disposal facilities.

A lithium cell battery is installed inside the device. Batteries must be disposed of properly by disposal facilities.

The packaging is made of environmentally friendly materials that can be used as secondary raw materials.

18.2.2 Disposal

Requirements

The device has been decontaminated.

Procedure

- ▶ Dispose of the device. Follow the disposal instructions on our website (www.sartorius.com). Inform the disposal facility that there is a lithium cell battery installed inside the device.
- ▶ Dispose of the packaging in accordance with local government regulations.
- ▶ Dispose of the consumables in accordance with local government regulations.

19 Conformity

The attached documents confirm compliance of the device with the directives or standards cited.

20 Legal information about licenses and trademarks

20.1 Millipore™

Millipore™ is a registered trademark of Merck KGaA, Darmstadt, Germany or its affiliates.

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in the form as delivered fulfils all the relevant provisions of the following European Directives – including any amendments valid at the time this declaration was signed – and meets the applicable requirements of the harmonized European Standards including any amendments valid at the time this declaration was signed – listed below:

2014/30/EU Elektromagnetische Verträglichkeit | *Electromagnetic compatibility*
EN 61326-1:2013

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The person authorised to compile the technical file:

Sartorius Lab Instruments GmbH & Co. KG
Electronics & Product Compliance
37070 Goettingen, Germany

Sartorius Lab Instruments GmbH & Co. KG
Goettingen, 2019-04-18

Dr. Reinhard Baumfalk
Head of Product Development
Lab Products and Services Division

Dr. Dieter Klausgrete
Senior Scientist Certification Management
Product Compliance Officer

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Project: 80003373

Date Issued: 2019-07-24

Issued To: Sartorius Lab Instruments GmbH & Co. KG
Otto-Brenner-Strasse 20
Goettingen, Niedersachsen, 37079
Germany

The products listed below are eligible to bear the CSA Mark shown with adjacent indicators 'C' and 'US' for Canada and US or with adjacent indicator 'US' for US only or without either indicator for Canada only.

Issued by: *Eric Fortier*
Eric Fortier

PRODUCTS

CLASS - C872106 - ELECTRICAL LABORATORY EQUIPMENT

CLASS - C872186 - ELECTRICAL LABORATORY EQUIPMENT Certified to US Standards

Peristaltic pumps, table top, portable, AC powered, Class I (grounded), models and rating as follow;
16419 and 16420, each rated: 100-240 V~, 50-60 Hz, 100 VA

Notes:

1. The above model is cord connected, Equipment Class I, Pollution Degree 2 and Overvoltage category II.
2. Mode of operation: Continuous
3. Environmental Conditions: +10 to +40 °C, 2000 m. max, 80% rH max.

Certificate: 80003373
Project: 80003373

Master Contract: 167555
Date Issued: 2019-07-24

CONDITIONS OF ACCEPTABILITY

- (1) If at any time there is a conflict between the system safety provisions and any relevant local (national or regional) requirements, the local requirements always take precedence.
- (2) For use in dry locations and indoor only.
- (3) No user replaceable battery. Internal non rechargeable lithium cell can only be replaced by Service personnel.
- (4) Appliance coupler on power cord is considered as a disconnect device, it must remain accessible.
- (5) Professional laboratory equipment; only trained personal is allowed to operate the unit.
- (6) Pump is only allowed to be operated once the tubing has been properly installed into the unit.

APPLICABLE REQUIREMENTS

- | | | |
|---|---|--|
| CAN/CSA-C22.2 No. 61010-1-12
+ UPD1:15 + UPD2:16 + AMD1:2018 | - | Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use, Part 1: General Requirements |
| ANSI/UL 61010-1-2012 +AMD1:2018 | - | Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use - Part 1: General Requirements |

Supplement to Certificate of Compliance

Certificate: 80003373

Master Contract: 167555

*The products listed, including the latest revision described below,
are eligible to be marked in accordance with the referenced Certificate.*

Product Certification History

Project	Date	Description
80003373	2019-07-24	CSA c/us certification of a sterilization pump, model Sterisart based on CB report from SIQ

Sartorius Lab Instruments GmbH & Co. KG
Otto-Brenner-Strasse 20
37079 Goettingen, Germany

Phone: +49.551.308.0
www.sartorius.com

The information and figures contained in these instructions correspond to the version date specified below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Masculine or feminine forms are used to facilitate legibility in these instructions and always simultaneously denote the other gender as well.

Copyright notice:

This instruction manual, including all of its components, is protected by copyright.

Any use beyond the limits of the copyright law is not permitted without our approval.

This applies in particular to reprinting, translation and editing irrespective of the type of media used.

© Sartorius Germany

Last updated:
06 | 2019