## Instructions For Use (IFU)

# Ultroser™ G

## Serum Substitute for Animal Cell Culture

### 1. Introduction

Ultroser G is a serum substitute used to replace FBS (Fetal Bovine Serum) in animal cell culture media. Ultroser G serum substitute is success fully used in the *in vitro* culture of many anchorage-dependent cells. It is also used in prenatal diagnostics for the culture of amniotic fluid cells for karyotype study. The biological activity is 5 times that of FBS. The product is supplied as a lyophilized powder corresponding to 20 mL of substitute.

### 2. Preparation and Solubilization

The lyophilized powder is solubilized in 20 mL sterile distilled pyrogen-free water (water for injection) as follows:

1. Open the vial by breaking the tip of the capsule,
2. Introduce the appropriate volume of water using:
  - A sterile pyrogen-free pipette after first having removed the stopper,
  - Or a sterile pyrogen-free syringe by injection through the stopper.

These two operations must be performed under sterile conditions.

3. Allow the lyophilisate to swell for 10 to 15 minutes. Solubilization is rapidly obtained by gently shaking the flask or by successive aspirations using a pipette or syringe.

The reconstituted solution should be clear and pink-coloured. For greater safety, the solution may be filtered through a sterilizing membrane of 0.2 µm pore diameter. This will not result in any loss of activity.

### 3. Concentration

After reconstitution, Ultroser G serum substitute has a biological activity 5 times greater than that of FBS. It should be used at concentrations between 0.5% and 4% (2.5% to 20% of FBS). A concentration of 2% (10% of FBS) is generally used. To obtain a standard concentration of 2%, mix 2% Ultroser G serum substitute with 98% base culture medium. For each cell type, various concentrations from 0.5% to 4% will have to be evaluated to determine the optimal growth conditions. The effects on cell growth are not always proportional to the concentration of Ultroser G serum substitute, nor in correlation with the results obtained with FBS in the same conditions of use.

### 4. Culture Medium

The culture medium may be the one generally used for culture with FBS (MEM, etc.). However, the use of an enriched medium such as HAM's F12, DMEM, IMDM, or RPMI is recommended.

### 5. Culture Support

The culture support is the same as that used with FBS. Ultroser G serum substitute contains adhesion factors. Therefore, it can be used for pre-incubation of the plastic dishes in the medium prior to adding the cells. However, for particular cell types which are difficult to grow, it may be necessary, to improve adhesion and differentiation, to coat the plastic dishes with one of the following products: poly-D-lysine, collagen, an extracellular matrix, or adhesion factors such as fibronectin.

### 6. Trypsinization

The trypsinization techniques for tissues and mono-layered cells are the same as those used with FBS. However, the trypsinization period should be as short as possible, to minimize cell damage. To remove any residual trypsin, it may be necessary to add a trypsin inhibitor such as soya bean extract or aprotinin, at a concentration of 0.5 to 1 mg/mL, before dispersing the cells. If the cells still remain aggregated to each other, the cell suspension may be passed through a 40 µm cell sieve.

### 7. Cell Transfer

Two techniques may be employed for sub-culture from a medium containing FBS:

- Direct transfer of the cells in the medium supplemented with Ultroser G serum substitute,
- Or after an adaptation period by progressively reducing the quantity of serum, in several passages:
  - 1<sup>st</sup> passage: 50% medium supplemented with serum (MSS) + 50% medium supplemented with Ultroser G serum substitute (MSU) in the optimum concentration
  - 2<sup>nd</sup> passage: 25% MSS + 75% MSU
  - 3<sup>rd</sup> passage: 10% MSS + 90% MSU
  - 4<sup>th</sup> passage: 100% MSU

This second technique allows better adaptation of the cells to the new culture medium. Before seeding the cells, check the pH of the final medium and, if necessary, adjust to 7,2 ± 0.1. Sodium bicarbonate or a synthetic buffer (HEPES) may be added to improve the buffering power of the medium.

### 8. Storage

Ultroser G serum substitute is available as a lyophilized powder. The product is shipped at ambient temperature and must then be stored between 2 and 8 °C. Once reconstituted, the solution must be used up quickly. Other wise it can be aliquoted and stored frozen at -20 °C for several weeks. Frequent freezing and thawing must be avoided.

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### 11. Symbols Used on the Label

	Manufacturer		Temperature limitation
	Catalog number		<i>In vitro</i> diagnostic medical device
	Batch code		Biological risks
	Use by		Consult instructions for use

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Specifications subject to change without notice. Copyright Sartorius Stedim Biotech GmbH. DIR: 2730265-000-00 Status: 03 | 2021
Sartorius Stedim Chromatography Resins SAS