



Lab Ultrafiltration and Purification Products

Simplifying Progress


SARTORIUS

Table of Contents


General Information	
Major Uses for Ultrafiltration	4
Sartorius Lab Ultrafiltration Devices	5
Process Alternatives	5
Membrane Performance Characteristics	6
Membrane Selection Guide	7
Protein and Macromolecule Concentration	9
Centrifugal Filtration	
▪ Vivaspin® 500	10
▪ Vivaspin® 2	12
▪ Centrisart® 1	15
▪ Vivaspin® Turbo 4	17
▪ Vivaspin® 6	19
▪ Vivaspin® 15R	22
▪ Vivaspin® Turbo 15 PES	24
▪ NEW! Vivaspin® Turbo 15 RC	26
▪ Vivaspin® 20	28
▪ Vivaclear	32
Pressure-Fugation	
▪ Vivacell 100	33
Tangential Flow Filtration	
▪ Vivaflow® 50	36
▪ Vivaflow® 50R	39
▪ Vivaflow® 200	41
Solvent Absorption	
▪ Vivapore® 5, 10 20	44
Ultrafiltration Membrane Discs	46

DNA Concentration	49
Vivacon® 500	50
Vivacon® 2	53
Protein Purification	57
Vivapure® Ion Exchange Protein Purification Products	58
Virus Purification and Concentration	61
Vivapure® Virus Purification and Concentration Kits	62
Adenovirus Purification with Vivapure® Adenopack Kits	63
▪ Vivapure® Adenopack 20	64
▪ Vivapure® Adenopack 100	65
▪ Vivapure® Adenopack 500	67
Lentivirus Purification with Vivapure® Lentiselect Kits	68
▪ Vivapure® Lentiselect 40	69
▪ Vivapure® Lentiselect 500	70
▪ Vivapure® Lentiselect 1000	71
Application Notes	73
1. Desalting and Buffer Exchange with Vivaspin® Centrifugal Concentrators	74
2. Treatment of Vivaspin® Concentrators for Improved Recovery of Low-concentrated Protein Samples	77
3. Scouting Protein Purification Conditions Using Vivapure® Centrifugal Ion Exchange Membrane Absorbers	80
4. Sartorius Ultrafiltration Products for the Concentration and Purification of Viruses – a Short Review	85
5. Sartorius Ultrafiltration Products in the Preparation of Biological Nanoparticles and Medical Nanocarriers	97
6. Vivaflow® and Vivaspin® Workflow in Protein Research Laboratories	103


Applications for Ultrafiltration



Solute concentration



Solute fractionation or clarification



Solute desalting or purification

Ultrafiltration is a convective process using anisotropic semi-permeable membranes to separate macromolecular species and solvents – primarily on the basis of size. It is particularly appropriate for the concentration of macromolecules and can also be used for purification or solvent exchange. Ultrafiltration is a gentle, non-denaturing method that is more efficient and flexible than alternative processes.

Typical applications for ultrafiltration

- Concentration | desalting of proteins, enzymes, DNA, monoclonal antibodies, immunoglobulins, extracellular vesicles, viruses and nanoparticles
- Bence Jones Protein concentration from urine samples prior to capillary electrophoresis
- Forensic DNA sample concentration prior to sequencing reaction
- Peptide fractionation in FASP (filter-aided sample preparation)
- Free drug | hormone assays
- Removal of primers from PCR amplified DNA
- Removal of labelled amino acids and nucleotides
- HPLC sample preparation
- Deproteinization of samples
- Recovery of biomolecules from cell culture supernatants | lysates
- Mammalian cell harvesting
- Cell washing, virus purification, cell debris removal and depyrogenation
- Environmental sample clarification | concentration

Solute concentration

Ultrafiltration membranes are used to increase the solute concentration of a desired biological or inorganic species. Macromolecules are retained by the membrane when they are significantly larger than the nominal pore size, while microsolute are removed convectively with the solvent.

Solute fractionation or clarification

Ultrafiltration is a cost effective method for separating samples into size-graded components providing that the desired fractions have at least a 10-fold difference in molecular weight. During filtration, the permeating solute remains at its initial concentration whilst the retained macromolecules will be enriched.


Solute desalting or purification

A solution may be purified from salts, non-aqueous solvents and low molecular weight materials. Multiple solvent exchanges, will progressively purify macromolecules from contaminating solutes. Microsolute are removed most efficiently by adding solvent to the solution being filtered at a rate equal to the speed of filtration. This is called diafiltration.


Solute concentration for diagnostics

Ultrafiltration devices can be used in the clinical setting for the concentration and separation of disease markers, such as Bence Jones Protein for multiple myeloma diagnostic sample prep, from clinical samples, such as blood serum, urine and cerebrospinal fluid. Use of devices for these applications require dedicated *in vitro* Diagnostic (IVD) registered devices. IVD devices are only available in registered countries, according to country specific regulations. Please contact Sartorius for more information on registered countries and availability.


Sartorius Lab Ultrafiltration Devices




Centrifugal concentrators



Pressure-fugation



Cross flow



Solvent absorption

Sartorius develops devices dedicated to optimising laboratory ultrafiltration processes with minimal time requirements while maximizing recovery, reliability and robustness.

In addition Sartorius are continually building technical and application

support tools to help users select the optimum device and process for their sample type.

Visit www.sartorius.com for more technical and application support material.

Alternative Ultrafiltration Methods

Sartorius offers a comprehensive choice of operating methods for sample ultrafiltration and concentration. Below is a guide to selecting the most suitable ultrafiltration method, depending on sample volume, equipment available, and the desired filtration speed and process control.

Centrifugal (0.1 to 100 mL starting volumes)

Centrifugation provides the vector to clear solvent and micromolecules through the ultrafiltration membrane and into a filtrate container positioned below. This is a gentle process that is characterized by quick set up and fast filtration speeds with most solutions. Twelve centrifugal devices are offered from the Vivaspin®, Centrisart®, Vivacell and Vivacon® families.

Pressure (5 to 100 mL starting volumes)

Pressurized air or an inert gas is used to provide the vector for filtration. For fastest filtration, Vivacell products are used with an orbital shaker, where agitation impedes macromolecules from polarizing on the membrane surface. Vivaspin® 20 and Vivacell 100 can be operated using gas pressure.

Pressure-Fugation (5 to 15 mL starting volumes)

A unique Sartorius method that combines gas pressure with centrifugation, with process times typically 30 to 50% faster than centrifugation alone. Vivaspin® 20 can be operated this way.

Crossflow | TFF (0.1 to 5 L starting volumes)

The solution to be processed is pumped under pressure across an ultrafiltration membrane and then returned to the original reservoir. The solution is progressively concentrated or purified as solvent and micromolecules pass through the membrane into a separate filtrate vessel. Vivaflow® cassettes are offered for this method.

Solvent Absorption (1 to 20 mL starting volumes)

An absorbent cellulose pad mounted behind the ultrafiltration membrane draws solvents and microsolute through the membrane. Retained macromolecules are concentrated into the bottom of the sample container. No additional equipment is required. Vivapore® is offered for this technique.

Membrane Performance Characteristics

Sartorius offers an extended range of membranes to cover the great majority of ultrafiltration requirements.

The following is a guide to selecting the most appropriate membranes according to their typical performance characteristics. Please note however, that membrane behavior and ultimate performance, largely depends on the specific characteristics of the solution being processed. Sartorius recommends that users experiment with alternative membrane materials to optimize their process performance.

Polyethersulfone (PES)

This is a low binding membrane that provides excellent performance with most solutions and exceptional recovery of negatively charged target molecules. Polyethersulfone membranes are usually preferred for their low fouling characteristics, exceptional flux and broad pH compatibility.

Regenerated Cellulose (RC)

The Sartorius regenerated cellulose membrane has been uniquely developed to ensure optimal performance

in the lab ultrafiltration devices.

This is a hydrophilic membrane suitable for general samples, with ultra-low protein adsorption and high chemical compatibility.

Regenerated cellulose is especially well suited to ultrafiltration of oligonucleotides and peptides.

Cellulose triacetate (CTA)

High hydrophilicity and very low non specific binding characterize this membrane. Cast without any support that could trap or bind passing microsolute, these membranes are preferred for sample cleaning and protein removal, and when high recoveries from the filtrate solution is of primary importance.

Hydrosart®

Demonstrating the same properties as regenerated cellulose, but with the added benefit of enhanced performance characteristics and extremely low protein binding. Hydrosart® is another membrane of choice for applications such as concentration and desalting of immunoglobulin fractions.

Membrane performance comparisons

Membrane	Frequently preferred for:
Polyethersulfone & Regenerated Cellulose 3 kDa MWCO 5 kDa MWCO 10 kDa MWCO 30 kDa MWCO 50 kDa MWCO 100 kDa MWCO	Concentration Desalting Buffer exchange Fractionation
Cellulose triacetate 5 kDa MWCO 10 kDa MWCO 20 kDa MWCO	Deproteinization Free bound drug studies Whenever the filtrate is being analyzed
Hydrosart® 2 kDa MWCO 5 kDa MWCO 10 kDa MWCO 30 kDa MWCO	Concentration Desalting Buffer exchange Fractionation Membrane evaluation for scale up

Membrane Selection Guide

The advanced designs and low adsorption materials that characterize Sartorius products, offer a unique combination of faster processing speeds and higher recovery of the concentrated sample. Providing that the appropriate device size and membrane cut-off is selected, Sartorius products will typically yield recoveries of the concentrated sample in excess of 90% when the starting sample contains over 0.1 mg/mL of the solute of interest. The majority of any loss is caused by non-specific binding both to the membrane surface and to exposed binding sites on the plastic of the sample container:

Adsorption to the membrane

Depending on sample characteristics relative to the membrane type used, solute adsorption on the membrane surface is typically 2-10 µg/cm². This can increase to 20-100 µg/cm² when the filtrate is of interest and the solute must pass through the whole internal structure of the membrane. Typically a higher cut-off membrane will bind more than a low molecular weight alternative.

Adsorption to the sample container

Although every effort is made to minimize this phenomenon by the selection of low binding materials and tool production to optical standards, some

solute will bind to the internal surface of the sample container. Whilst the relative adsorption will be proportionately less important than on the membrane, due to the higher total surface area, this can be the major source of yield loss.

Process optimization

When highest recoveries are crucial, in particular when working with solute quantities in the microgram range, Sartorius recommends that users consider the following:

- Select the smallest device that suits the sample volume.
- Take advantage of the extra speed of Sartorius products by refilling a smaller device repeatedly.
- Select the lowest MWCO membrane that suits the application.
- Reduce pressure or centrifugal force to approximately half of the recommended maximum.
- Avoid over concentration. The smaller the final concentrate volume, the more difficult it is to achieve complete recovery.
- If feasible, after sample retrieval, rinse the device with one or more drops of buffer.
- Pretreat the device overnight with a passivation solution such as 5% SDS, Tween 20, or Triton X-100, then rinse thoroughly before use.

Membrane selection guide (recommended MWCO)

Application	< 5 kDa	10 kDa	30 kDa	50 kDa	100 kDa	> 300 kDa
Bacteria					■	■
DNA fragments		■	■	■	■	
Enzymes	■	■				
Extracellular Vesicle					■	■
Growth factors	■	■				
Immunoglobulins			■	■	■	
mAb			■	■	■	
Nucleic Acids	■	■	■	■	■	
Oligonucleotides	■					
Peptides	■					
Viruses			■	■	■	
Yeast		■	■	■	■	■

For highest recovery, select a membrane MWCO which is a maximum one third to half of the molecular weight of the solute to be retained



SARTORIUS

VIVAFLOW 50

SARTORIUS

VIVAFLOW 50

SARTORIUS

Protein and Macromolecule Concentration

Table of Contents

Vivaspin® 500	10
Vivaspin® 2	12
Centrisart® 1	15
Vivaspin® Turbo 4	17
Vivaspin® 6	19
Vivaspin® 15R	22
Vivaspin® Turbo 15 PES	24
Vivaspin® Turbo 15 RC	26
Vivaspin® 20	28
Vivaclear	32
Vivacell 100	33
Vivaflow® 50	36
Vivaflow® 50R	39
Vivaflow® 200	41
Vivapore® 5, 10 20	44
Ultrafiltration Membrane Discs	46

Vivaspin® 500



100 to 500 µL samples

Vivaspin® 500 centrifugal filter units offer a simple, one step procedure for sample preparation. They can effectively be used in fixed angle rotors accepting 2.2 mL centrifuge tubes.

The legacy patented vertical membrane design and thin channel filtration chamber (US 5,647,990), minimizes membrane fouling and provides fast concentrations – even with particle-loaded solutions.

Technical specifications

Concentrator capacity

Swing bucket rotor	do not use
Fixed angle rotor	500 µL

Dimensions

Total length	50 mm
Diameter	11 mm
Active membrane area	0.5 cm ²
Hold-up volume, membrane and support	< 5 µL
Dead-stop volume	5 µL

Materials of construction

Body	Polycarbonate
Filtrate vessel	Polypropylene
Concentrator cap	Polycarbonate
Membrane	Polyethersulfone


Equipment required

Centrifuge

Rotor type	Fixed angle
Minimum rotor angle	40°
Rotor cavity	To fit 2.2 mL (11 mm) conical bottom tubes
Maximum speed	12,000 g

Concentrate recovery

Pipette type	Fixed or variable volume
Recommended tip	Thin gel loader type




Performance characteristics

	Time to concentrate up to 30x at 20°C and solute recovery	
Rotor	Fixed angle	
Centrifugal force	12,000 <i>g</i>	
Start volume	500 μ L	
	Time	Recovery
Aprotinin 0.25 mg/mL (6.5 kDa) 3 kDa MWCO PES	30 min	96%
BSA 1.0 mg/mL (66 kDa) 5 kDa MWCO PES	15 min	96%
10 kDa MWCO PES	5 min	96%
30 kDa MWCO PES	5 min	96%
IgG 0.25 mg/mL (160 kDa) 30 kDa MWCO PES	10 min	96%
50 kDa MWCO PES	10 min	96%
100 kDa MWCO PES	10 min	96%

Ordering information

Vivaspin® 500 PES	Pack size	Prod. no.
3 kDa MWCO	25	VS0191
3 kDa MWCO	100	VS0192
5 kDa MWCO	25	VS0111
5 kDa MWCO	100	VS0112
10 kDa MWCO	25	VS0101
10 kDa MWCO	100	VS0102
30 kDa MWCO	25	VS0121
30 kDa MWCO	100	VS0122
50 kDa MWCO	25	VS0131
50 kDa MWCO	100	VS0132
100 kDa MWCO	25	VS0141
100 kDa MWCO	100	VS0142
300 kDa MWCO	25	VS0151
300 kDa MWCO	100	VS0152
1,000 kDa MWCO	25	VS0161
1,000 kDa MWCO	100	VS0162
0.2 μ m	25	VS0171
0.2 μ m	100	VS0172

Vivaspin® 2



0.4 to 3 mL samples

The Vivaspin® 2 bridges the gap between the 500 µL and 4 mL centrifugal concentrators. This device combines the speed of the classic Vivaspin® products with low internal surface and membrane area for superior recoveries from very dilute solutions.

Also unique to the Vivaspin® 2 is the choice of directly pipetting the concentrate from the dead-stop pocket built into the bottom of the concentrator, or alternatively reverse spinning into the concentrator recovery cap. Both methods result in near total concentrate recoveries.


Available with a choice of polyether-sulfone, cellulose triacetate or Hydrosart® membranes, Vivaspin® 2 offers the highest flexibility for process optimization.

Technical specifications


Concentrator capacity	
Swing bucket rotor	3 mL
Fixed angle rotor	2 mL
Dimensions	
Total length	126 mm
Diameter	17 mm
Active membrane area	1.2 cm ²
Hold-up volume, membrane	< 10 µL
Dead-stop volume	8 µL
Materials of construction	
Body	Polycarbonate
Filtrate vessel	Polycarbonate
Concentrator cap	Polycarbonate
Membrane	PES, CTA, HY

Equipment required

Centrifuge		
Rotor type	Swing bucket	Fixed angle
Minimum rotor angle	-	25°
Rotor cavity	To fit 15 mL (17 mm) conical bottom tubes	To fit 15 mL (17 mm) conical bottom tubes
Maximum speed	4,000 g	8,000 g
Concentrate recovery		
Pipette type	Fixed or variable volume	Fixed or variable volume
Recommended tip	Thin gel loader type	Thin gel loader type



PES, CTA, or Hydrosart® membranes; Fits rotors accepting standard 15 mL tubes



Reverse spin concentrate retrieval

Performance characteristics

	Time to concentrate up to 30× at 20°C and solute recovery	
Rotor	Fixed angle	
Centrifugal force	5,000 <i>g</i>	
Start volume	2 mL	
	Time	Recovery
Aprotinin 0.25 mg/mL (6.5 kDa)		
3 kDa MWCO PES	50 min	96%
BSA 1.0 mg/mL (66 kDa)		
5 kDa MWCO PES	12 min	98%
5 kDa MWCO Hydrosart®	22 min	98%
10 kDa MWCO PES	8 min	98%
10 kDa MWCO CTA	10 min	96%
10 kDa MWCO Hydrosart®	12 min	98%
20 kDa MWCO CTA	5 min	96%
30 kDa MWCO PES	8 min	97%
30 kDa MWCO Hydrosart®	5 min	97%
IgG 0.25 mg/mL (160 kDa)		
20 kDa MWCO CTA	6 min	97%
30 kDa MWCO PES	10 min	96%
50 kDa MWCO PES	10 min	96%
100 kDa MWCO PES	8 min	95%


Ordering tips

- Choose a membrane pore size at least 50% smaller than the size of the molecule to be retained.
- Usually choose PES membranes for fastest concentrations.
- Usually choose CTA for protein removal | filtrate recovery.
- Usually choose Hydrosart® membranes for highest recovery of immunoglobulins.

Ordering information

Vivaspin® 2 PES	Pack size	Prod. no.
3 kDa MWCO	25	VS0291
3 kDa MWCO	100	VS0292
5 kDa MWCO	25	VS0211
5 kDa MWCO	100	VS0212
10 kDa MWCO	25	VS0201
10 kDa MWCO	100	VS0202
30 kDa MWCO	25	VS0221
30 kDa MWCO	100	VS0222
50 kDa MWCO	25	VS0231
50 kDa MWCO	100	VS0232
100 kDa MWCO	25	VS0241
100 kDa MWCO	100	VS0242
300 kDa MWCO	25	VS0251
300 kDa MWCO	100	VS0252
1,000 kDa MWCO	25	VS0261
1,000 kDa MWCO	100	VS0262
0.2 µm	25	VS0271
0.2 µm	100	VS0272
Vivaspin® 2 CTA	Pack size	Prod. no.
10 kDa MWCO	25	VS02V1
10 kDa MWCO	100	VS02V2
20 kDa MWCO	25	VS02X1
20 kDa MWCO	100	VS02X2
Vivaspin® 2 Hydrosart®	Pack size	Prod. no.
2 kDa MWCO	25	VS02H91
2 kDa MWCO	100	VS02H92
5 kDa MWCO	25	VS02H11
5 kDa MWCO	100	VS02H12
10 kDa MWCO	25	VS02H01
10 kDa MWCO	100	VS02H02
30 kDa MWCO	25	VS02H21
30 kDa MWCO	100	VS02H22

Centrisart® 1



0.5–2.5 mL samples

Centrisart® 1 is a ready-to-use unit for small-volume, centrifugal ultrafiltration to separate proteins from low molecular weight substances in biological samples.

Centrisart® 1 features a unique design that enables ultrafiltration in the direction opposite to centrifugal force. This is so effective in preventing premature blockage of the filter that even whole blood samples can be deproteinized.

The ultrafiltrate is collected in the floating filtrate tube, where it is readily accessible without disassembly.

Centrisart® 1 is ideal for the following applications:

- Drug binding studies
- Isolation of metabolites from serum
- Protein removal from blood samples
- Cleaning of liposomes
- Virus removal

Technical specifications

Concentrator capacity

Swing bucket rotor	2.5 mL
Fixed angle rotor	2.5 mL

Dimensions

Total length	93 mm
Diameter	14 mm
Active membrane area	0.79 cm ²
Hold-up volume, membrane	< 5 µL
Dead-stop volume	100 µL

Materials of construction

Centrifuge tube	Polystyrene
Filtrate tube	Cellulose propionate
Concentrator cap	Polyethylene
Membrane	CTA, PES

Equipment required


Centrifuge

Rotor type	Swing bucket	Fixed angle
Minimum rotor angle	-	25°
Rotor cavity	To fit 15 mL (17 mm) conical flat bottom tubes	To fit 15 mL (17 mm) conical flat bottom tubes
Maximum speed	2,500 g	2,000 g


Concentrate recovery

Pipette type	Fixed or variable volume	Fixed or variable volume
Recommended tip	Thin gel loader type	Thin gel loader type


Easy-to-use




Remove filtrate tube, pour in sample



Replace filtrate tube



Centrifuge



Recover the filtrate...

...or use forceps to remove the filtrate tube and access the concentrate

Performance characteristics

	Time to filter 50% of sample volume	Time to filter 90% of sample volume	Passage of sample species volume
Centrifugal force	2,000 g	2,000 g	2,000 g
Start volume	2.5 mL	2.5 mL	2.5 mL
BSA 1.0 mg/mL (66 kDa) 5 kDa MWCO	300 min	–	0%
10 kDa MWCO	35 min	80 min	2%
20 kDa MWCO	9 min	20 min	2%
IgG 0.25 mg/mL (160 kDa) 100 kDa MWCO	13 min	35 min	3%
Blue Dextran 0.1 mg/mL (2,000 kDa) 300 kDa MWCO	9 min	25 min	28%

Ordering information

Centrisart® 1 CTA	Pack size	Prod. no.
5 kDa MWCO CTA	12	13229-E
10 kDa MWCO CTA	12	13239-E
20 kDa MWCO CTA	12	13249-E
Centrisart® 1 PES		
100 kDa MWCO PES (IVD device)	12	13269-E*
300 kDa MWCO PES	12	13279-E

References


P. Nebinger and Koel (1993).
Determination of acyclovir by ultrafiltration and high-performance liquid chromatography.
J. Chromatography **619**, 342-344

F. da Fonseca-Wollheim, K.-G. Heinze, K. Lomsky and H. Schreiner (1988).
Serum ultrafiltration for the elimination of endogenous interfering substances in creatinine determination.
J. Clin. Chem. Clin. Biochem. **26**, 523-525

R. H. Christenson, S. D. Studenberg, S. Beck-Davis and F. A. Sedor (1987).
Digoxin-like immunoreactivity eliminated from serum by centrifugal ultrafiltration before fluorescence polarization immunoassay of digoxin.
Clinical Chemistry **33**, 606-608

* IVD devices are available only in countries where registered according to local IVD regulations

Vivaspin® Turbo 4 PES



2 to 4 mL samples

Vivaspin® Turbo 4 PES offers the fastest sample concentration with the highest recoveries. This device can handle up to 4 mL sample volumes in swing bucket and fixed angle rotors that accept 15 mL conical bottom centrifuge tubes.

The optimized design and sleek internal profile ensure maximum process speeds all the way down to the last few microliters, resulting in more than 100-fold concentration.

UV joining technology provides a smooth transition between membrane and housing, allowing collection of the entire concentrated sample from the unique angular dead-stop pocket.

Technical specifications

Concentrator capacity

Swing bucket rotor	4 mL
Fixed angle rotor	4 mL

Dimensions

Total length	122.5 mm
Diameter	17 mm
Active membrane area	3.2 cm ²
Hold-up volume, membrane	< 10 µL
Dead-stop volume, swing bucket	40 µL
Dead-stop volume, fixed angle	30 µL

Materials of construction

Body	Styrene butadiene copolymer
Filtrate vessel	Polypropylene
Concentrator cap	Polypropylene
Membrane	Polyethersulfone

Equipment required

Centrifuge

Rotor type	Swing bucket	Fixed angle
Minimum rotor angle	-	25°
Rotor cavity	To fit 15 mL (17 mm) conical bottom tubes	To fit 15 mL (17 mm) conical bottom tubes
Maximum speed	4,000 g	7,500 g

Concentrate recovery

Pipette type	Fixed or variable volume	Fixed or variable volume
Recommended tip	Thin gel loader type	Thin gel loader type

Visit us at

www.sartorius.com/

VivaspinTurbo4

for further information.

Here you can find instructions on how to use Vivaspin® Turbo 4 PES for:

- Desalting and buffer exchange
- Preparation of biological nanoparticles and medical nanocarriers
- Concentration and purification of viruses
- Urine protein concentration
- Separation of proteins and metabolites for disease detection

Performance Characteristics

	Time to concentrate up to 30× at 20°C and solute recovery			
	Swing bucket		Fixed angle (25°)	
Rotor	4,000 g		7,500 g	
Centrifugal force*	4 mL		4 mL	
Start volume	4 mL		4 mL	
	Time	Recovery	Time	Recovery
Cytochrome c (12.4 kDa)				
3 kDa MWCO PES	60 min	98%	80 min	96%
5 kDa MWCO PES	40 min	95%	50 min	94%
Lysozyme (14.3 kDa)				
3 kDa MWCO PES	65 min	95%	70 min	93%
5 kDa MWCO PES	50 min	94%	60 min	92%
α-Chymotrypsin (25 kDa)				
10 kDa MWCO PES	10 min	95%	8 min	95%
BSA (66 kDa)				
10 kDa MWCO PES	10 min	98%	7 min	97%
30 kDa MWCO PES	8 min	96%	6 min	97%
IgG (160 kDa)				
30 kDa MWCO PES	18 min	94%	13 min	92%
50 kDa MWCO PES	16 min	93%	12 min	90%
100 kDa MWCO PES	17 min	94%	13 min	92%


Ordering information

Vivaspin® Turbo 4 PES	Pack size	Prod. no.
3 kDa MWCO	25	VS04T91
3 kDa MWCO	100	VS04T92
5 kDa MWCO	25	VS04T11
5 kDa MWCO	100	VS04T12
10 kDa MWCO	25	VS04T01
10 kDa MWCO	100	VS04T02
10 kDa MWCO (IVD device)	25	VS04T01IVD**
10 kDa MWCO (IVD device)	100	VS04T02IVD**
30 kDa MWCO	25	VS04T21
30 kDa MWCO	100	VS04T22
50 kDa MWCO	25	VS04T31
50 kDa MWCO	100	VS04T32
100 kDa MWCO	25	VS04T41
100 kDa MWCO	100	VS04T42

* 3,000 g (swing bucket) or 5,000 g (fixed angle) centrifugal force for 100 kDa MWCO devices.

** IVD devices are available only in countries where registered according to local IVD regulations

Vivaspin® 6



2 to 6 mL samples

Vivaspin® 6 concentrators have been developed to offer increased volume flexibility and performance.

Vivaspin® 6 can process an impressive 6 mL in either swing bucket or fixed angle rotors accepting standard 15 mL conical bottom centrifuge tubes.

Featuring twin vertical membranes for unparalleled filtration speeds and more than 100-fold concentration, the retentate volume is easily estimated from the printed graduations on the side of the concentrator. The modified dead-stop pocket further simplifies direct pipette retrieval of the final concentrate.

Technical specifications

Concentrator capacity

Swing bucket rotor	6 mL
Fixed angle rotor	6 mL

Dimensions

Total length	122 mm
Diameter	17 mm
Active membrane area	2.5 cm ²
Hold-up volume, membrane	< 10 µL
Dead-stop volume	30 µL

Materials of construction

Body	Polycarbonate
Filtrate vessel	Polycarbonate
Concentrator cap	Polypropylene
Membrane	Polyethersulfone

Equipment required

Centrifuge

Rotor type	Swing bucket	Fixed angle
Minimum rotor angle	-	25°
Rotor cavity	To fit 15 mL (17 mm) conical bottom tubes	To fit 15 mL (17 mm) conical bottom tubes
Maximum speed	4,000 g	8,000 g
Maximum speed, ≥100 kDa MWCO	4,000 g	6,000 g

Concentrate recovery

Pipette type	Fixed or variable volume	Fixed or variable volume
Recommended tip	Thin gel loader type	Thin gel loader type

Performance characteristics

	Time to concentrate up to 30× at 20°C and solute recovery			
	Swing bucket		Fixed angle (25°)	
Rotor	Swing bucket		Fixed angle (25°)	
Centrifugal force	3,000 g		7,500 g	
Start volume	6 mL		6 mL	
	Time	Recovery	Time	Recovery
Cytochrome c 0.25 mg/mL (12.4 kDa) 5 kDa MWCO PES	-	-	90 min	97%
BSA 1.0 mg/mL (66 kDa) 5 kDa MWCO PES	20 min	98%	12 min	98%
10 kDa MWCO PES	13 min	98%	10 min	98%
30 kDa MWCO PES	12 min	98%	9 min	97%
IgG 0.25 mg/mL (160 kDa) 30 kDa MWCO PES	18 min	96%	15 min	95%
50 kDa MWCO PES	17 min	96%	14 min	95%
100 kDa MWCO PES	15 min	91%	12 min	91%
Latex beads 0.004% in DMEM + 10% FCS (55 nm) 300 kDa MWCO PES	-	-	25 min	99%
Latex beads 0.004% in DMEM + 10% FCS (240 nm) 1,000 kDa MWCO PES	-	-	4 min	99%
Yeast 1.0 mg/mL (<i>S. Cerevisiae</i>) 0.2 µm PES	4 min	97%	3 min	97%

Ordering information

Vivaspin® 6 PES	Pack size	Prod. no.
3 kDa MWCO	25	VS0691
3 kDa MWCO	100	VS0692
5 kDa MWCO	25	VS0611
5 kDa MWCO	100	VS0612
10 kDa MWCO	25	VS0601
10 kDa MWCO	100	VS0602
10 kDa MWCO (IVD device)	25	VS0601IVD**
10 kDa MWCO (IVD device)	100	VS0602IVD**
30 kDa MWCO	25	VS0621
30 kDa MWCO	100	VS0622
50 kDa MWCO	25	VS0631
50 kDa MWCO	100	VS0632
100 kDa MWCO	25	VS0641
100 kDa MWCO	100	VS0642
300 kDa MWCO	25	VS0651
300 kDa MWCO	100	VS0652
1,000 kDa MWCO	25	VS0661
1,000 kDa MWCO	100	VS0662
0.2 µm	25	VS0671
0.2 µm	100	VS0672

** IVD devices are available only in countries where registered according to local IVD regulations

Vivaspin® 15R



2 to 15 mL samples

Vivaspin® 15R is designed for initial sample volumes up to 15 mL and features a modified regenerated cellulose membrane; Hydrosart®. This membrane is ideal where extremely high recovery with very low adsorption is needed. An example of this application includes desalting and concentration of immunoglobulin fractions.

- Ultimate recoveries (95 – 98%)
- Extremely short concentration time (30-fold in 15 minutes)
- Simple and convenient handling
- Easy scale-up to 0.1 to 5 L with Vivaflow® 50R or 200 with Hydrosart® membranes
- Very low hold-up volume (< 20 µL)

Technical specifications

Concentrator capacity

Swing bucket rotor	15 mL
Fixed angle rotor	12.5 mL

Dimensions

Total length	116 mm
Diameter	30 mm
Active membrane area	3.9 cm ²
Hold-up volume, membrane	< 20 µL
Dead-stop volume	30 µL

Materials of construction

Body	Polycarbonate
Filtrate vessel	Polypropylene
Concentrator cap	Polycarbonate
Membrane	Hydrosart®


Equipment required


Centrifuge

Rotor type	Swing bucket	Fixed angle
Minimum rotor angle	-	25°
Rotor cavity	To fit 50 mL (30 mm) conical bottom tubes	To fit 50 mL (30 mm) conical bottom tubes
Maximum speed	3,000 g	6,000 g


Concentrate recovery

Pipette type	Fixed or variable volume	Fixed or variable volume
Recommended tip	Thin gel loader type	Thin gel loader type





Spin



Recover

Performance characteristics


	Time to concentrate up to 30× at 20°C and solute recovery			
	Swing bucket		Fixed angle (25°)	
Rotor	3,000 g		6,000 g	
Centrifugal force	15 mL		12.5 mL	
Start volume	Time		Recovery	
Aprotinin 0.1 mg/mL* (6.5 kDa) 5 kDa MWCO	47 min	95%	45 min	95%
Cytochrome c 0.25 mg/mL* (12.4 kDa) 5 kDa MWCO	45 min	96%	45 min	96%
	25 min	94%	18 min	94%
α-chymotrypsin 0.25 mg/mL* (25 kDa) 5 kDa MWCO	50 min	98%	45 min	98%
	25 min	98%	18 min	98%
Ovalbumin 1.0 mg/mL* (45 kDa) 10 kDa MWCO	20 min	98%	14 min	98%
	15 min	94%	12 min	94%
BSA 1.0 mg/mL* (66 kDa) 30 kDa MWCO	18 min	98%	15 min	98%
	IgG 0.1 mg/mL in DMEM (160 kDa) 30 kDa MWCO	30 min	98%	25 min

Ordering information

Vivaspin® 15R Hydrosart®	Pack size	Prod. no.
2 kDa MWCO	12	VS15RH91
2 kDa MWCO	48	VS15RH92
5 kDa MWCO	12	VS15RH11
5 kDa MWCO	48	VS15RH12
10 kDa MWCO	12	VS15RH01
10 kDa MWCO	48	VS15RH02
30 kDa MWCO	12	VS15RH21
30 kDa MWCO	48	VS15RH22

* Proteins other than IgG made up in 50 mM potassium phosphate, 150 mM sodium chloride, pH 7.4

Vivaspin® Turbo 15 PES



4 to 15 mL samples

Vivaspin® Turbo 15 PES enables the fastest sample concentration and highest recoveries. This device can handle samples up to 15 mL in rotors accepting 50 mL centrifuge tubes. The optimized design and sleek internal profile ensure maximum process speeds all the way down to the last few microlitres, resulting in more than 100-fold concentration.

UV joining technology provides a smooth transition between membrane and housing, allowing collection of the entire concentrated sample from the unique, angular dead-stop pocket.

Stable polyethersulfone membranes are suited to a wide pH range and especially recommended for high recovery of negatively charged target molecules. Now complemented with a regenerated cellulose option, Vivaspin® Turbo offers the best membrane, whatever the sample.

Technical specifications

Concentrator capacity

Swing bucket rotor	15 mL
Fixed angle rotor (25°)	11 mL

Dimensions

Total length, concentrator insert assembled	77 118 mm
Diameter, concentrator insert	27 mm
Active membrane area	7.2 cm ²
Hold-up volume, membrane	<10 µL
Dead-stop volume, swing bucket fixed angle	100 60 µL

Materials of construction

Body	Styrene butadiene copolymer
Filtrate vessel	Polypropylene
Concentrator cap	Polypropylene
Membrane	Polyethersulfone

Equipment required

Centrifuge

	Swing bucket	Fixed angle
Rotor type	Swing bucket	Fixed angle
Minimum rotor angle	-	25°
Rotor cavity	To fit 15 mL (17 mm) conical bottom tubes	To fit 15 mL (17 mm) conical bottom tubes
Maximum speed	4,000 g	4,000 g
Maximum speed, ≥100 kDa MWCO	2,000 g	2,000 g

Concentrate recovery

Pipette type	Fixed or variable volume	Fixed or variable volume
Recommended tip	Thin gel loader type	Thin gel loader type

Performance characteristics


	Time to concentrate up to 20× at 20 °C and solute recovery			
	Swing bucket		Fixed angle (25°)	
Rotor	Swing bucket		Fixed angle (25°)	
Centrifugal speed	4,000 g		4,000 g	
Start volume	15 mL		11 mL	
	Time	Recovery	Time	Recovery
Cytochrome c 0.25 mg/mL (12.4 kDa) 5 kDa MWCO PES	30 min	98%	50 min	98%
Lysozyme 0.25 mg/mL (14.3 kDa) 5 kDa MWCO PES	33 min	96%	50 min	96%
α-Chymotrypsin 1.0 mg/mL (25 kDa) 10 kDa MWCO PES	10 min	95%	10 min	95%
BSA 1.0 mg/mL (66 kDa) 10 kDa MWCO PES	10 min	99%	10 min	99%
30 kDa MWCO PES	8 min	98%	10 min	98%
IgG 1.0 mg/mL (160 kDa) 30 kDa MWCO PES	23 min	95%	17 min	95%

Ordering information

Vivaspin® Turbo 15 PES	Pack size	Prod. no.
3 kDa MWCO	12	VS15T91
3 kDa MWCO	48	VS15T92
5 kDa MWCO	12	VS15T11
5 kDa MWCO	48	VS15T12
10 kDa MWCO	12	VS15T01
10 kDa MWCO	48	VS15T02
10 kDa MWCO (IVD device)	12	VS15T01IVD*
10 kDa MWCO (IVD device)	48	VS15T02IVD*
30 kDa MWCO	12	VS15T21
30 kDa MWCO	48	VS15T22
50 kDa MWCO	12	VS15T31
50 kDa MWCO	48	VS15T32
100 kDa MWCO	12	VS15T41
100 kDa MWCO	48	VS15T42

* IVD devices are available only in countries where registered according to local IVD regulations

NEW! Vivaspin® Turbo 15 RC



4 to 15 mL samples

Vivaspin® Turbo 15 RC enables the fastest sample concentration and highest recoveries. This device can handle samples up to 15 mL in rotors accepting 50 mL centrifuge tubes. The optimized design and sleek internal profile ensure maximum process speeds all the way down to the last few microlitres, resulting in more than 100-fold concentration.

Solvent-free, heat weld technology provides a smooth transition between

membrane and housing, allowing complete concentrate recovery from the unique, angular dead-stop pocket.

Regenerated cellulose membranes developed specifically for Sartorius lab ultrafiltration devices are suited to general samples, with ultra-low adsorption and high chemical compatibility, and especially recommended for oligonucleotides and peptides. Complemented with a polyethersulfone option, Vivaspin® Turbo offers the best membrane, whatever the sample.

Technical specifications

Concentrator capacity

Swing bucket rotor	15 mL
Fixed angle rotor (25°)	11 mL

Dimensions

Total length, concentrator insert assembled	77 118 mm
Diameter, concentrator insert	27 mm
Active membrane area	8.1 cm ²
Hold-up volume, membrane	<10 µL
Dead-stop volume, swing bucket fixed angle	100 60 µL

Materials of construction

Body	Styrene butadiene copolymer
Filtrate vessel	Polypropylene
Concentrator cap	Polypropylene
Membrane	Regenerated Cellulose

Equipment required

Centrifuge

Rotor type	Swing bucket	Fixed angle
Minimum rotor angle	-	25°
Rotor cavity	To fit 15 mL (17 mm) conical bottom tubes	To fit 15 mL (17 mm) conical bottom tubes
Maximum speed	4,000 g	6,000 g
Maximum speed, ≥100 kDa MWCO	3,000 g	-

Concentrate recovery

Pipette type	Fixed or variable volume	Fixed or variable volume
Recommended tip	Thin gel loader type	Thin gel loader type

Performance characteristics

	Time to concentrate up to 30× at 20 °C and solute recovery			
	Swing bucket		Fixed angle (25°)	
Rotor	Swing bucket		Fixed angle (25°)	
Centrifugal speed	4,000 g*		6,000 g	
Start volume	15 mL		11 mL	
	Time	Recovery	Time	Recovery
Cytochrome c 0.25 mg/mL (12.4 kDa) 5 kDa MWCO RC	23 min	94%	37 min	92%
Lysozyme 0.25 mg/mL (14.3 kDa) 5 kDa MWCO RC	23 min	94%	37 min	89%
α-Chymotrypsin 1.0 mg/mL (25 kDa) 10 kDa MWCO RC	7 min	93%	9 min	92%
BSA 1.0 mg/mL (66 kDa) 10 kDa MWCO RC	8 min	94%	10 min	98%
30 kDa MWCO RC	4 min	96%	4 min	93%
IgG 1.0 mg/mL (160 kDa) 50 kDa MWCO RC	17 min	95%	11 min	96%
100 kDa MWCO RC	18 min	89%	12 min	89%

Ordering information

Vivaspin® Turbo 15 RC	Pack size	Prod. no.
5 kDa MWCO	12	VS15TR11
5 kDa MWCO	48	VS15TR12
10 kDa MWCO	12	VS15TR01
10 kDa MWCO	48	VS15TR02
30 kDa MWCO	12	VS15TR21
30 kDa MWCO	48	VS15TR22
50 kDa MWCO	12	VS15TR31
50 kDa MWCO	48	VS15TR32
100 kDa MWCO	12	VS15TR41
100 kDa MWCO	48	VS15TR42

* 3,000 g for 100 kDa MWCO devices

Vivaspin® 20



5 to 20 mL samples

Vivaspin® 20 centrifugal concentrators have been developed to offer increased volume flexibility and performance.

Vivaspin® 20 handles up to 20 mL in swing bucket centrifuges and 14 mL in 25° fixed angle rotors accepting 50 mL centrifuge tubes.

Featuring twin vertical membranes for unparalleled filtration speeds the Vivaspin® 20 can achieve 100× plus concentrations.


Remaining volume is easy to read off the printed scale on the side of the concentrator and the modified dead stop pocket further simplifies direct pipette recovery of the final concentrate.

More process flexibility

Vivaspin® 20 is available with unique accessories and operating methods that are designed to provide more process flexibility and further time saving.


Gas pressure ultrafiltration

When an appropriate centrifuge is unavailable, or for single sample processing, Vivaspin® 20 can be filled with up to 15 mL and then pressurised for bench top concentration. For even faster processing, gas pressure can be combined with centrifugal force. “Pressure-fugation” is particularly suitable for difficult or viscous samples such as serum, or when using a low process temperature which reduces filtration speed, and generally when minimum process time is essential.



Technical specifications

Concentrator capacity	
Swing bucket rotor	20 mL
Fixed angle rotor	14 mL
With pressure head	15 mL
Dimensions	
Total length	116 mm 125 mm with pressure head
Width	30 mm
Active membrane area	6.0 cm ²
Hold-up volume, membrane	< 20 µL
Dead-stop volume	50 µL
Materials of construction	
Body	Polycarbonate
Filtrate vessel	Polycarbonate
Concentrator cap	Polypropylene
Pressure head	Acetal aluminium
Membrane	Polyethersulfone



Using the Vivaspin® 20 pressure cap

Constant volume diafiltration


In this procedure following concentration, a diafiltration cup is filled with buffer and then spun one time to achieve 98% salt removal. This compares to the need for two spins to achieve the same result with the traditional refill and re-spin procedure.

The improved performance is due to the constant washing action of the buffer solution from the diafiltration cup as it replaces solvent and salts as they pass through the ultrafiltration membrane.

Equipment required

Centrifuge		
Rotor type	Swing bucket	Fixed angle
Minimum rotor angle	-	25°
Rotor cavity	To fit 50 mL (30 mm) conical bottom tubes	To fit 50 mL (30 mm) conical bottom tubes
Maximum speed	4,000 g*	6,000 g*
Maximum speed, ≥100 kDa MWCO	3,000 g	6,000 g
Pressure		
Pressure accessories	VCA002, VCA005 and VCA200	
Maximum pressure	5 bar (75 psi)	
Rotor type for pressure-fuge	Swing bucket	
Maximum speed, pressure-fuge	3,000 g	
Maximum speed, pressure-fuge, ≥100 kDa MWCO	2,000 g	
Concentrate recovery		
Pipette type	Fixed or variable volume	Fixed or variable volume
Recommended tip	Thin gel loader type	Thin gel loader type

* When pressurized, maximum speed is 3,000 g (up to 50 kDa MWCO) or 2,000 g (≥100 kDa MWCO)



Desalting of concentrated sample

Performance characteristics

	Time to concentrate up to 30× at 20°C and solute recovery							
	Centrifuge		Centrifuge		Bench top		Press-fuge	
Mode	Swing bucket		25° Fixed angle		Pressure		Swing bucket	
Rotor	3,000 g		6,000 g		4 bar		3,000 g + 4 bar	
Centrifugal force pressure	20 mL		14 mL		10 mL		10 mL	
Start volume	Min.	Rec.	Min.	Rec.	Min.	Rec.	Min.	Rec.
Cytochrome c 0.25 mg/mL (12.4 kDa)								
3 kDa MWCO PES	110	97%	180	96%	60	96%	-	-
BSA 1.0 mg/mL (66 kDa)								
5 kDa MWCO PES	23	99%	29	99%	50	98%	14	98%
10 kDa MWCO PES	16	98%	17	98%	32	97%	8	97%
30 kDa MWCO PES	13	98%	15	98%	32	97%	8	97%
IgG 0.25 mg/mL (160 kDa)								
30 kDa MWCO PES	27	97%	20	95%	46	94%	13	97%
50 kDa MWCO PES	27	96%	22	95%	46	93%	13	96%
100 kDa MWCO PES	25	91%	20	90%	42	88%	12	94%
Latex beads 0.004% in DMEM +10% FCS (55 nm)								
300 kDa MWCO PES	20	99%	35	99%	10	99%	-	-
Latex beads 0.004% in DMEM +10% FCS (240 nm)								
1,000 kDa MWCO PES	4	99%	12	99%	4	99%	-	-
Yeast 1.0 mg/mL (<i>S. Cerevisiae</i>) 0.2 µm PES								
0.2 µm PES	15	95%	5	95%	20	95%	2	95%

Ordering information

Vivaspin® 20 Polyethersulfone	Pack size	Prod. no.
3 kDa MWCO	12	VS2091
3 kDa MWCO	48	VS2092
5 kDa MWCO	12	VS2011
5 kDa MWCO	48	VS2012
10 kDa MWCO	12	VS2001
10 kDa MWCO	48	VS2002
10 kDa MWCO (IVD device)	12	VS2001IVD*
10 kDa MWCO (IVD device)	48	VS2002IVD*
30 kDa MWCO	12	VS2021
30 kDa MWCO	48	VS2022
50 kDa MWCO	12	VS2031
50 kDa MWCO	48	VS2032
100 kDa MWCO	12	VS2041
100 kDa MWCO	48	VS2042
300 kDa MWCO	12	VS2051
300 kDa MWCO	48	VS2052
1,000 kDa MWCO	12	VS2061
1,000 kDa MWCO	48	VS2062
0.2 µm	12	VS2071
0.2 µm	48	VS2072

Vivaspin® 20 accessories

Air pressure controller (APC) fitted with pressure gauge, regulator, over-pressure safety valve and female coupling. APC is supplied with extension line (4 mm pneumatic tubing, 1 m) with male and female couplings, and inlet tubing (6 mm pneumatic tubing, 1 m)	1	VCA002
Charge valve for pressure head VCA200	1	VCA005
Female coupling	1	VCA010
Male coupling	1	VCA011
Replacement extension line (4 mm pneumatic tubing, 3 m)	1	VCA012
Vivaspin® 20 pressure head	1	VCA200
Diafiltration cups	12	VSA005

* IVD device article codes available only in countries where registered according to local IVD regulations

Vivaclear Centrifugal Filters



Vivaclear centrifugal filters are disposable microfiltration devices for the fast and reliable clarification or filtration of biological samples in the range 100 to 500 μL . They can be used in fixed angle rotors accepting 2.2 mL centrifuge tubes.

Applications

- Clarification of samples before loading onto Vivapure[®] protein purification spin columns
- Removal of particles and particulates
- Filtration of plasma and serum
- Removal of cells or cell debris

Product Features

- High-flux polyethersulfone membrane
- 0.8 μm pore size
- Low hold-up volume (<5 μL)
- Fast and reproducible performance


Technical specifications Vivaclear Centrifugal Filters

Rotor	40–45° Fixed angle rotor 500 μL
Pore size	0.8 μm
Dimensions	
Total length	43 mm
Diameter, filtrate tube	11 mm
Active membrane area	0.34 cm^2
Hold-up volume, membrane and support	< 5 μL
Maximum RCF	2,000 g
Materials of construction	
Body	Polypropylene
Filtrate tube	Polypropylene
Membrane	Polyethersulfone

Ordering information

Vivaclear Mini	Pack size	Prod. no.
0.8 μm PES	100	VK01P042

Vivacell 100



20 to 98 mL samples

Vivacell 100 bridges the gap between centrifugal concentrators and cross-flow cassettes. These devices feature vertical membranes for high speed processing of even high particle loaded samples. In addition, a unique choice between centrifugal, pressure or pressure-shake operating methods provides unrivaled process flexibility.

Fitting swing bucket rotors accepting 250 mL bottles, Vivacell 100 offers the highest sample capacity available in a centrifugal device – up to an astonishing 90 mL.

Vivacell 100 units can also be used for single or extremely sensitive samples of up to 98 mL when pressurized and left on the bench. For temperature-sensitive samples, they can be placed into a refrigerator when pressurized. Pressurization is made easy by use of quick-release connectors and can be combined with orbital shaking for even faster sample concentration. In whichever mode Vivacell 100 is used, the vertical membrane design inhibits membrane fouling while the integrated dead-stop impedes concentration to dryness and loss of sample.

Technical specifications

Concentrator capacity

Swing bucket rotor	90 mL
With pressure head	98 mL

Dimensions


Total length	123 mm centrifugal 197 mm with pressure head
Width	62 mm
Active membrane area	23.5 cm ²
Hold-up volume of membrane	< 250 µL
Dead-stop volume	350 µL

Operating requirements


Rotor type	Swing bucket
Rotor cavity	To fit 250 mL (62 mm) centrifuge bottles (maximum cavity depth 105 mm)
Maximum speed	2,000 g
Maximum pressure	5 bar (75 psi)

Materials of construction

Body	Polycarbonate
Filtrate vessel	Polycarbonate
Concentrator cap	Polypropylene
Pressure head	Acetal
Pressure head seal	Thermoplastic Elastomer Vulcanizate
Membrane	Polyethersulfone




Device fits standard 250 mL rotors




Centrifuge

- Process convenience
- Low shear, no foaming
- Less visual control



Pressure

- Simplicity and highest process control
- Ideal for refrigerated use
- Slower concentrations



Pressure-shake

- Speed and process control
- Ideal for single samples

Performance characteristics

90 mL start volume	Time to concentrate up to 30× at 20°C			Solute recovery
	Swing bucket, 2,000 g	Pressure, 4 bar (60 psi)		
		Static	Orbital shake	
BSA 1.0 mg/mL (66 kDa)				
5 kDa MWCO PES	22 min	75 min	25 min	96%
10 kDa MWCO PES	16 min	60 min	20 min	96%
30 kDa MWCO PES	16 min	60 min	20 min	94%
IgG 0.25 mg/mL (160 kDa)				
50 kDa MWCO PES	20 min	70 min	30 min	94%
100 kDa MWCO PES	20 min	85 min	30 min	90%
Latex beads 0.004% in DMEM + 10% FCS (55 nm)				
300 kDa MWCO PES	35 min	-	120 min	99%
Latex beads 0.004% in DMEM + 10% FCS (240 nm)				
1,000 kDa MWCO* PES	4 min	5 min	4 min	99%

*2,000 g in centrifuge, 2 bar (29 psi) pressure



Ordering information

Vivacell 100 PES with PP cap	Pack size	Prod. no.
5 kDa MWCO	2	VC1011
5 kDa MWCO	10	VC1012
10 kDa MWCO	2	VC1001
10 kDa MWCO	10	VC1002
30 kDa MWCO	2	VC1021
30 kDa MWCO	10	VC1022
50 kDa MWCO	2	VC1031
50 kDa MWCO	10	VC1032
100 kDa MWCO	2	VC1041
100 kDa MWCO	10	VC1042
300 kDa MWCO	2	VC1051
300 kDa MWCO	10	VC1052
1,000 kDa MWCO	2	VC1061
1,000 kDa MWCO	10	VC1062
0.2 µm	2	VC1071
0.2 µm	10	VC1072


Vivacell 100 accessories

Air pressure controller (APC) fitted with pressure gauge, regulator, over-pressure safety valve and female coupling. APC is supplied with extension line (4 mm pneumatic tubing, 1 m) with male and female couplings, and inlet tubing (6 mm pneumatic tubing, 1 m)	1	VCA002
Female coupling	1	VCA010
Male coupling	1	VCA011
Replacement extension line (4 mm pneumatic tubing, 3 m)	1	VCA012
Replacement pressure head seals	10	VCA014
Vivacell 100 pressure head with seals (5×)	1	VCA800

Vivaflow® 50



Multiple cassettes



Single cassette

0.1 to 5 L

The novel Vivaflow® 50 system provides a standard of ease of use, performance, flexibility and economy which is unrivalled by any laboratory or pilot scale filtration system on the market.

Unique features

- Thin channel flip-flow recirculation path provides high cross flow velocities with minimum pump requirements.
- No need for pressure holders.
- Crystal clear for simple control of remaining hold up and membrane status.
- Unique Interlocking modules with series connectors for easy scale up.
- Disposable | single use.

Unique performance

- A single 50 cm² module will typically reduce 500 mL to less than 15 mL in under 50 minutes.
- Less than 10 mL minimum system recirculation for highest concentrations.
- Less than 500 µL non recoverable hold up volume.
- Near total recoveries achievable with a single 10 mL rinse.

Unique “flip-flow” thin channel flow path results in high turbulence and linear velocity for exceptional flux even at high concentrations

Technical specifications Vivaflow® 50

Dimensions

Overall L W H	25 107 84 mm
Channel W H	15 mm 0.3 mm
Active membrane area	50 cm ²
Minimum recirculation volume	< 10 mL
Hold-up volume, cassette	1.5 mL
Non recoverable hold-up	< 0.5 mL

Operating conditions

Pump flow rate	200 – 400 mL/min
Maximum pressure	3 bar (45 psi)
Maximum temperature	60°C

Materials of construction

Main housing	Polycarbonate
Flow channel	TPX (PMP)
Membrane support	TPX (PMP)
Seals and O rings	Silicone
Pressure indicator	Polypropylene, SS spring
Flow restrictor	Polypropylene
Fittings	Nylon
Tubing	PVC (medical grade)



Performance characteristics

	Time to concentrate up to 20x at 3 bar inlet pressure, 20°C			
	Single device 250 mL start volume	Three devices 1 L start volume	Solute recovery	
			Direct	10 mL rinse
BSA 1.0 mg/mL (66 kDa)				
5 kDa MWCO PES	34 min	49 min	96%	>99%
10 kDa MWCO PES	22 min	32 min	94%	>99%
30 kDa MWCO PES	22 min	32 min	92%	99%
50 kDa MWCO PES	20 min	29 min	92%	98%
IgG 1.0 mg/mL (160 kDa)				
100 kDa MWCO PES	43 min	62 min	92%	98%
100 kDa MWCO RC	40 min	58 min	92%	98%
Yeast 1.0 mg/mL (<i>S. Cerevisiae</i>)				
0.2 µm PES	33 min	47 min	92%	98%


Ordering information

Vivaflow® 50 modules include filtrate tube, size 16 peristaltic tubing, flow restrictor and fittings	Pack size	Prod. no.
3 kDa MWCO PES	2	VF05P9
5 kDa MWCO PES	2	VF05P1
10 kDa MWCO PES	2	VF05P0
30 kDa MWCO PES	2	VF05P2
50 kDa MWCO PES	2	VF05P3
100 kDa MWCO PES	2	VF05P4
1,000 kDa MWCO PES	2	VF05P6
0.2 µm PES	2	VF05P7
100 kDa MWCO RC	2	VF05C4
Vivaflow® 50 complete system		
Pump (240 V), Easy-Load pump head (size 16), tubing, 500 mL sample diafiltration reservoir, cassette stand, pressure indicator, T-connectors, series interconnectors	1	VFS502
Pump (115 V), Easy-Load pump head (size 16), tubing, 500 mL sample diafiltration reservoir, cassette stand, pressure indicator, T-connectors, series interconnectors	1	VFS504
Vivaflow® 50 tubing and fittings		
Size 16 PVC pump tubing (3 metres, 3.2 × 1.6 mm)	1	VFA004
Flow restrictor set (2 each of 0.4, 0.6, 0.8 mm)	6	VFA009
T-connectors for running 2 stacks	2	VFA030
Series interconnectors	6	VFA031
Female luer fittings	10	VFA032
VF50 tubing kit (2 × 1 m size 16 PVC tubing with inlet fittings, 2 × 50 cm size 16 PVC tubing with 0.6 mm flow restrictors, 1 × series interconnector)	1	VFA034
Flow restrictor 0.6 mm	6	VFA035


Vivaflow® 50 accessories

Masterflex Economy Drive variable speed peristaltic pump (240 V)	VFP001
Masterflex Economy Drive variable speed peristaltic pump (115 V)	VFP002
500 mL sample diafiltration reservoir	VFA006
Masterflex Easy-Load pump head - size 16	VFA012
Vivaflow® 50 stand	VFA016
Pressure indicator (1-3 bar)	VFA020

Vivaflow® 50R



Vivaflow® 50R - Single cassette



Vivaflow® 50R - Two cassettes

0.1 to 1 L samples

Concentrate 100 mL to under 20 mL in just a few minutes or concentrate one liter 50 times in less than 60 minutes. Alternatively, speed up your process by using two Vivaflow® 50R cassettes in parallel and concentrate 1 liters in under 30 min.

Vivaflow® 50R is a plug-and-play laboratory crossflow cassette for concentrating up to 1 L aqueous samples. The active membrane area per device is 50 cm².

One unit comes with all the necessary accessories for running the device with a laboratory pump and a size 16 pump head. For speeding up concentration, two cassettes can be run simultaneously.

- Fast and easy protein sample concentration
- Reusable
- Concentrates volumes from 0.1 L to 1 L
- Optimal for concentration of culture supernatants and viruses
- The most compact crossflow cassette with a premium Hydrosart® membrane

Technical specifications Vivaflow® 50R

Dimensions

Overall L W H	24 100 100 mm
Channel W H	7.5 0.4 mm
Active membrane area	50 cm ²
Minimum recirculation volume	10 mL
Hold-up volume, cassette	1.7 mL
Non-recoverable hold-up	< 0.5 mL

Operating conditions

Pump flow rate	200 - 400 mL/min
Maximum pressure	4 bar (60 psi)
Maximum temperature	60°C

Materials of construction

Main housing	Acrylic
Flow channel	Acrylic
Membrane support	Polypropylene
Seals and O-rings	Silicone
Pressure indicator	Polypropylene, SS spring
Flow restrictor	Polypropylene
Fittings	Nylon
Tubing	PVC (medical grade)

Visit us at

www.sartorius.com/

Vivaflow50R

for further information.

Here you can find

instructions on how to

use Vivaflow® 50R for:

- Preparation of biological nanoparticles and medical nanocarriers
- Concentration and purification of viruses

Performance characteristics

	Time to concentrate up to 20× at 3.0 bar inlet 2.5 bar outlet pressure, 20°C			
	Start volume 250 mL	Average flux mL/min	Recovery	
			Direct	25 mL rinse
Lysozyme 0.25 mg/mL (14.3 kDa)				
5 kDa MWCO Hydrosart®	70	3.4	96%	98%
10 kDa MWCO Hydrosart®	23	10.3	94%	96%
BSA 1.0 mg/mL (66 kDa)				
10 kDa MWCO Hydrosart®	24	9.9	98%	>99
30 kDa MWCO Hydrosart®	15	15.8	97%	>99
IgG 1.0 mg/mL (160 kDa)				
100 kDa MWCO Hydrosart®	46	5.2	97%	>99
Start volume 1 L (one Vivaflow® 50R at 3 bar), BSA 1.0 mg/mL				
10 kDa MWCO Hydrosart®	95	10.0	98%	>99
Start volume 1 L (two Vivaflow® 50R in parallel at 3 bar), BSA 1.0 mg/mL				
10 kDa MWCO Hydrosart®	48	19.8	98%	>99

Ordering information


Vivaflow® 50R*	Pack size	Prod. no.
5 kDa MWCO Hydrosart®	1	VF05H1
10 kDa MWCO Hydrosart®	1	VF05H0
30 kDa MWCO Hydrosart®	1	VF05H2
100 kDa MWCO Hydrosart®	1	VF05H4
Vivaflow® 50R complete system		
Pump (230 V), Easy-Load pump head (size 16), tubing, 500 mL sample diafiltration reservoir	1	VFS202
Pump (115 V), Easy-Load pump head (size 16), tubing, 500 mL sample diafiltration reservoir	1	VFS204
Vivaflow® 50R tubing and fittings		
Size 16 PVC pump tubing and Luer fittings (3 m, 3.2 × 1.6 mm)	1	VFA004
T-connectors for running 2 stacks	2	VFA030
Flow restrictor set (2 each of 0.4 mm, 0.6 mm, 0.8 mm)	6	VFA009
Female luer fittings - size 16	10	VFA032
Flow restrictors 0.6 mm	6	VFA035
Female luer fittings - size 15	10	VFA036

Vivaflow® 50R accessories

Masterflex Economy Drive variable speed peristaltic pump (230 V)	1	VFP001
Masterflex Economy Drive variable speed peristaltic pump (115 V)	1	VFP002
500 mL sample diafiltration reservoir	1	VFA006
Masterflex Easy-Load pump head - size 16	1	VFA012

* Vivaflow® 50R modules include pressure indicator, flow restrictor and size 16 pvc peristaltic tubing and fittings.

Vivaflow® 200



Vivaflow® 200 set-up for diafiltration

0.5 to 5 L

Concentrate 250 mL to under 20 mL in just a few minutes or concentrate one litre 50 times in less than 30 minutes. Alternatively, use two Vivaflow® 200 cassettes in parallel and concentrate 5 litres in under 75 minutes.

Near total sample recoveries can be expected with most solutions.

The economical standard package comes complete with tubing, pressure indicator, flow restrictor and high pressure pump tubing. All you need is a peristaltic pump capable of handling 6.4 mm OD (size 16) tubing. Should your pump head require larger tubing, link your own peristaltic tube up to the standard product, using the interconnector provided.

Two modules in parallel will concentrate 5 litres in under 75 minutes.

Technical specifications Vivaflow® 200

Dimensions

Overall L | H | W 38 | 126 | 138 mm

Channel W | H 10 mm | 0.4 mm

Active membrane area 200 cm²

Minimum recirculation volume < 20 mL

Hold up volume, cassette 5.3 mL

Non-recoverable hold-up < 2 mL

Operating conditions

Pump flow 200–400 mL/min

Maximum pressure 4 bar (60 psi)

Maximum temperature 60°C

Materials of construction

Main housing Acrylic

Flow channel Acrylic

Membrane support Polypropylene

Seals and O rings Silicone

Pressure indicator Polypropylene, SS spring


Flow restrictor Polypropylene

Fittings Nylon


Tubing PVC (medical grade)

Performance characteristics

	Time to concentrate up to 20× at 3 bar inlet pressure, 20°C			
	1 litre start volume	Average flux mL/min	Recovery	
			Direct	25 mL rinse
Lysozyme 0.25 mg/mL (14.3 kDa)				
2 kDa MWCO Hydrosart®	160	6	97%	>99%
3 kDa MWCO PES	180	5	97%	>99%
BSA 1.0 mg/mL (66 kDa)				
5 kDa MWCO PES	29	33	98%	>99%
5 kDa MWCO Hydrosart®	70	14	98%	>99%
10 kDa MWCO PES	23	41	96%	>99%
10 kDa MWCO Hydrosart®	35	27	98%	>99%
30 kDa MWCO PES	25	38	96%	99%
30 kDa MWCO Hydrosart®	20	48	96%	>99%
50 kDa MWCO PES	22	43	96%	98%
IgG 1.0 mg/mL (160 kDa)				
100 kDa MWCO PES	54	18	96%	99%
Yeast 1.0 mg/mL (<i>S. Cerevisiae</i>)				
0.2 µm PES	11	86	92%	98%
Dilute solute concentration, start volume 1 litre at 3 bar, 10 kDa MWCO PES				
BSA 0.001 mg/mL	18	52	90%	98%
BSA 0.01 mg/mL	20	47	92%	98%
BSA 0.1 mg/mL	21	45	94%	99%
Start volume 5 litres (two VF200 in parallel at 3 bar) 10 kDa MWCO PES				
BSA 1.0 mg/mL (66 kDa)	67	70	97%	>99%



Operation – Single cassette



Operation – Two cassettes

Ordering information

Vivaflow® 200 modules include pressure indicator, flow restrictor and size 16 pvc peristaltic tubing and fittings	Pack size	Prod. no.
3 kDa MWCO PES	1	VF20P9
5 kDa MWCO PES	1	VF20P1
10 kDa MWCO PES	1	VF20P0
30 kDa MWCO PES	1	VF20P2
50 kDa MWCO PES	1	VF20P3
100 kDa MWCO PES	1	VF20P4
0.2 µm PES	1	VF20P7
2 kDa MWCO Hydrosart®	1	VF20H9
5 kDa MWCO Hydrosart®	1	VF20H1
10 kDa MWCO Hydrosart®	1	VF20H0
30 kDa MWCO Hydrosart®	1	VF20H2
100 kDa MWCO Hydrosart®	1	VF20H4
Vivaflow® 200 complete system		
Pump (240 V), Easy-Load pump head (size 16), tubing, 500 mL sample diafiltration reservoir	1	VFS202
Pump (115 V), Easy-Load pump head (size 16), tubing, 500 mL sample diafiltration reservoir	1	VFS204
Vivaflow® 200 tubing and fittings		
Size 15 PVC pump tubing and Luer fittings (3 m, 4.8 × 2.6 mm)	1	VFA003
Size 16 PVC pump tubing and Luer fittings (3 m, 3.2 × 1.6 mm)	1	VFA004
Y-connector (size 15 to 2 × size 16)	1	VFA005
Flow restrictor set (2 × 0.4, 0.6, 0.8 mm)	6	VFA009
Female luer fittings size 16	10	VFA032
Flow restrictors 0.6 mm	6	VFA035
Female luer fittings size 15	10	VFA036

Vivaflow® 200 accessories

Masterflex Economy Drive variable speed peristaltic pump (240 V)	VFPO01
Masterflex Economy Drive variable speed peristaltic pump (115 V)	VFPO02
500 mL sample diafiltration reservoir	VFA006
Masterflex Easy-Load head – size 16	VFA012
Masterflex Easy-Load pump head – size 15	VFA013

Vivapore® Solvent Absorption Concentrators

1 to 20 mL samples


With no need for additional equipment, pressure or vacuum, solvent absorption is the most economic and user friendly concentration technique available to the clinician and research scientist.

Just fill the unit with the solution to be concentrated, wait for the desired concentration level to be achieved and then pipette the concentrated sample from the bottom of the reservoir.


Vivapore® is ideal for general purpose laboratory concentration or purification prior to further analysis. It is particularly suited for labile solutions that can denature with alternative shear or pressure inducing methods or that require processing in a cold room environment.

Vivapore® concentrators extend the solvent absorption technique to a totally new level of performance, application potential and ease of use.

Vivapore® solvent absorption concentrators are IVD registered devices and only available in countries where they are registered according to local regulations. Please contact Sartorius for more information on registered countries and availability.



Vivapore® 5



Vivapore® 10 | 20

Technical specifications

Dimensions	Vivapore® 5	Vivapore® 10 20
Sample capacity	1 - 5 mL	2 - 10 mL 20 mL*
Overall W H	42 82 mm	46 100 mm
Active membrane area	20 cm ²	28 cm ²
Dead-stop volume	50 µL	50 µL
Maximum concentration factor	100×	400×
Materials of construction		
Concentrator	SAN	SAN
Expansion reservoir (optional)	-	SAN
Membrane	Polyethersulfone	Polyethersulfone

* up to 20 mL with optional 10 mL expansion reservoir (Prod. No. VPA006)

Performance characteristics

Product	Time to concentrate up to 10× [min.]			Concentrate recovery		
	VP5	VP10 20	VP10 20**	VP5	VP10 20	VP10 20*
Start volume	5 mL	10 mL	20 mL	5 mL	10 mL	20 mL
Cytochrome c (12.6 kDa) [mg/mL]	0.25	0.25	0.1	0.25	0.25	0.1
7.5 kDa MWCO PES	35	75	150	90%	90%	92%
BSA (66 kDa)						
7.5 kDa MWCO PES	30	55	115	92%	92%	92%
IgG (160 kDa MW)						
7.5 kDa MWCO PES	40	70	160	75%	77%	78%

Product	Time to concentrate up to 50× [min.]			Concentrate recovery		
	VP5	VP10 20	VP10 20*	VP5	VP10 20	VP10 20*
Start volume	5 mL	10 mL	20 mL	5 mL	10 mL	20 mL
Cytochrome c (12.6 kDa)	65	70	160	91%	88%	90%
7.5 kDa MWCO PES						
BSA (66 kDa)	45	50	105	90%	90%	92%
7.5 kDa MWCO PES						
IgG (160 kDa)	50	65	140	53%	65%	74%
7.5 kDa MWCO PES						

Ordering information

Vivapore® 5	Pack size	Prod. no.
7.5 kDa MWCO PES, including stand and recovery pipettes	4	VP0503**
7.5 kDa MWCO PES, including stand and recovery pipettes	30	VP0501**
7.5 kDa MWCO PES	100	VP0502**
Vivapore® 10 20		
7.5 kDa MWCO PES, including stand and recovery pipettes	4	VP2003**
7.5 kDa MWCO PES including stand and recovery pipettes	30	VP2001**
7.5 kDa MWCO PES	100	VP2002**


Vivapore® accessories

Disposable stand for 4 units	6	VPA002
Plastic recovery pipettes for Vivapore® 5	100	VPA007
10 mL expansion reservoir for Vivapore® 10 20	10	VPA006
Plastic recovery pipettes for Vivapore® 10 20	100	VPA005

* with expansion reservoir

** IVD devices are available only in countries where registered according to local IVD regulations

Ultrafiltration Membrane Discs



Polyethersulfone (PES)

This is a general purpose membrane that provides excellent performance with most solutions when high recoveries in the retentate are of primary importance. Polyethersulfone membranes exhibit no hydrophobic or hydrophilic interactions and are usually preferred for their low fouling characteristics, exceptional flux and broad pH compatibility.

Cellulose Triacetate (CTA)

High hydrophilicity and very low non-specific binding characterize this membrane. Cast without any mem-

brane support that could trap or bind passing microsolute, these membranes are preferred for sample cleaning, protein removal, and when high recoveries from the filtrate are of primary importance.

Regenerated Cellulose (RC)

These membranes are also highly hydrophilic and are often preferred for their higher protein recovery when processing very dilute solutions. Resistance to autoclaving, ease of cleaning and extended chemical compatibility also characterize this membrane material.

Typical Performance and Specifications

	Polyethersulfone, Type 146	Cellulose Triacetate, Type 145	Regenerated Cellulose, Type 144
Thickness	120 μm	120 μm	180 μm
pH range	1 - 14	4 - 8	1 - 13
Flux with water			
10 kDa MWCO	0.2 mL/min/cm ²	0.11 mL/min/cm ²	0.08 mL/min/cm ²
Protein retention			
Cytochrome c	95%	90%	99%

Ordering information

PES Membrane Discs, Type 146	Diameter	Pack size	Prod. no.
1 kDa MWCO	47 mm	10	14609--47-----D
5 kDa MWCO	25 mm	10	14629--25-----D
	47 mm	10	14629--47-----D
	63 mm	10	14629--63-----D
	76 mm	10	14629--76-----D
10 kDa MWCO	25 mm	10	14639--25-----D
	47 mm	10	14639--47-----D
	63 mm	10	14639--63-----D
	76 mm	10	14639--76-----D
30 kDa MWCO	47 mm	10	14659--47-----D
	63 mm	10	14659--63-----D
50 kDa MWCO	47 mm	10	14650--47-----D
100 kDa MWCO	47 mm	10	14668--47-----D
	63 mm	10	14668--63-----D
300 kDa MWCO	47 mm	10	14679--47-----D
CTA Membrane Discs, Type 145			
5 kDa MWCO	47 mm	10	14529--47-----D
10 kDa MWCO	47 mm	10	14539--47-----D
	50 mm	10	14539--50-----D
20 kDa MWCO	43 mm	10	14549--43-----D
	47 mm	10	14549--47-----D
	47 mm	10	14549--47-----N
RC Membrane Discs, Type 144			
5 kDa MWCO	25 mm	10	14429--25-----D
	47 mm	10	14429--47-----D
	63 mm	10	14429--63-----D
	76 mm	10	14429--76-----D
10 kDa MWCO	25 mm	10	14439--25-----D
	47 mm	10	14439--47-----D
	63 mm	10	14439--63-----D
	76 mm	10	14439--76-----D
30 kDa MWCO	47 mm	10	14459--47-----D
	63 mm	10	14459--63-----D



DNA Concentration

Table of Contents

Vivacon® 500	50
Vivacon® 2	53

Vivacon® 500



100 to 500 µL samples

Vivacon® 500 centrifugal concentrators offer the optimal solution for DNA and protein concentration and buffer exchange applications. For optimal performance with very dilute samples, e.g. forensic samples, Vivacon® 500 is equipped with the patented regenerated cellulose membrane, Hydrosart®.

High recoveries and excellent reproducibilities are paired with convenience offered by molecular weight cut-off printed on individual devices.



The possibility of a reverse-spin after sample processing assures complete concentrate recovery which is especially important when working with low sample concentrations.

Vivacon® 500 PCR Grade

When using DNA amplification technologies, any traces of DNA originating from the equipment have to be eliminated.

Vivacon® 500 PCR Grade units are treated with ethylene oxide (ETO) in a validated process in order to deactivate all traces of DNA that might interfere with subsequent amplification procedures.

Ref.: K. Shaw et al., Int. J. Legal Med. (2008) 122: 29–33



Reverse spin

Feature	Benefit
Reverse-spin enabled	Complete and highly reproducible sample recovery
Low binding material	High recoveries of low sample concentrations

Technical Specifications Vivacon® 500

Concentrator capacity	
Fixed angle rotor	0.5 mL
Dimensions	
Total length, concentration	45 mm
Total length, reverse-spin	47.5 mm
Width	12.4 mm
Active membrane area	0.32 cm ²
Hold-up volume, of membrane and support	< 5 µL
Dead-stop volume (40° rotor)	5 µL
Materials of construction	
Body	Polycarbonate
Filtrate vessel	Polypropylene
Membrane	Hydrosart®

Conversion table for MWCO to Nucleotide Cut-Off

Membrane	MWCO	Double-Stranded Nucleotide Cut-Off
Hydrosart®	2 kDa	> 10 bp
Hydrosart®	10 kDa	> 30 bp
Hydrosart®	30 kDa	> 50 bp
Hydrosart®	50 kDa	> 300 bp
Hydrosart®	100 kDa	> 600 bp
Cellulose Acetate	125 kDa	> 650 bp

Performance characteristics for DNA

Start volume 0.5 mL, sample concentration 50 ng/mL

MWCO	Sample size	Time to concentrate up to 30× at 20°C	Concentrate recovery	RCF
2 kDa	10 bp	60 min	93%	7,500 g
10 kDa	30 bp	25 min	94%	7,500 g
30 kDa	50 bp	18 min	88%	5,000 g
50 kDa	300 bp	18 min	91%	5,000 g
100 kDa	600 bp	10 min	87%	3,000 g
125 kDa	650 bp	12 min	85%	2,000 g
125 kDa	900 bp	9 min	94%	3,000 g

Performance characteristics for proteins

Start volume 0.5 mL, sample and concentration of proteins as specified in table

MWCO	Sample	Time to concentrate up to 30× at 20°C	Concentrate recovery	RCF
2 kDa	0.25 mg/mL cytochrome c	30 min	95%	14,000 g
10 kDa	0.25 mg/mL cytochrome c	15 min	92%	14,000 g
30 kDa	1.0 mg/mL BSA	10 min	95%	14,000 g
50 kDa	1.0 mg/mL BSA	10 min	92%	14,000 g
100 kDa	1.0 mg/mL bovine IgG	11 min	90%	8,000 g
125 kDa	1.0 mg/mL bovine IgG	10 min	81%	8,000 g



Ordering information

Vivacon® 500	Pack size	Prod. No.
2 kDa MWCO	25	VN01H91
2 kDa MWCO	100	VN01H92
10 kDa MWCO	25	VN01H01
10 kDa MWCO	100	VN01H02
30 kDa MWCO	25	VN01H21
30 kDa MWCO	100	VN01H22
30 kDa MWCO	500	VN01H23
50 kDa MWCO	25	VN01H31
50 kDa MWCO	100	VN01H32
50 kDa MWCO	500	VN01H33
100 kDa MWCO	25	VN01H41
100 kDa MWCO	100	VN01H42
100 kDa MWCO	500	VN01H43
125 kDa MWCO	25	VN01H81
125 kDa MWCO	100	VN01H82
125 kDa MWCO	500	VN01H83
Vivacon® 500 PCR Grade		
30 kDa MWCO	100	VN01H22ETO
30 kDa MWCO	500	VN01H23ETO
100 kDa MWCO	100	VN01H42ETO
125 kDa MWCO	100	VN01H82ETO
125 kDa MWCO	500	VN01H83ETO

Vivacon® 500 accessories

Additional collection tubes	100	VNCT01
-----------------------------	-----	--------

Vivacon® 2



Reverse spin

0.4 – 2 mL samples

Vivacon® 2 centrifugal concentrators offer the optimal solution for DNA and protein concentration and buffer exchange applications. For optimal performance with very dilute samples, e.g. forensic samples, Vivacon® 2 is equipped with the patented regenerated cellulose membrane Hydrosart®.

High recoveries and excellent reproducibilities are paired with convenience offered by volume graduation and molecular weight cut-off printed on individual devices.

The possibility of a re-spin after sample processing assures complete concentrate recovery which is especially important when working with low sample concentrations.

Vivacon® 2 PCR Grade

Vivacon® 2 PCR Grade units are treated with ethylene oxide (ETO) in a validated process in order to deactivate all traces of DNA that might interfere with subsequent amplification procedures.

Feature	Benefit
Re-spin possibility	Complete and highly reproducible sample recovery
Low binding material	High recoveries of low sample concentration
Easy to remove re-spin cap	Convenient sample handling
Graduation printed on	Optimal process control

Technical Specifications Vivacon® 2

Concentrator capacity	
Fixed angle rotor	2 mL
Dimensions	
Total length, concentration	125 mm
Total length, reverse spin	115 mm
Diameter	16 mm
Active membrane area	0.95 cm ²
Hold-up volume, membrane and support	10 µL
Dead-stop volume (25° rotor)	55 µL
Materials of construction	
Body	Polycarbonate
Filtrate vessel	Polypropylene
Reverse-spin recovery vial	Polypropylene
Recovery vial cap	Polypropylene
Membrane	Hydrosart®

Conversion table for MWCO to Nucleotide Cut-Off

Membrane	MWCO	Double-Stranded Nucleotide Cut-Off
Hydrosart®	2 kDa	> 10 bp
Hydrosart®	10 kDa	> 30 bp
Hydrosart®	30 kDa	> 50 bp
Hydrosart®	50 kDa	> 300 bp
Hydrosart®	100 kDa	> 600 bp
Cellulose acetate	125 kDa	> 650 bp

Performance characteristics for DNA

Volume 2 mL, sample concentration 50 ng/mL, start volume: 2 mL

MWCO	Sample size	Time to concentrate up to 30× at 20°C	Concentrate recovery	RCF
2 kDa	10 bp	120 min	92%	7,500 g
10 kDa	30 bp	60 min	94%	5,000 g
30 kDa	50 bp	60 min	95%	2,500 g
50 kDa	300 bp	45 min	96%	2,500 g
100 kDa	600 bp	30 min	93%	2,500 g
125 kDa	650 bp	30 min	88%	2,500 g
125 kDa	900 bp	30 min	89%	2,500 g

Performance characteristics for proteins

Start volume 2 mL, sample and concentration of proteins as specified in table

MWCO	Test molecule	Time to concentrate up to 30× at 20°C	Concentrate recovery	RCF
2 kDa	0.25 mg/mL cytochrome c	120 min	95%	7,500 g
10 kDa	0.25 mg/mL cytochrome c	90 min	96%	5,000 g
30 kDa	1.0 mg/mL BSA	40 min	96%	5,000 g
50 kDa	1.0 mg/mL BSA	30 min	94%	5,000 g
100 kDa	1.0 mg/mL bovine IgG	30 min	92%	5,000 g
125 kDa	1.0 mg/mL bovine IgG	27 min	81%	5,000 g

Ordering information

Vivacon® 2	Pack size	Prod. No.
2 kDa MWCO	25	VN02H91
2 kDa MWCO	100	VN02H92
10 kDa MWCO	25	VN02H01
10 kDa MWCO	100	VN02H02
30 kDa MWCO	25	VN02H21
30 kDa MWCO	100	VN02H22
30 kDa MWCO	500	VN02H23
50 kDa MWCO	25	VN02H31
50 kDa MWCO	100	VN02H32
100 kDa MWCO	25	VN02H41
100 kDa MWCO	100	VN02H42
100 kDa MWCO	500	VN02H43
125 kDa MWCO	25	VN02H81
125 kDa MWCO	100	VN02H82
125 kDa MWCO	500	VN02H83
Vivacon® 2 PCR Grade		
30 kDa MWCO	100	VN02H22ETO
50 kDa MWCO	100	VN02H32ETO
100 kDa MWCO	100	VN02H42ETO
100 kDa MWCO	500	VN02H43ETO
125 kDa MWCO	500	VN02H83ETO




Protein Purification

Table of Contents

Vivapure® Ion Exchange Protein Purification Products

58

Vivapure® Ion Exchange Protein Purification Products



SEM comparing chromatography beads (right) with the Sartobind® membrane adsorber, which features 50x larger pore sizes.


Fast and easy-to-use spin columns

Vivapure® Ion Exchange (IEX) devices incorporate Sartobind® membrane adsorber technology as the chromatography matrix. The ready-to-use spin column format makes protein purification as easy as filtration. With no risk of running dry, Vivapure® replaces time-consuming and expensive resin-based chromatography in many protein purification workflows.

The rapid bind-wash-elute protocol is especially ideal in screening applications, where multiple samples or purification conditions can be conveniently processed in parallel.

The microporous structure of these membrane adsorbers has a pore size $> 3 \mu\text{m}$, which is orders of magnitude larger than conventional chromatography resins. This allows molecules to be transported to the ligands immobilized on the membrane adsorber by convective flow, overcoming the diffusion limitations of chromatography resins, and leading to very high flow rates. The large pore sizes also prevent gel filtration effects and minimize non-specific binding.

With Vivapure®, there is no need for column packing, saving time and ensuring reproducibility. Furthermore, Sartobind® membrane adsorber technology is available in process scale formats, making Vivapure® an indispensable tool for process development prior to purification scale-up.




Fast and easy protein purification with Vivapure® spin columns




Vivapure® Mini H



Vivapure® Maxi H



Vivapure® Mini H



Vivapure® Maxi H

Membrane availability

Functional groups	Ion exchanger type
Sulphonic acid (S)	Strong acidic cation exchanger: R-CH ₂ -SO ₃ ⁻ Na ⁺
Quaternary ammonium (Q)	Strong basic anion exchanger: R-CH ₂ -N ⁺ -(CH ₃) ₃ Cl ⁻
Diethylamine (D)	Weak basic anion exchanger: R-CH ₂ -NH ⁺ -(CH ₂ H ₅) ₂

Performance characteristics

Vivapure® spin columns	Protein binding capacity*	Max. volume, swing bucket	Max. volume, fixed angle
Vivapure® Mini H	4 mg	-	0.4 mL
Vivapure® Maxi H	60 – 80 mg	19 mL	10.5 mL

Typical applications


- Fractionation of protein mixtures prior to 1D or 2D-PAGE
- Scouting purification conditions for new protein targets
- Removal of endotoxins from monoclonal antibodies
- Preparation of heme moiety from heme containing protein prior to functional analysis
- General protein purification and polishing
- Detergent removal from protein solutions
- Purification of antibodies from serum, ascites or cell culture supernatant
- Intermediate sample purification prior to further HPLC | FPLC
- Purification of membrane-bound proteins

Detailed application notes are available on our website: www.sartorius.com

Ordering information

Vivapure® Mini Ion Exchange Spin Columns (up to 0.5 mL)	Spin Columns	Centrifuge Tubes	Prod. no.
Vivapure® Mini H starter kit (8 of each S and Q)	16	32	VS-IX01SQ16
Vivapure® D Mini H	24	48	VS-IX01DH24
Vivapure® Q Mini H	24	48	VS-IX01QH24
Vivapure® S Mini H	24	48	VS-IX01SH24
Vivapure® Maxi Ion Exchange Spin Columns (up to 20 mL)			
Vivapure® D Maxi H	8	16	VS-IX20DH08
Vivapure® Q Maxi H	8	16	VS-IX20QH08
Vivapure® S Maxi H	8	16	VS-IX20SH08

* Binding capacities established using 1 mg/mL BSA in 25 mM Tris-HCl (pH 8) for Vivapure® Q and D devices, or 1 mg/mL cytochrome c in 25 mM sodium acetate buffer (pH 5.5) for Vivapure® S devices. Actual capacities depend on the target molecule and selected buffer conditions.



Virus Purification and Concentration

Table of Contents

Vivapure® Virus Purification and Concentration Kits	62
Adenovirus Purification with Vivapure® Adenopack Kits	63
Vivapure® Adenopack 20	64
Vivapure® Adenopack 100	65
Vivapure® Adenopack 500	67
Lentivirus Purification with Vivapure® Lentiselect Kit	68
Vivapure® Lentiselect 40	69
Vivapure® Lentiselect 500	70
Vivapure® Lentiselect 1000	71

Vivapure® Virus Purification and Concentration Kits



Recombinant virus vectors are the preferred method for a wide range of gene delivery applications. Especially adenovirus type 5 and VSV-G pseudotyped lentivirus are two frequently utilized viral vectors for in vitro and in vivo applications.

Recombinant adenovirus vectors

Recombinant adenovirus vectors are versatile tools in research and therapeutic applications for gene transfer and protein expression in cell lines that have low transfection efficiency with liposomes.

After entering cells, the virus remains epichromosomal (i.e. does not integrate into the host chromosome, leaving the host genome unaffected). The delivery of RNAi into cells is becoming a major application for adenovirus vectors.

Lentivirus vectors

Lentivirus vectors are frequently used in gene transfer studies, due to their ability of gene transfer and integration into dividing and non-dividing cells. The pseudotyped envelope with vesicular stomatitis virus envelope G (VSV-G) protein broadens their target cell range. Lentiviral vectors have been shown to deliver genes into cell types (e.g. neurons, lymphocytes and macrophages) which other retrovirus vectors could not be used for. The lentivirus vector is increasingly used to integrate siRNA efficiently in a wide variety of cell lines and primary cells, both in vitro and in vivo.

Rapid virus purification by Membrane Chromatography

The Sartobind® ion exchange membrane adsorber technology used in Adenopack and Lentiselect is unique in its capability to efficiently and rapidly capture and recover large virus particles. When compared to chromatography media, membrane adsorbers provide large 3000 nm pores allowing unrestricted access and recovery of virus from the charged adsorber surface. Convective flow through the syringe filter devices provides high-speed separations not possible with traditional chromatography, cesium chloride density gradients and ultracentrifugation methods.

Our membrane adsorbers with porous matrices, high capacities, low differential pressures, high flow rates and low unspecific adsorption show an excellent performance in small scale virus purification. Additionally, they are also scalable and confirm to cGMP facilities to large volume, high performance separation, reducing the processing time by a factor of 10 in the final process.

Adenovirus Purification with Vivapure® Adenopack Kits

Adenopack 20 | 100 | 500

The Adenopack adenovirus purification and concentration kits offer researchers who need to recover up to 3×10^{13} purified recombinant adenovirus particles for invitro transfection a fast, safe and easy to use solution. The kits include all reagents and devices necessary for clarification, purification and concentration of adenovirus type 5 from HEK293 cell cultures in only two hours. These straight forward kits replace time-consuming and labor-intensive 48 hour CsCl density gradients.

Adenopack kits are offered as Adenopack 20, Adenopack 100 and Adenopack 500, for the purification and concentration of adenovirus type 5 from 20 mL to 500 mL cell culture, leading to 1×10^{11} - 3×10^{13} purified viral particles. For each sample volume, the most convenient handling method is offered for ultimate convenience.

To this end, preparations using Adenopack 20 are pursued in spin column format in a centrifuge, Adenopack 100 is a manually operated kit in syringe filter format*, and Adenopack 500 is a pump driven kit.

Adenopack advantages

Fast and easy virus purification

- Purification completed in 2 hours
- Convenient, over 10 × faster alternative to CsCl density gradient

Quantitative yields

- In contrast to CsCl density gradient, the complete cell culture is used for virus purification and not only the viral pellet

Flexible product range

- Applicable from initial construct screening to large scale virus production

Complete Kit

- Including filtration devices, Adenopack units for virus purification, Vivaspin® and all buffers

Low endotoxin levels

- High cell viability and infection rates due to endotoxin levels of <0.025 EU/mL

Purification results from preparations with Ad5 GFP-constructs

Purification method	Process time	Eluate	Recovery***	Viral Particles
Adenopack 20 20 mL culture	1 hour	1 mL	65 – 70%	$1 \times 10^{11-12}$
Adenopack 100 60 mL culture	1–2 hours	1 mL	65%	$1–3 \times 10^{12}$
Adenopack 100 200 mL culture	2 hours	1 mL	80%	1×10^{13}
Adenopack 500 500 mL culture	2 hours	1 mL	80%	$1–3 \times 10^{13}$
500 mL CsCl density gradient	24–48 hours	1–2 mL**	60–70%	$1 \times 10^{11-12}$

* Vivapure® Adenopack 100 can optionally be operated with a laboratory pump and an infusion pump, for which protocols are provided on our web page www.sartorius.com. Additionally, the tubes and adaptors needed for these operation modes can be ordered.

** after dialysis

*** before buffer exchange

Vivapure® Adenopack 20

The optimal kit for construct screening



Vivapure® Adenopack 20 is the downscale kit in the Adenopack series, purifying up to 1×10^{12} adenovirus type 5 particles from 20 mL cell culture. Especially when testing new constructs, parallel and fast purifications of different adenoviruses are essential. This kit allows the rapid, simple and affordable spin column based purification of 6 different samples in

parallel and bridges a gap in the CsCl density gradient method – for the first time adenovirus type 5 can efficiently be purified from less than 100 mL cell culture volume!

Typical performance

For a normal yielding vector, 1×15 cm culture plate purified using this method yields up to 1×10^{12} viral particles.


Vivapure® Adenopack 20 contents and ordering information

Vivapure® Adenopack 20	VS-AVPQ020
Vivapure® Adenopack 20 RT*	VS-AVPQ022
Adenopack Maxi spin columns	6
Vivaclear Maxi 0.45 µm PES	6
Empty 50 mL tubes	6
Loading Buffer (10×)	25 mL
Washing Buffer (10×)	30 mL
Elution Buffer	20 mL
Benzonase® (12.5 U/µL)	120 µL
Vivaspin® 20, 100 kDa MWCO	6
Instructions	1 each for Kit and Vivaspin®

Technical data

Kit specifications

Sample size	20 mL of cell culture
Number of purifications	6×20 mL
Virus particles (VP) per mL	Typically up to $1 \times 10^{11} - 10^{12}$
VP IU	50-100
Processing time	Typically 1 hour
Endotoxin level	< 0.025 EU/mL



* Adenopack 20 RT does not include Benzonase®

** Benzonase® Nuclease is manufactured by Merck KGaA, Darmstadt, Germany and is covered by US Patent 5,173,418 and EP Patent 0,229,866. Nycomed Pharma A/S (Denmark) claims worldwide patent rights to Benzonase® Nuclease, which are licensed exclusively to Merck KGaA, Darmstadt, Germany. Benzonase® is a registered trademark of Merck KGaA, Darmstadt, Germany.

Vivapure® Adenopack 100

Fast purification of up to 1×10^{13} viral particles



Vivapure® Adenopack 100 is optimally suited for adenovirus purification from up to 200 mL cell culture for in vitro transfection. This flexible kit contains two Adenopack 100 units, which can be either used in tandem for the purification of up to 200 mL cell culture for recovering 1×10^{13} viral particles or individually for purifying $1-3 \times 10^{12}$ viral particles from up to 60 mL cell culture. The purification is pursued manually with a syringe optimally attached to a retort stand.

However, for even more convenience, protocols are provided for optionally running the virus purification with a peristaltic pump or with an infusion pump, in addition to detailed instructions for a manual operation supplied with the kit. The accessories needed for the operation with a pump are supplied as individual products.

Typical performance

For a normal yielding vector, 10×15 cm culture plate purified using this method yields up to 1×10^{13} viral particles.

Vivapure® Adenopack 100 contents and ordering information

Vivapure® Adenopack 100	VS-AVPQ101
Vivapure® Adenopack 100 RT*	VS-AVPQ102
Adenopack 100 units	2
Minisart® NML Plus	4
20 mL syringe	4
Tubing set and one way valve	2
10 mL syringe (elution)	2
Loading Buffer (10×)	1 × 25 mL
Washing Buffer	1 × 120 mL
Elution Buffer	1 × 20 mL
Benzonase® 12.5 U/μL	200 μL
Vivaspin® 20 concentrator	4
Instructions	1 each for Kit and Vivaspin®

Adenopack 100 accessories


Pump tubing set for Vivapure® Adenopack 100	VS-AVPA001
---	------------

* Adenopack 100 RT does not include Benzonase®

Technical data

Kit specifications

Sample size	20–200 mL of cell culture
Number of purifications	2 × 20–60 mL 1 × 200 mL
Virus particles (VP) per mL	Typically up to 1×10^{13}
VP IU	20–50
Processing time	Typically 2 hours
Endotoxin level	< 0.025 EU/mL



Vivapure® Adenopack 500

Pump driven kit for larger volumes



Vivapure® Adenopack 500 is the direct upscale kit to the Adenopack 100, for adenovirus purification. In only 2 hours up to 3×10^{13} adenovirus particles are purified and concentrated from 500 mL cell culture. This completely ready-to-use kit is conveniently operated by a laboratory pump, offering optimal flow control and minimal

hands-on time. This easy to use product replaces lengthy and inefficient cesium chloride density gradient methods.

Typical performance

For a normal yielding vector, 25×15 cm culture plate purified using this method yields up to 3×10^{13} viral particles.


Vivapure® Adenopack 500 contents and ordering information

Vivapure® Adenopack 500	VS-AVPQ501
Vivapure® Adenopack 500 RT*	VS-AVPQ502
Adenopack 500 unit	1
Sartopore® 2 150	1
Tubing set and one way valve	2
10 mL syringe	1
Loading Buffer (10×)	60 mL
Washing Buffer (10×)	30 mL
Elution Buffer	20 mL
Benzonase® 12.5 U/μL	500 μL
Vivaspin® 20 concentrator	2
Instructions	1 each for Kit and Vivaspin®

Technical data

Kit specifications

Sample size	500 mL of cell culture
Number of purifications	1 × 500 mL
Virus particles (VP) per mL	Typically up to 3×10^{13}
VP IU	20-50
Processing time	Typically 2 hours
Endotoxin level	< 0.025 EU/mL



* Adenopack 500 RT does not include Benzonase®

Lentivirus Purification with Vivapure® Lentiselect Kit

Lentiselect 40 | 500 | 1000

The Lentiselect lentivirus purification and concentration kits offer researchers who need to recover up to 5×10^9 infective lentivirus particles per mL for invitro transfection or animal studies a fast and easy to use solution.

These straight forward kits replace time-consuming ultracentrifugation protocols, which typically take approximately one day for large sample volumes, thus reducing the purification time to only a few hours.

Lentiselect kits are offered as Lentiselect 40, Lentiselect 500 and Lentiselect 1000 for the purification and concentration of VSV-G pseudotyped lentivirus from 40 mL to 1 L cell culture, leading to $8 \times 10^8 - 1 \times 10^{10}$ purified infective particles. For each sample volume, the most convenient handling method is offered. To this end, 40 mL sample volumes are processed manually with Lentiselect 40, while Lentiselect 500 and 1000 are pump driven kits.

Lentiselect advantages

Fast and easy virus purification

- Purification completed in under one to six hours, depending on sample volume
- Kit as easy to use as filtration

No need for expensive instruments

- Lentivirus purification with Lentiselect is independent of equipment such as ultracentrifuges

High virus purity

- Achieve pure virus due to a chromatography purification for your experiments instead of a crude and variable cell culture supernatant pellet

Optimal for multiple virus construct screening

With Lentiselect 40, four purification runs can be conducted in parallel with one kit

Complete Kits

- Including Lentiselect units for virus purification, Vivaspins for concentration | buffer exchange and all buffers and syringes necessary

Low endotoxin levels

- High cell viability and infection rates due to endotoxin levels of <0.025 EU/mL

Purification results from preparations with VSV-G pseudotyped lentivirus constructs

Purification method	Process time	Eluate	Viral Particles/mL	Recovery	Infective Viral Particles
Lentiselect 40 40 mL sample	45 min	200 μ L*	4×10^9	50%	8×10^8
Lentiselect 500 500 mL sample	3 hours	1 mL*	3×10^9	35%	$2 - 5 \times 10^9$
Lentiselect 1000 1 L sample	6 hours	2 mL*	5×10^9	35%	1×10^{10}
Ultracentrifugation 500 mL sample	10 - 11 hours	500 μ L	6×10^9	25%	3×10^9

* After desalting | buffer exchange

Vivapure[®] Lentiselect 40

Fast purification of up to 8×10^8 viral particles



Vivapure[®] Lentiselect 40 is optimally suited for lentivirus purification for up to 40 mL cell culture and contains all components necessary for 4 purifications. Up to 8×10^8 viral particles are recovered in less than one hour. In contrast to traditional ultracentrifugation methods, virus purification with Vivapure[®] Lentiselect is fast and simple, without the need for expensive equipment like an ultra-

centrifuge. Additionally, this chromatographic procedure leads to pure virus samples in contrast to the crude ultracentrifuge pellet, resulting in higher reproducibility and increased gene transfer efficiency.

Typical performance

For a normal yielding vector, 2×15 cm culture plate purified using this method yield up to 8×10^8 particles.


Vivapure[®] Lentiselect 40 contents and ordering information

Vivapure [®] Lentiselect 40	VS-LVPQ040
Lentiselect unit	4
50 mL syringe	4
10 mL syringe	4
Tube set with one-way valve	4
Loading buffer (10 ×)	30 mL
Washing buffer	150 mL
Elution buffer	20 mL
Vivaspin [®] 20, 100 kDa MWCO	8
Instructions	1 each for Kit and Vivaspin [®]

Technical data

Kit specifications

Sample size	40 mL cell culture
Number of purifications	4 × 40 mL
Infective particles (P) per mL	Typically up to 3×10^9
VP IU	5–15
Processing time	Typically 45 minutes
Endotoxin level	< 0.025 EU/mL



Vivapure® Lentiselect 500

Fast purification of up to $2-5 \times 10^9$ infective particles per mL from 500 mL cell culture



Vivapure® Lentiselect 500 is optimally suited for VSV-G pseudotyped lentivirus purification from up to 500 mL cell culture and contains all reagents and devices necessary for purifying up to $2-5 \times 10^9$ infective particles.

The whole purification procedure is simply operated by a laboratory pump, which minimizes hands-on time. Unlike conventional purification methods as ultracentrifugation,

Vivapure® Lentiselect 500 offers a fast and simple solution for purifying VSV-G pseudotyped lentiviruses making expensive purification equipment like ultracentrifuges redundant.

Typical performance

For a normal yielding vector, 500 mL cell culture purified using this method yield up to $2-5 \times 10^9$ infective particles per millilitre (total volume 1 mL).


Vivapure® Lentiselect 500 contents and ordering information

Vivapure® Lentiselect 500	VS-LVPQ500
Lentiselect unit	1
Sartopore® 2 150	1
50 mL syringe	1
Tube set with one-way valve	1
Loading buffer (10 ×)	30 mL
Washing buffer	170 mL
Elution buffer	30 mL
Vivaspin® 20, 100 kDa MWCO	2
Operating manual	1 each for Kit and Vivaspin®

Technical data

Kit Specifications

Sample volume	500 mL cell culture
Number of purifications	1 × 500 mL
Infective particles (IP) per mL	Typically up to $2-5 \times 10^9$ *
Processing time	Typically up to 3 hours
Endotoxin level	< 0.025 EU/mL



*1 mL final elution sample

Vivapure® Lentiselect 1000

Pump driven kit for larger sample volumes



Vivapure® Lentiselect 1000 is the direct scale up kit to Lentiselect 500, for VSV-G pseudotyped lentivirus purification. The rapid 6 hour protocol results in a recovery of $4-5 \times 10^9$ infective particles per mL (total volume 2 mL) from 1 L cell culture supernatant.

This kit is to be operated by a laboratory pump and contains all necessary buffers and ultrafiltration devices for

optimal convenience. The traditional time consuming ultracentrifugation method is replaced by this fast and simple Vivapure® Lentiselect 1000 kit.

Typical performance

For a normal yielding vector, 1 L cell culture purified using this method yield up to $4-5 \times 10^9$ infective particles per millilitre (total volume 2 mL).


Vivapure® Lentiselect 1000 contents and ordering information

Vivapure® Lentiselect 1000	VS-LVPQ1000
Lentiselect unit	2
Sartopore® 2 150	1
50 mL syringe	1
Tube set with one-way valve	1
Loading buffer (10 ×)	30 mL
Washing buffer	170 mL
Elution buffer	60 mL
Vivaspin® 20, 100 kDa MWCO	2
Operating manual	1 each for Kit and Vivaspin®

Technical data

Kit specifications

Sample volume	1 L cell culture
Number of purifications	1 × 1 L
Infective particles (IP) per mL	Typically up to $4-5 \times 10^9$ *
Processing time	Typically up to 6 hours
Endotoxin level	< 0.025 EU/mL



* 2 mL final elution sample



Application Notes

Table of Contents

1. Desalting and Buffer Exchange with Vivaspin® Centrifugal Concentrators	74
2. Treatment of Vivaspin® Concentrators for Improved Recovery of Low-concentrated Protein Samples	77
3. Scouting Protein Purification Conditions Using Vivapure® Centrifugal Ion Exchange Membrane Absorbers	80
4. Concentration and Purification of Viruses by using Ultrafiltration, Incl. Coronavirus – a Short Review	85
5. Sartorius Ultrafiltration Products in the Preparation of Biological Nanoparticles and Medical Nanocarriers	97
6. Vivaflow® and Vivaspin® Workflow in Protein Research Laboratories	103

1. Desalting and Buffer Exchange with Vivaspin® Centrifugal Concentrators

Introduction

Vivaspin® centrifugal concentrators, with patented vertical membrane technology, combine fast filtration with high recovery of target proteins. This makes Vivaspin® the technology of choice for desalting or buffer exchange, avoiding lengthy dialysis steps.

While proteins are retained by an ultra-filtration membrane, salts can pass freely through, independent of protein concentration or membrane MWCO. In consequence, the composition of the buffer in the flow-through and retentate is unchanged after protein concentration. By diluting the concentrate back to the original volume, the salt concentration is lowered. The concentrate can be diluted with water or salt-free buffer if simple desalting is required; however, it is also possible to dilute the concentrate with a new buffer, thereby exchanging the buffering substance entirely. For example, a 10 mL protein sample containing 500 mM salt, if concentrated 100-fold still contains 500 mM salt. If this concentrate is then diluted 100× with water or salt-free buffer, the protein concentration returns to the original level, while the salt concentration is reduced 100× to only 5 mM (i.e. a 99% reduction in salt concentration).

The protein sample can then be concentrated again to the desired level, or the buffer exchange can be repeated to reduce the salt concentration even further before a final concentration of the protein. This process is called diafiltration. For proteins with a tendency to precipitate at higher concentrations, it is possible to perform several diafiltration steps in sequence, with the protein concentrated each time to only 5 or 10×. For example, if a precipitous protein sample is concentrated to 5× then diluted back to the original volume, and this process is repeated a further two times, this still results in a >99% reduction in salt concentration, without over-concentrating the protein.

Methods

Select an appropriate MWCO for your sample. For maximum recovery, select a MWCO $\frac{1}{3}$ to $\frac{1}{2}$ the molecular weight of the molecule of interest.

1. Add protein sample up to the maximum fill volume of the concentrator (as stated in the device operating instructions). If the sample volume is lower than the maximum device volume, it can be diluted to the maximum fill volume before the first centrifugation step. This will increase the salt removal rate.
2. Centrifuge for the recommended amount of time at an appropriate spin speed (see device operating instructions).
3. Empty filtrate container and refill the concentrator with an appropriate exchange solvent.*
4. Centrifuge again as before.
5. Recover the concentrated, desalted sample from the bottom of the concentrate pocket with a pipette.




Figure 1: Step-by-step method for desalting and concentration

* Filtrate volumes should be retained until the concentrated sample has been analyzed.

Results

Vivaspin® 20

MWCO	5 kDa		30 kDa		50 kDa		100 kDa	
	Cytochrome C 0.25 mg/mL		BSA 1 mg/mL		BSA 1 mg/mL		IgG 1 mg/mL	
	Protein Recovery	NaCL Removal	Protein Recovery	NaCL Removal	Protein Recovery	NaCL Removal	Protein Recovery	NaCL Removal
Spin 1	100%	99%	97%	99%	97%	99%	90%	98%
Spin 2	96%	100%	92%	100%	93%	100%	87%	100%

Four Vivaspin® 20 devices of each MWCO were tested with 20 mL samples. Each sample contained 500 mM NaCl. To perform diafiltration, devices were centrifuged at 4,000 g for 45 min (5 kDa MWCO) or 30 min (>5 kDa MWCOs). After the first and second spins, the retentate samples were

brought up to 20 mL with ultrapure water from an Arium® system (Sartorius). OD readings were taken at 410 nm for Cytochrome C or 280 nm for BSA and IgG samples. Salt concentrations were measured using a Qcond 2200 conductivity measuring instrument.

Vivaspin® 6

MWCO	5 kDa		30 kDa		50 kDa		100 kDa	
	Cytochrome C 0.25 mg/mL		BSA 1 mg/mL		BSA 1 mg/mL		IgG 1 mg/mL	
	Protein Recovery	NaCL Removal	Protein Recovery	NaCL Removal	Protein Recovery	NaCL Removal	Protein Recovery	NaCL Removal
Spin 1	98%	99%	92%	99%	93%	99%	92%	98%
Spin 2	85%	100%	86%	100%	83%	100%	89%	100%

Four Vivaspin® 6 devices of each MWCO were tested with 6 mL samples. Each sample contained 500 mM NaCl. To perform diafiltration, devices were centrifuged at 4,000 g for 45 min (5 kDa MWCO) or 30 min (>5 kDa MWCOs). After the first and second spins, the retentate samples were

brought up to 6 mL with ultrapure water from an Arium® system (Sartorius). OD readings were taken at 410 nm for Cytochrome C or 280 nm for BSA and IgG samples. Salt concentrations were measured using a Qcond 2200 conductivity measuring instrument.

Conclusions

As the results show, the efficient design of Vivaspin® devices allowed >95% of the salt to be removed during the first centrifugation step. Only one subsequent centrifugation step was needed to increase the typical salt

removal to 99% with >92% recovery of the target protein. Diafiltration using ultrafiltration devices such as Vivaspin® 6 and 20 represents a faster and more efficient solution to desalting and buffer exchange, than conventional techniques such as dialysis.

2. Treatment of Vivaspin® Concentrators for Improved Recovery of Low-Concentrated Protein Samples



Introduction

With appropriate device size and membrane cut-off selected, Vivaspin® products will typically yield recoveries for the concentrated sample > 90% when the starting sample contains over 0.1 mg/mL protein of interest. Depending on sample characteristics relative to the membrane type used, solute (protein) adsorption on the membrane surface is typically very low (2 – 10 µg/cm²) and in practice not detectable.

This can increase to 20 – 100 µg/cm² when the filtrate is of interest and the sample must pass through the whole internal structure of the membrane. Whilst the relative adsorption to the plastic of the sample container will be proportionately less important than on the membrane, due to the higher total surface area, this can be also be a source of yield loss. Typically, a higher cut-off membrane will bind more than a low molecular weight alternative.

Whenever possible, the smallest MWCO and device size applicable should be chosen. Swinging bucket rotors are preferred to fixed angle rotors. This reduces the surface area of the concentrator that will be exposed to the solution during centrifugation.

An important factor not to be neglected is the thorough recovery of the retentate. Make sure to carefully remove all traces of solution from the sample container and, if feasible, rinse the device after recovering the sample with one or more drops of buffer and then recover again.

The intention of the following “passivation” procedure is to improve recovery of protein samples in the nano- to microgram concentration range by pretreating the device (membrane & plastic). For this purpose a range of solutions are suggested in Table 1.

Table 1: Passivation Solutions

Type	Concentration
Powdered milk	1% in Arium® water
BSA	1% in PBS
Tween 20	5% in Arium® water
SDS	5% in Arium® water
Triton X-100	5% in Arium® water
PEG 3000	5% in Arium® water

Passivation procedure for Vivaspin® ultrafiltration concentrators

A) Passivation procedure

1. Wash the concentrators once by filling with Arium® water and spin the liquid through according to the respective protocol.
2. Remove residual water thoroughly by pipetting.
Caution: Take care not to damage the membrane with the pipette tip.
3. Fill concentrators with the blocking solution of choice as given in Table 1.
4. Incubate the filled concentrators at room temperature for at least 2 hours (overnight is also possible except for Triton X-100 which is not recommended for overnight incubation).
5. Pour out the blocking solution.
6. Rinse the device 3 – 4 × very thoroughly with Arium® water and finally spin through.
7. The “passivated” devices are now ready for use. We recommend comparing different passivation reagents with an untreated device.

Note

It is necessary to rinse the device thoroughly before each washspin to ensure that traces of passivation compound are removed from the deadstop. Use the device immediately for protein concentration or store it at 4°C filled with Arium® water, to prevent the membrane from drying.

B) Evaluation of passivation effects (exemplary with BSA)

1. Prepare a 10 µg/mL BSA stock solution e.g. by diluting 90 µL of the 4 mg/mL stock solution in 36 ml 0.1 M sodium borate pH 9.3. Mix well.
2. Fill Vivaspin® 2 devices with 2 ml of this 10 µg/mL BSA solution and close with cap provided.
3. Spin the device in a swing-out rotor at 4,000 × g until the volume is to app. 100 µL.
4. Recover the concentrate and make back up to 2 mL with 0.1 M sodium borate pH 9.3
5. Determine recovered protein concentrations e.g. according to Bradford or BCA assays.

Results and discussion

As an example, the effect of milk powder was analysed. It could be shown (Figure 1) that the protein recovery of a 10 µg/mL BSA solution could be increased from around 70 to 90%. If milk powder is not interfering with sample purity and quality, it is a good starting point to improve recovery of diluted sample solutions.

Protein recovery (10 µg/mL BSA) with Vivaspin® PES 10 kDa after passivation

In another example, detergents were analysed with only 250 and 500 ng BSA (Figure 2) BSA recovery declined to 50 – 30% in untreated devices as the protein concentration was reduced. Significant improvement to 60 – 90% recovery could be demonstrated when using the passivation strategy. Often, Triton X-100 seemed to work though the optimal reagent has to be selected for the respective protein and its hydrophilic | -phobic characteristics.




Figure 1: Protein recovery (10 µg/mL BSA) with Vivaspin® PES 10 kDa after passivation




Figure 2: Protein recovery (250 and 500 ng BSA) with Vivaspin® 2 PES 10 kDa after passivation

Summary

Passivation is an appropriate method to achieve increasing sample recovery when using very dilute samples. In addition to skimmed milk, other proteins (BSA), detergents and compounds are possible. However, it should be noted that this is a general procedure, not specific for any particular application. Depending on the hydrophilic | -phobic character of the protein non-specific binding may be more or less of a problem and the suggested passivation solutions may lead to different results. Even with the Hydrosart membrane, which is recommended for dilute samples, passivation of the

device will reduce losses on the plastic surface. One very important thing to remember is that the blocking agent is potentially introduced into the sample. It should be assured that this will not interfere with downstream analysis.

For example, proteins must not be used for passivation if a pure protein is intended to be concentrated for x-ray crystallography, as even the smallest traces would interfere with the diffraction pattern. Other subsequent analyses methods include activity testing, gel electrophoresis or labelling are less problematic.

3. Scouting Protein Purification Conditions Using Vivapure® Centrifugal Ion Exchange Membrane Absorbers

Introduction

For separation and purification of proteins from biological samples, different characteristics of the target protein e.g. its size, charge, hydrophobicity or specifically engineered tags, are exploited.

With ion exchange chromatography, separation is achieved on the basis of charge differences between biomolecules. This makes it a versatile method often used for pre-fractionation or purification of a target protein from crude protein mixtures. To optimize the purification procedure for an individual target, several binding and elution conditions have to be tested on cation and anion exchange matrices.

In contrast to traditional column chromatography methods, Vivapure® IEX centrifugal columns allow scouting of several chromatography conditions in parallel, leading quickly to different fractions which can be further analyzed for enriched or even already purified target protein.

Here, we demonstrate the performance of Vivapure® IEX Mini spin columns for evaluation of optimal purification conditions of cloned SH2 domains from an *E. coli* lysate in a two step procedure. This protocol can generally be employed for identifying a purification method based on ion exchange chromatography for a given target protein, as it is fast and only uses small amounts of the sample.

In the first step of this study, binding conditions were evaluated by loading the sample on Vivapure® Q and S columns at various pH values, eluting bound proteins with a high salt concentration buffer and analyzing all fractions for the target protein.

The results from this experiment provided the optimal binding pH and the best ion exchange chemistry for the purification of SH2 domain.

In a second step, the best elution method was evaluated by applying increasing salt concentrations to columns which were shown to bind the target protein in step one, leading to a complete purification protocol in less than one hour.

Materials and Methods

Buffers tested to determine the optimum pH and salt concentration for binding and elution in ion exchange purification of SH2 domain

Buffer A:	25 mM Citrate, pH 4
Buffer B:	25 mM Potassium phosphate, pH 6
Buffer C:	25 mM HEPES, pH 8
Buffer D:	25 mM Sodium bicarbonate, pH 10
Buffer E:	25 mM Citrate, pH 4, supplemented with 1 M NaCl.
Buffer F:	25 mM Potassium phosphate, pH 6, supplemented with 0.2 M, 0.4 mM, 0.6 mM, 0.8 mM and 1 M NaCl, respectively.
Buffer G:	25 mM HEPES, pH 8, supplemented with 1 M NaCl
Buffer H:	25 mM Sodium bicarbonate, pH 10, supplemented with 1 M NaCl

Scouting Binding Conditions

300 mL LB media was inoculated with 4 mL of an overnight culture and incubated at 37°C, shaking at 150 rpm until an OD₆₀₀ of 1.0 was reached. IPTG was added to a final concentration of 1 mM and the culture incubated for a further 4 h with shaking at 150 rpm. Cells were harvested by centrifugation at 4,000 g for 30 min at 4°C. The pellet was resuspended in 35 mL PBS (150 mM KPi, pH 7.3) and cells were lysed by addition of lysozyme to a final concentration of 0.1 mg/mL and incubation for 1 h at 37°C. Insoluble particles and cell debris were removed by centrifugation at 10,000 g for 30 min at 4°C.

4 × 200 µL aliquots of the cell lysate were diluted with 1.8 mL binding buffer A to D, to adjust each sample to the respective pH being tested. To avoid clogging of the membranes in the Vivapure® Mini spin columns, samples were clarified by passage through 0.45 µm CA Minisart NML syringe filters (Sartorius).

4 × Q and 4 × S Vivapure® Mini spin columns were labeled 4, 6, 8 and 10, corresponding to the pH of the buffer to be used. To each spin column, 400 µL of the corresponding binding buffer was added and spun for 5 minutes at 2,000 g (45° fixed angle rotor).

400 µL of the clarified samples adjusted to pH 4, 6, 8 or 10 were applied to each of the correspondingly equilibrated Vivapure® Q and S spin columns. Columns were spun for 5 min at 2,000 g.

Afterwards, Vivapure® Mini spin columns were reloaded with 400 µL sample and spun again for 5 min at 2,000 g. Loosely bound proteins were washed away with the application of 400 µL of the respective binding buffer to each of the columns and spinning for 5 min at 2,000 g. Flow-through and wash fractions were collected for subsequent detection of the target protein.

200 µL of elution buffer E, F, G or H, were applied to the washed columns and spun for 3 min at 2,000 g. Eluates were saved for subsequent analysis.

4 µL of flow-through, wash, and eluate fractions from each column were analyzed by reducing SDS-PAGE followed by silver staining.

Optimizing Elution Conditions

Taking account of the results of the first experiment (Scouting Binding Conditions) 200 µL cell lysate was diluted with 1.8 mL binding buffer B (25 mM KPi, pH 6). To avoid clogging of the membrane in the Vivapure® Mini spin column, the pH adjusted sample was clarified by passage through a 0.45 µm CA Minisart NML syringe filter (Sartorius).

400 µL binding buffer B was applied to one Vivapure® S Mini spin column and spun for 5 minutes at 2,000 g.

400 µL of the clarified sample was applied to the equilibrated Vivapure® S column and spun for 5 min at 2,000 g. Afterwards, the Vivapure® S Mini spin column was reloaded with 400 µL sample and spun again for 5 min at 2,000 g.

Loosely bound proteins were washed away by application of 400 μ L binding buffer to the column and spinning for 5 min at 2,000 g. Flow-through and wash fractions were saved for analysis.

To elute the target protein, 100 μ L elution buffer F, supplemented with 0.2 M NaCl was applied to the Vivapure[®] S Mini spin column and spun for 3 min at 2,000 g. The eluate was collected. For the next elution step, 100 μ L of elution buffer F, supplemented with 0.4 M NaCl was applied and again spun for 3 min at 2,000 g. Elution was continued with 0.2 M NaCl increments until a final salt concentration of 1 M was reached, saving the eluates from each step.

4 μ L of flow-through, wash, and eluate fractions from each column were analyzed by reducing SDS-PAGE followed by silver staining.


Results

Dilution of the *E. coli* lysate with binding buffer A (25 mM Citrate, pH 4) lead to complete precipitation of sample proteins. Thus, pH 4 could not be tested in this experiment. As can be seen on the SDS gel in (Figure 1), the target protein was present in the eluates from Vivapure[®] Q Mini spin columns at all pH values tested, together with most of the *E. coli* proteins (Lanes Q "e"). In contrast, using the Vivapure[®] S Mini spin column, at all pH values tested, most *E. coli* proteins did not bind to the membrane and were found in the flow-through (Lane S "f"), thus resulting in purer target protein in all eluate fractions (Lane S "e").

Differences could be detected in the binding efficiency of the target protein. At pH 8, traces of the target protein were already found in the flow-through, with slightly higher amounts at pH 10 (Lane S "e"). At pH 6, the most efficient binding of the target protein to the S membrane was observed.

The purification conditions determined for Vivapure[®] S with potassium phosphate buffer (pH 6) were further optimized to determine the ideal salt concentration for SH2 domain elution. The target protein started to elute with 200 mM NaCl, however the main fraction eluted with 400 mM NaCl. Traces of the target protein were also found in the next elution step with 600 mM NaCl, but this might be due to the low elution volume.


Figure 1



Sample	Sample Volume (μL)	Volume Loaded Onto Gel (μL)
M = Broad range marker		
s = Sample before application	800	4
f = Flow-through	800	4
w = Wash fraction	400	4
e = Elution with 1 M NaCl	200	4

Note: Scouting for optimal binding conditions of a SH2 domain expressed in *E. coli*. 12% reducing SDS gel, silver stained, shows the sample before purification (s), flow-through (f), wash (w) and eluate (e) fractions (1 M NaCl) from Vivapure® Q and S Mini spin columns, at the various pH values tested.

Figure 2



Sample	Process Volume (μL)	Volume Loaded Onto Gel (μL)
M = Broad range marker		
s = Sample before application	800	16
f = Flow-through	800	16
w = Wash fraction	400	16
e1 = 25 mM KPi, pH 6, 200 mM NaCl	100	8
e2 = 25 mM KPi, pH 6, 400 mM NaCl	100	8
e3 = 25 mM KPi, pH 6, 600 mM NaCl	100	8
e4 = 25 mM KPi, pH 6, 800 mM NaCl	100	8
e5 = 25 mM KPi, pH 6, 1 M NaCl	100	8

Note: Optimizing elution conditions for a SH2 domain expressed in *E. coli*, using Vivapure® S Mini spin column at pH 6. 12% reducing SDS gel, silver stained, shows the sample before purification (s), flow through (f), wash (w) and eluate (e1-5) fractions.

Conclusion

A two-step procedure was used to rapidly scout optimal purification conditions for a target protein (a SH2 domain from *E. coli* lysate) with ion exchange chromatography. In the first step, the most suitable ion exchanger and buffer pH for binding the target protein was verified. In the second step, the elution condition was optimized, building on the results gained in step one. With the scouting procedure described here, it was possible to quickly and conveniently purify the target protein to homogeneity.

The results obtained in this experiment can be used to various ends, e.g:

- Polishing a specific protein after purification with another chromatographic technique
- Quickly establishing a FPLC method for a new protein
- Identification of the optimal purification method prior to scale up with Vivapure® IEX Maxi spin columns.

For these purposes Vivapure® IEX Mini and Maxi spin columns and Sartobind® membrane adsorber units with FPLC connectors are available.

4. Concentration and Purification of Viruses by Using Ultrafiltration, Incl. Coronavirus – A Short Review

Introduction

Evolutionary, viruses developed various mechanisms to interact and manipulate the genetic material of their target cells. Based on this, modern molecular biology utilizes viruses in a constantly growing number of applications.¹ They range from controlled genetic transfection of cells to a variety of different basic studies in medical science.² In medical studies the strategic focus is on recombinant vaccines and on the development of potential vectors for gene therapy.^{3,4}

Besides the great relevance of viruses for medical applications, the assessment of virus type and content is important for the risk assessment of food and drinking water.⁵ Also, the classification of virus content is often of high relevance for the quality control of aquatic biotopes.⁶

During the preparation, handling, or analysis of viruses or virus-like particles (VLPs), a concentration and/or purification step is frequently required.⁵ Typical viruses have a size within the range of about 20 nm up to several hundred nanometers.⁷ Therefore they are ideally suited for the retention on ultrafiltration membrane systems and such ultrafilters are widely used in basic virus research. The specifications of such ultrafiltration devices depend on the particular type of virus and the purpose of the subsequent application.

This short review highlights methods for the purification of various mammalian viruses for basic medical research. Also, the concentration of pathogenic viruses from water and food samples and the purification of marine bacteriophages (virioplankton) are highlighted. It will also give guidance for the selection of an ideal performing device with the optimum molecular weight cut-off (MWCO) for the user specified ultrafiltration process.

Concentration of mammalian viruses in medical research

In medical research viruses and VLPs are of major interest, particularly for investigations on infectious viral diseases and for the development of vaccines or antiviral drugs. Moreover, certain VLPs can manipulate genetic material in a directed manner and are used broadly in the development of genetic therapy approaches. Additionally, viral vectors are well established as a transfection method for gene transfer to cell lines e. g. to manipulate mammalian cells in vivo and in vitro.

An overview of thematically linked publications using Sartorius ultrafiltration devices for the purification and concentration of viruses and VLPs in the medical context is given in Table 1. Among other applications, Vivaspin® devices were employed for the concentration of adeno-associated virus (AAV) and lentiviral vectors after purification via ion exchange chromatography,⁸⁻¹⁰ on blood sera to prepare blank samples from hepatitis C virus (HCV)-positive blood sera,¹¹ for the development of a vaccine against human immunodeficiency virus (HIV) and of an antiviral drug against Chikungunya virus.^{12,13}

Table 1: Summarized examples of applications with Vivaspin® and Vivaflow® for of viruses in medical research

Goal of Research (Type of virus, Host Organism)	Purpose of Filtration (Buffer System)	Sartorius Ultrafiltration Device (MWCO)	Subsequent Step	Ref.
Gene therapy (Adenovirus type 5, VLP, human)	Diafiltration (20 mM Tris saline buffer)	Vivaflow® (100 kDa)	Storage, chromatography on Sartobind® STIC membrane absorber (FPLC)	14
Reduction of HCV-induced fibrosis (Hepatitis C Virus; human)	Removal of HCV from human blood serum (Blood serum)	Vivaspin® (30 kDa)	Preparation of negative control (from positive sample) for immunofluores- cence assay, fibrosis induction assays	11
Development of a viral entry inhibitor for HIV (HIV, human)	Removal of protein fraction from virus (PBS)	Vivaspin® 20 (1,000 kDa)	Virus inactivation	12
Gene therapy for cancer treatment (adeno-associated virus; rAAV-2, human)	Concentration and purification after expression, Buffer exchange after His tag (FreeStyle 293 Expression Medium (Gibco), serum-free)	Vivaspin® 20 (1,000 kDa)	Titer, ELISA, cell binding assay, apoptosis cell cycle assay	8
System for controlled gene expression in mice brain (Adeno- associated virus, mice)	Concentration of eluate after anion exchange chromatog- raphy (elution buffer)	Vivaspin® 20 (100 kDa)	Transduction of mice neurons	9
Efficient gene transfer into the CNS (Lentivirus, human)	Concentration after ion exchange chromatography (PBS)	Vivaspin® (100 kDa)	Quantification via real- time PCR and end-point dilution. Transduction of murine neuronal and glial cells in vivo	10
Identification of effective chikungunya antiviral drugs (Chikungunya-Virus, human)	Concentration	Vivaspin® 20 (100 kDa)	Quantification by TCID ₅₀	13
Gene therapy of achromatopsia in mice (Recombinant adeno-associated virus, human virus used in mice)	Concentration (Anion exchange chromatography elution buffer)	Vivaspin® 4 (10 kDa)	Titer determination by dot-blot analysis, subretinal injections	15

Concentration of viruses from drinking water and food samples

The guidelines for drinking-water quality by the world health organization describe safety plans to reduce potential risks from various virus infections.¹⁶ It states that, due to the increased resistance of viruses to disinfection methods, an absence of bacterial contamination after disinfection cannot be used as a reliable indicator of the presence | absence of pathogenic viral species in drinking water supplies. Considering this, ultrafiltration can play a vital role in detecting such viral contaminations for the research on drinking water quality and food safety.

For an ultrafiltration step, the water sample does not have to be pre-conditioned and its efficacy in concentrating the virus is virtually independent

of the chemical properties and the structure of the virus.¹⁷ Thus ultrafiltration is very well suited to isolate and concentrate virus particles from water samples and is a valuable aid during the assessment of water quality. Most of the viruses which are found in water and also food samples are of fecal origin. Screening for these viruses is crucial to prevent infections. The most frequent ones are hepatitis A, hepatitis E and norovirus.¹⁸ Ultrafiltration has been described as the most appropriate method for the recovery of hepatitis A virus from vegetables and other food items.¹⁹ Detection of infectious viruses is mainly done by propagation in cell culture (plaque assay) or the detection of the viral genomes by molecular amplification techniques such as quantitative reverse transcriptase polymerase chain reaction (RT-PCR).²⁰

Table 2: Summarized examples of ultrafiltration application with Vivaspin® and Vivaflow® with viruses from drinking water and food samples

Goal of Research (Type of virus, Host Organism)	Purpose of Filtration (Buffer System)	Sartorius Ultrafiltration Device (MWCO)	Subsequent Step	Ref.
Method for the detection of norovirus genogroup I (Norovirus, human)	Concentration (PBS processed food samples)	Vivaspin® (5 kDa)	RNA extraction for real-time RT-PCR	22
Analysis of viral content in groundwater (A set of pathogenic viruses, potentially human)	Concentration of drinking water sample (Drinking water)	Vivaflow® 200 (10 kDa)	Qualitative analysis (enterovirus) by RT-nested PCR and microtiter neutralization test	21
Comparative Analysis of Viral Concentration Methods (Hepatitis A virus, human)	Concentration (0.25 M threonine, 0.3 M NaCl, pH 9.5)	Vivaspin® 20 (100 kDa)	RNA extraction for real-time RT-PCR	19
Analysis of regional outbreak of gastroenteritis due to drinking water contamination (Norovirus, Astrovirus, Rotavirus, Enterovirus, Hepatitis A virus; human)	Concentration (50 mmol/L glycine buffer, 1% beef extract)	Vivaspin® 2	Nucleic acid extraction	23

Concentration of viruses and bacteriophages from marine biological samples

In marine biology, the concentration and subsequent analysis of marine bacteriophages (virioplankton) is of major interest. They outnumber the bacterioplankton (their host organisms) by an order of magnitude and thus have an important influence on the whole marine biosphere.²⁴

As described by Wyn-Jones & Sellwood (ref. 17) ultrafiltration can be used to concentrate virus particles in water samples without any prior pretreatment of the sample and it is also practically independent from the

chemical and structural properties of the viruses. Thus, it finds wide use for the analysis of aquatic viruses. For instance, Schroeder et al. (ref. 26) were able to determine the diversity and monitor population dynamics of viruses that infect *Emiliana huxleyi*, a globally important form of photosynthetic plankton. In this study a reusable Vivaflow® 50 unit equipped with a polyethersulfone (PES) membrane with MWCO of 50 kDa was used to concentrate viruses in sea water samples prior to storage and analysis. For further examples of virus concentration from marine biological samples see table 3.

Table 3: Summarized examples of ultrafiltration applications with Sartorius Vivaflow® and Vivaspin® of samples from marine biology

Goal of Research (Type of virus, Host Organism)	Purpose of Filtration (Buffer System)	Sartorius Ultra-filtration Device (MWCO)	Subsequent Step	Ref.
Assessment of viroplankton diversity (Viroplankton, Plankton)	0.2 µm filtration for clarification, filtrate subjected to 3 kDa filter for concentration (Sea water)	Vivaflow® 200 (0.2 µm and 30 kDa)	Subsequent analysis by DNA separation on Agarose gel	25
Classification of virus (MpRNAV-01B, <i>Micromonas pusilla</i>)	Vivaflow® 200: harvest and concentration of whole cell lysate; Vivaspin®: washing (removal of CsCl)	Vivaflow® 200, Vivaspin® (30 kDa)	Classification of new virus: genome, proteins, stability, etc.	28
Assessment of genetic diversity in viroplankton (<i>Emiliana huxleyi</i> Bloom virus, Eukaryotic phytoplankton - alga)	After 0.45 µm filtration, concentration 1 L to 20 mL (Sea water)	Vivaflow® 50 (50 kDa)	PCR and Denaturing gradient gel electrophoresis	26
Investigation of gene expression during infection (<i>Emiliana huxleyi</i> virus strain 86, Eukaryotic phytoplankton - alga)	Concentration from 5 L to 20 mL (f/2 medium)	Vivaflow® 50 (50 kDa)	CsCl-gradient	27
Study on host genome integration (<i>virophage mavirus</i> , <i>Cafeteria roenbergensis</i>)	Clarification with 0.2 µm filter and concentration with 100 kDa filter (<i>Cafeteria roenbergensis</i> , f/2 medium)	Vivaflow® 200 (0.2 µm and 100 kDa)	CsCl gradients, electron microscopy	29

Concentration of Coronavirus for general research and protein research (spike protein)

Coronaviruses are spherical, enveloped, RNA based viruses that are typically 80–120 nm in diameter, but in many cases have a diameter outside of this range. Coronavirus genomes are the largest of all RNA viruses which offers a relatively large area of study. Correlatingly the potential for future mutations in this large genome may lead to future human diseases that may evolve into epidemics and pandemics, such as the previous Middle East Respiratory Syndrome (MERS-CoV), and Severe Acute Respiratory Syndrome 1 (SARS-CoV-1) and 2 (SARS-CoV-2). Hence further research into the replication, transmission, genome and structure will continue with greater investment of time and funding in the years to come.

A key component to the infection cycle is the coronavirus spike (S) protein, that mediates entry into host cells, through both attachment and membrane fusion. As such, it is a primary target for the development of novel antiviral drugs and vaccines.

The concentration and purification of both the virions and the spike proteins from cell culture and supernatants is often a key requirement to isolate the respective target, prior to structural, functional analysis and binding assays, etc.

Table 4 highlights several applications where Vivaspin® centrifugal concentrators, or Vivaflow® tangential flow filtration cassettes have been used for the concentration or Coronavirus proteins, including the spike protein. References are also provided to direct readers for detailed reading.

Table 5 provides examples of concentration of intact virions, or Coronavirus virus like particles (VLPs), with the same devices.

Table 4: Summarized examples of ultrafiltration applications with Sartorius Vivaflow® and Vivaspin® of coronavirus protein samples.

Goal of Research (Type of virus, Host Organism)	Purpose of Filtration (Buffer System)	Sartorius Ultrafiltration Device (MWCO)	Subsequent Step	Ref.
Neutralization of a SARS-CoV-2 antibody to a functionally conserved receptor binding domain (RBD) on the trimeric spike (S) protein	Buffer exchange of a SARS-CoV-2 RBD protein	Vivaspin® 20, PES (10 kDa)	Protein concentration by UV/Vis and binding affinity by Streptavidin BLI (Sartorius Octet)	30
Investigation of neutralising antibody response on a SARS-CoV-2 spike glycoprotein RBD-SpyVLP (virus-like particle) platform	Concentration of SpyTag-RBD protein construct	Vivaspin® 20, PES (10 kDa)	Purification by SEC	31
Investigation of exosome based vaccines containing coronavirus spike (S) protein, for SARS-CoV-1	Concentration of solubilized spike protein in supernatant	Vivaspin®, PES (10 kDa)	Western blot analysis	32
Analyze of the ability to redirect the functionality of the Mouse Hepatitis coronavirus spike (S) protein to infect human cancer cells	Concentration of cellular receptor protein constructs	Vivaspin®, PES	Western blot analysis	33
Structure determination of coronavirus SARS-CoV-1 non-structural protein 1 (nsp1)	Concentration of coronavirus nsp1 during purification process	Vivaspin®, PES	Crystalization screening	34
Structure determination of the ADRP domain of Feline Coronavirus (FCoV) non-structural protein 3 (nsp3)	Concentration of coronavirus nsp3 during purification process	Vivaspin®, PES (10 kDa)	Crystalization screening	35
Investigation into the role of three transmembrane proteases in the activation of SARS-CoV-1 spike (S) protein	Concentration of VLPs from HEK 293T cell culture supernatant	Vivaspin®, PES	Cell-cell fusion assay	36
Cryo-electron microscopy of Human Coronavirus HCoV-NL63 spike glycoprotein trimer that is a potential target for neutralizing antibodies during infection	Concentration of recombinant HCoV-NL63 viruses from clarified Drosophila S2 cell culture supernatant	Vivaflow®, PES (10 kDa)	Affinity purification	37

Table 5: Summarized examples of ultrafiltration applications with Sartorius Vivaflow® and Vivaspin® of coronavirus viron and VLP samples

Goal of Research (Type of virus, Host Organism)	Purpose of Filtration (Buffer System)	Sartorius Ultrafiltration Device (MWCO)	Subsequent Step	Ref.
Characterisation of phenotypic changes in virus isolates, such as MERS-CoV, that could relate to pandemic potential	Concentration of MERS-CoV virus isolates	Vivaspin®, PES (100 kDa)	Quantification using plaque titration Viral RNA sequencing analysis	38
Investigation of antiviral potential of Echinacea purpurea (Echinaforce®) against human coronaviruses; SARS-CoV and MERS-CoV	Concentration of MERS-CoV and SARS-CoV virus dilutions	Vivaspin® 20, PES	Limiting dilution assay (TCID ₅₀)	39
Investigation into inactivation of SARS-CoV-2 through heating and chemical protocols	Concentration and separation of deactivated SARS-CoV-2 from lysis buffer	Vivaspin® 500, PES	Inoculation onto Vero-E6 monolayer	40
Investigation of viral and cellular determinants governing hCoV-EMC entry into host cells	Concentration of SARS-CoV and hCoV-EMC virus like particles (VLPs)	Vivaspin®, PES	Western blot analysis	41, 42

Concentration and capture of virions and | or viral RNA in wastewater

In humans and birds Coronaviruses may inflict mild to fatal respiratory tract infections, but in other animal groups a range of other diseases may also occur, such as hepatitis and neurological illness⁴². SARS-CoV-2 is the most recent among a string of Coronavirus epidemics, which early indications suggest that due to its high infectivity, rates of asymptomatic infection, significant incubation time, our relatively limited knowledge of transmission dynamics and overall lack of global pandemic preparation, has evolved into a true global pandemic and has caused significant impact on global health, society and economy.

The severity of this pandemic is driving increased research and funding in all associated areas. One area is on the tracking and epidemiological studies of SARS-CoV-2 infections. One area of focus is in the use of regional wastewater systems, where the compartmentalisation of these systems offers distinct tracking in real time, without the lag for symptom appearance and clinical diagnosis⁴³. In addition, the data collected can be used as a supplemental and low-cost surveillance indicator on the circulation of the virus in a community without the need to screen individuals. Further, it contributes to the tracking of infection prevalence, by adding another epidemic indicator⁴⁴.

RT-PCR is the standard method to test for SARS-CoV-2, but samples typically require concentration and removal of non-Coronavirus material prior to testing to ensure optimal results. Ultrafiltration is a successful method for this⁴³, and some examples have been given in Table 6.

Table 5: Summarized examples of ultrafiltration applications with Sartorius Vivaflow® and Vivaspin® of virus and viral RNA in wastewater samples

Goal of Research (Type of virus, Host Organism)	Purpose of Filtration (Buffer System)	Sartorius Ultrafiltration Device (MWCO)	Subsequent Step	Ref.
Measurement of SARS-CoV-2 RNA in sewage	Concentration of viral RNA	Vivaspin®, PES (50 kDa)	Viral RNA extraction and purification RT-qPCR quantification	43, 44, 45
Benchmarking virus concentration methods for quantification of SARS-CoV-2 in raw wastewater	Concentration of viral RNA	Vivacell, PES (10 kDa)	Viral RNA extraction and purification RT-qPCR quantification	46
Evaluation of two methods to concentrate SARS-CoV-2 from untreated wastewater	Concentration of viral RNA from 40 ml (total) to 700-1000 µl	Vivaspin® (10 kDa)	Viral RNA extraction and RT-qPCR and ddPCR quantification	47
Virus detection in full scale membrane bioreactor (MBR) plant by virus concentration monitoring, inc. Norovirus, Sapovirus and Rotavirus	Concentration of viral particles in effluent	Vivaflow® 50, PES	PEG precipitation Viral RNA quantification	48
Evaluation of membrane bioreactor wastewater virus removal, inc. Norovirus, Sapovirus, Adenovirus	Concentration of effluent from 1 L to 40 mL	Vivaflow® 50, PES	Nucleic acid extraction RT-PCR quantification	49
Evaluation of membrane bioreactor wastewater Norovirus removal	Concentration of viruses in effluent	Vivaflow® 50, PES	Nucleic acid extraction RT-PCR quantification	50
Evaluation of the association between number of hepatitis E cases in the community and concentration in local sewage	Concentration of viruses in effluent	Vivaflow® 50, PES	Nucleic acid extraction RT-qPCR quantification	51

Conclusion

The purification of virus by ultrafiltration is virtually independent of the chemical properties and the structure of the virus particles. As viruses have a size within the range of about 20 nm up to several hundred nanometers, they are typically several orders of magnitude bigger than even the biggest protein complexes.⁷ Therefore, most viruses are unfailingly retained on membranes with large MWCOs of up to 1,000 kDa. The exact specifications of the ideal ultrafiltration membranes depend on the purpose of the subsequent application.

Ultrafiltration for the concentration of Coronavirus species plays an important role in a range of workflows. Perhaps due to the size distribution of viruses and VLPs, the exact MWCO used is not standard across each study. Although typically, for 80–120 nm particles the 100 kDa MWCO would provide the optimal balance between recovery, removal of interfering substances, speed and shear stresses. Whereas for the recovery of RNA material, lower MWCOs (10–50 kDa) are recommended to capture a greater range of RNA chain lengths. However, until further standardisation is confirmed for each application, it is prudent to test specific devices before implementing into procedures.

During the preparation of viral vectors for medical studies, a buffer exchange after column purification can be performed with various MWCOs of all sizes.^{8,9,10,15} To separate virus particles from small proteins, a 1,000 kDa cut off has been shown to work.¹² For the complete removal of HCV from blood serum a 30 kDa MWCO has been utilized.¹¹ When the assessment of whole virus content is crucial (e.g. food, drinking water or marine water samples) smaller MWCOs (5–100 kDa) are used to ensure full recovery of virus particles.^{19,21,22,25–29}

Abbreviations

AAV	Adeno-associated virus
CNS	Central nervous system
	DNA Deoxyribonucleic acid
CoV	Coronavirus
ELISA	Enzyme-linked immunosorbent assay
FPLC	Fast protein liquid chromatography
fCoV	Feline Coronavirus
hCoV	Human Coronavirus
HCV	Hepatitis C virus
HIV	Human immunodeficiency virus
kDa	Kilodalton (1000 g per mole)
M	Molarity (mole per litre)
MERS	Middle east respiratory syndrome
mol	Mole
MWCO	Molecular weight cut-off
nsp	Nonstructural protein
PBS	Phosphate buffered saline
PCR	Polymerase chain reaction
PEG	Polyethylene Glycol
PES	Polyethersulfone
RNA	Ribonucleic acid
SARS	Severe acute respiratory syndrome
RBD	Receptor binding domain
BLI	Bio-Layer Interferometry
RT-PCR	Reverse transcriptase-polymerase chain reaction
ddPCR	Droplet digital polymerase chain reaction
TCIDP50	50% Tissue culture infective dose
VLP	Virus-like particle

References

- Vannucci, L., Lai, M., Chiuppesi, F., Ceccherini-nelli, L. & Pistello, M. Viral vectors : a look back and ahead on gene transfer technology. *New Microb.* 36, 1-22 (2013).
- Luo, D. & Saltzman, W. M. Synthetic DNA delivery systems. *Nat. Biotechnol.* 8, 33-37 (2000).
- Ura, T., Okuda, K. & Shimada, M. Developments in Viral Vector-Based Vaccines. *Vaccines* 2, 624-41 (2014).
- Mingozi, F. & High, K. A. Therapeutic in vivo gene transfer for genetic disease using AAV: progress and challenges. *Nat Rev Genet* 12, 341-355 (2011).
- Soule, H., Genoulaz, O., Gratacap-Cavallier, B. Chevallier, P., Liu, J.-X. & Seigneurin, J.-M. Ultrafiltration and reverse transcription-polymerase chain reaction: an efficient process for poliovirus, rotavirus and hepatitis A virus detection in water. *Water Res.* 34, 1063-1067 (2000).
- Bergh, O., BOrsheim, K. Y., Bratbak, G. & Heldal, M. High abundance of viruses found in aquatic environments. *Nature* 340, 467-468 (1989).
- Hulo, C. et al. ViralZone: A knowledge resource to understand virus diversity. *Nucleic Acids Res.* 39, 576- 582 (2011).
- Hagen, S. et al. Modular adeno-associated virus (rAAV) vectors used for cellular virus-directed enzyme prodrug therapy. *Sci. Rep.* 4, 3759 (2014).
- Schindler, S. E. et al. Photo-activatable Cre recombinase regulates gene expression in vivo. *Sci. Rep.* 5, 13627 (2015).
- Scherr, M. et al. Efficient gene transfer into the CNS by lentiviral vectors purified by anion exchange chromatography. *Gene Ther.* 9, 1708-1714 (2002).
- Granato, M. et al. HCV derived from sera of HCV- infected patients induces pro-fibrotic effects in human primary fibroblasts by activating GLI2. *Sci. Rep.* 6, 30649 (2016).
- Martin, L. et al. Rational design of a CD4 mimic that inhibits HIV-1 entry and exposes cryptic neutralization epitopes. *Nat. Biotechnol.* 21, 71-76 (2003).
- Karlas, A. et al. A human genome-wide loss-of-function screen identifies effective chikungunya antiviral drugs. *Nat. Commun.* 7, 11320 (2016).
- Nestola, P. et al. Rational development of two flowthrough purification strategies for adenovirus type 5 and retro virus-like particles. *J. Chromatogr. A* 1426, 91-101 (2015).
- Carvalho, L. S. et al. Long-term and age-dependent restoration of visual function in a mouse model of CNGB3-associated achromatopsia following gene therapy. *Hum. Mol. Genet.* 20, 3161-3175 (2011).
- Guidelines for drinking-water quality - 4th ed. World Health Organization 2011.
- Wyn-Jones, a P. & Sellwood, J. Enteric viruses in the aquatic environment. *J. Appl. Microbiol.* 91, 945-962 (2001).
- Botzenhart, K. Viren im Trinkwasser. *Bundesgesundheitsblatt - Gesundheitsforsch. - Gesundheitsschutz* 50, 296-301 (2007).
- Lee, K. B., Lee, H., Ha, S. D., Cheon, D. S. & Choi, C. Comparative analysis of viral concentration methods for detecting the HAV genome using real-time RT-PCR amplification. *Food Env. Virol.* 4, 68-72 (2012).
- Bosch, A. et al. Analytical Methods for Virus Detection in Water and Food. *Food Anal. Methods* 4, 4-12 (2011).
- Masciopinto, C. et al. Unsafe tap water in households supplied from groundwater in the Salento Region of Southern Italy. *J. Water Health* 5, 129-148 (2007).
- Dreier, J., Störmer, M., Mäde, D., Burkhardt, S. & Kleesiek, K. Enhanced reverse transcription-PCR assay for detection of norovirus genogroup I. *J. Clin. Microbiol.* 44, 2714-2720 (2006).
- Maunula, L. et al. Enteric Viruses in a Large Waterborne Outbreak of Acute Gastroenteritis in Finland. *Food Environ. Virol.* 1, 31-36 (2009).
- Wommack, K. E. & Colwell, R. R. Virioplankton: viruses in aquatic ecosystems. *Microbiol. Mol. Biol. Rev.* 64, 69-114 (2000).
- Parada, V., Baudoux, A.-C., Sintes, E., Weinbauer, M. G. & Herndl, G. J. Dynamics and diversity of newly produced virioplankton in the North Sea. *ISME J.* 2, 924-936 (2008).
- Schroeder, D. C., Oke, J., Hall, M., Malin, G. & Wilson, W. H. Virus Succession Observed during an Emiliana huxleyi Bloom. *Virus. Appl. Environ. Microbiol.* 69, 2484- 2490 (2003).
- Allen, M. J. et al. Locus-Specific Gene Expression Pattern Suggests a Unique Propagation Strategy for a Giant Algal Virus. *J. Virol.* 80, 7699-7705 (2006).
- Brussaard, C. P. D., Noordeloos, A. A. M., Sandaa, R. A., Heldal, M. & Bratbak, G. Discovery of a dsRNA virus infecting the marine photosynthetic protist *Micromonas pusilla*. *Virology* 319, 280-291 (2004).
- Fischer, M. G. & Hackl, T. Host genome integration and giant virus-induced reactivation of the virophage mavirus. *Nature* 540, 288-291 (2016).
- Liu, H. et al. Cross-Neutralization of a SARS-CoV-2 Antibody to a Functionally Conserved Site Is Mediated by Avidity. *Immunity* 53, 1272-1280 (2020)

31. Tan, T.K. et al. A COVID-19 vaccine candidate using SpyCatcher multimerization of the SARS-CoV-2 spike protein receptor-binding domain induces potent neutralizing antibody responses. *Nature Communications* 12:542 (2021)
32. Kuate, S. et al. Exosomal vaccines containing the S protein of the SARS coronavirus induce high levels of neutralizing antibodies. *Virology* 362, 26-37 (2007).
33. Würdinger, T. et al. Soluble receptor-mediated targeting of mouse hepatitis coronavirus to the human epidermal growth factor receptor. *J. Virology* 79, 15314-15322 (2005)
34. Jansson, A. M. Structure of Alphacoronavirus transmissible gastroenteritis virus nsp1 has implications for coronavirus nsp1 function and evolution. *J Virology* 87, 2949-2955 (2013),
35. Justyna, A. et al. Structure of the X (ADRP) domain of nsp3 from feline coronavirus. *Bio. Crystallography* 65, 1292-1300 (2009).
36. Bertram, S, et al. Cleavage and activation of the severe acute respiratory syndrome coronavirus spike protein by human airway trypsin-like protease. *J Virology* 85, 13363-13372 (2011).
37. Walls, A. C. et al. Glycan shield and epitope masking of a coronavirus spike protein observed by cryo-electron microscopy. *Nature Str & M. Biology* 23 (2016).
38. Schroeder, S. et al. Functional comparison of MERES-coronavirus lineages reveals increased zoonotic potential of the recombinant lineage 5. Pending publication.
39. Signer, J. In vitro antiviral activity of Echinaforce®, an Echinacea purpurea preparation, against common cold coronavirus 229E and highly pathogenic MERS-CoV and SARS-CoV. Pending publication
40. Pastorino, B. et al. Evaluation of heating and chemical protocols for inactivating SARS-CoV-2. Pending publication
41. Gierer, S. et al. The spike protein of the emerging betacoronavirus EMC uses a novel coronavirus receptor for entry, can be activated by TMPRSS2, and is targeted by neutralising antibodies. *J. Virology* 87, 5502-5511 (2013).
42. Masters, P. S. The Molecular Biology of Coronaviruses. *Adv. Virus Research* 66, 193-292 (2006)
43. Larsen, D. A, Wigginton K R. Tracking COVID-19 with Wastewater. *Nature Biotechnology* 38, 1151-1153 (2020).
44. Trottier, J. et al. Post-lockdown detection of SARS-CoV-2 RNA in the in the wastewater of Montpellier, France. *One Health* 10 (2020)
45. Hokkaido University. SARS-CoV-2 RNA Detected in Untreated Wastewater from Louisiana. *ScienceDaily*. (26 Aug 2020)
46. Jafferalli, M. H. et al. Benchmarking virus concentration methods for quantification of SARS-CoV-2 in raw wastewater. *Science of the Total Environment* 10, 755. (2021)
47. Dumke, R. et al. Evaluation of Two Methods to Concentrate SARS-CoV-2 from Untreated Wastewater. *Pathogens* 195 (2021)
48. Takayuki, M. et al. Virus type-specific removal in a full-scale membrane bioreactor treatment process. *Food and Env. Virology* 10, 176-186 (2017).
49. Sima, L. C. et al. Calicivirus removal in a membrane bioreactor wastewater treatment plant. *Applied and Inv. Microbiology* 77, 5170-5177 (2011).
50. Schaeffer, J. et al. Improving the efficacy of sewage treatment decreases norovirus contamination in oysters. *Int. J. of Food Microbiology* 286, 1-5 (2018).
51. Takayuki, M. et al. Detection of hepatitis E virus in sewage after an outbreak on a French island. *Food and Env. Virology* 8, 194-100 (2016).

5. Sartorius Ultrafiltration Products in the Preparation of Biological Nanoparticles and Medical Nanocarriers



Introduction

Paul Ehrlich was inspired by the idea of the “magic bullet”^{*} when he for the first time described in theory toxic drugs assembled to so-called “Nanocarriers” in 1908.¹ Today, Nanocarriers have found multiple applications in modern medicine and biotechnology. A key application for these special nanomaterials is a targeted delivery of drugs where they act as transport modules (i. e. as nanoparticles, vesicles, or micelles) for the active ingredient.^{2,3,4,5} This is assumed to be more effective and less toxic to the (human) organism compared to traditionally administered drug substances.⁶ Besides drug delivery, various further fields using Nanocarriers evolved during the last decades; e. g. magnetic resonance imaging or stem cell gene therapy with metal-based nanoparticles,^{7,8} or optical imaging with quantum dots.⁹

Nanocarriers can be categorized by their starting material (i. e. metal-, lipid-, polymer-, and protein-based) and by their formation after preparation (i. e. vesicles, particles and micelles). In general, the preparation of a nanoparticle suspension or a vesicle dispersion in an aqueous medium consists of three steps: a) assembly of the Nanocarriers (for example, by injections, film hydration, or reverse phase evaporation), b) purification (for example, by chromatography, dialysis or ultrafiltration), and c) concentration (for example, by ultrafiltration or evaporation).

This short review provides examples of recent literature dealing with the preparation of Nanocarriers. Particular focus is laid on the concentration and purification steps which were performed via ultrafiltration with Sartorius VivaSpin[®] or Vivaflow[®] devices with different pore sizes (respectively

molecular weight cut-off, MWCO). The VivaSpin[®] portfolio spans a volume range from 0.1 to 20 mL, whereas the Vivaflow[®] system covers volumes from 0.1 to 5 liters. Thus, Sartorius offers an unrivaled wide range of processable sample volumes, membrane materials and MWCOs to meet the different requirements of their intended use. Challenges in this context are buffer exchange after synthesis, desalting and washing,^{10,11} exclusion of solubilized compounds,^{12,13,14} or aggregates.¹⁵

Purification is essential to obtain isotonic conditions for in vivo applications, to prevent aggregation or agglomeration and to remove free toxic drugs, ligands, or other substrates potentially triggering side effects. Concentration steps are essential to adjust the amount of pharmaceutical active ingredient in the drug and achieve the anticipated therapeutic or diagnostic effect.

During purification, the separation of free substances (starting material) from the desired Nanocarriers via size-exclusion chromatography (SEC) leads to an unavoidable dilution and to the necessity of a subsequent concentration step. In contrast, dialysis purifies without significant dilution but a concentration step can still be mandatory, if higher Nanocarrier concentrations are necessary. Both separation methods require quite extensive, costly and timeconsuming manual handling. This drawback is overcome with the ultrafiltration utilized by centrifugation in VivaSpin[®] or with a peristaltic pump for the Vivaflow[®] system. This technique is less expensive and quickly performed with very little manual input. Noteworthy is that purification and concentration steps are performed simultaneously.¹⁶

^{*} In German “Zauberkegel”, opera “Freischütz” by Carl Maria von Weber

After the Nanocarrier is purified, the determination of drug loading (conjugation or encapsulation efficiency) is commonly performed. The conjugation or encapsulation efficiency is one of the reference values to describe and characterize Nanocarriers. Other important properties are the zeta potential and the size distribution determined via photon correlation spectroscopy (PCS), high-resolution transmission electron microscopy (HRTEM) imaging, or dynamic light scattering (DLS).

Prior to performing these different characterizations, a successful purification and concentration of the suspension or dispersion is essential.

In the following tables you can find an overview of publications using ultrafiltration steps for the purification and concentration of different kinds of Nanocarriers. Table 2 provides guidance on which devices and MWCOs to use.

Table 1 summarizes examples of Nanocarrier ultrafiltration applications with Sartorius Vivaspin® or Vivaflow®:

Nanocarrier: Nanoparticle, Vesicle, Micelle	Size distribution obtained via (HR)TEM or DLS, Z-Average via PCS and others-if reported	Application	Ref.
Nanoparticles from metal, metal oxides and functionalized metals			
Iron oxides nanoparticles with cisplatinbearing polymer coating	SD: 4.5 ± 0.9 nm via X-Ray-Diffraction Analysis	Magnetic resonance imaging	7
Functionalized iron oxide nanoparticles	SD: 38 and 40 nm via DLS	Stem cell gene therapy and tracking	8
Gold nanoparticles	SD: 0.8 – 10.4 nm via Atomic Force Microscopy	Antimicrobial activity	17
Protein coated gold nanoparticles	SD: 15 and 80 nm via TEM	Drug delivery	18
Functionalized gold nanoparticles	Core-SD: 2 nm via TEM	Targeted imaging tool and antigen delivery	19
Functionalized gadolinium-based nanoparticles	Z-Average: 1.1 ± 0.6 nm and 4 – 14 nm	Diagnostic and therapeutic application	20, 21
Functionalized nanocrystals	10 to 20 nm	Quantum dots for imaging	9
Nanoparticles from polymers, functionalized polymers and polymersomes			
Polymer based Nanoparticles		Drug delivery	22
Curdlan coated polymer nanoparticles	Z-Average: 280 – 480 nm depending on the composition	Macrophage stimulant activity and drug delivery	23
Docetaxel-carboxymethylcellulose Polymer Nanoparticles	Z-Average: 118 ± 1.8 nm	Anti-cancer efficacy studies	4
Functionalized Polymersomes	Z-Average: 185 nm	Surface functionalization studies	3

Lipid Nanoparticles and Liposomes			
Liposomes and micelles	Z-Average: 100 nm for Liposomes and 15 nm for micelles	Ischemia-reperfusion injury	25
Solid lipid Nanoparticles	Z-Average: 100 – 120 nm depending on the used lipid	Drug delivery (Brain Targeting)	26
Bacterial outer membrane vesicles	SD: 124 nm via TRPS	Tunable resistive pulse sensing (TRPS) Analysis	27
Bacterial outer membrane vesicles		Basic research	28
Bacterial outer membrane vesicles	SD: 95 nm	Basic research	29
Bacterial outer membrane vesicles	SD: 50 – 150 nm via TEM	Basic research	30
Liposomes		Drug delivery	2
Liposomes		Encapsulated hydrophilic drugs (Drug delivery)	31
Micelles			
Micelles		Drug delivery	4
Hydrophobic drug micelles based on polymers	SD via DLS: 39 – 165 nm depending on compound in use	Drug delivery	14
Protein Nanoparticles			
Protein Nanoparticles	SD: 20 – 40 nm via DLS	Drug carrier studies	32

SD = Size distribution

Table 2 lists example Sartorius devices and typical MWCOs used for each nanocarrier ultrafiltration application:


Nanocarrier: Nanoparticle, Vesicle, Micelle	Sartorius Ultrafiltration Device	MWCO	Ultrafiltration purpose	Ref.
Nanoparticles from metal, metal oxides and functionalized metals				
Iron oxides nanoparticles with cisplatinbearing polymer coating	Vivaspin® 20	100 kDa	Purification and concentration	7
Functionalized iron oxide nanoparticles	Vivaspin® 20	100 kDa	Washing step	8
Gold nanoparticles	Vivaspin® 20	5 kDa	Purification step	17
Protein coated gold nanoparticles	Vivaspin® 6	10 kDa	Separation of Nanoparticles Dyes and washing	18
Functionalized gold nanoparticles	Vivaspin®	10 kDa	Purification step	19
Functionalized gadolinium-based nanoparticles	Vivaspin®	5 kDa and 10 kDa	Purification and concentration	20, 21
Functionalized nanocrystals	Vivaspin®	300 kDa and 50 kDa	Separation of quantum dots-antibody conjugates from starting material (prior to enumeration)	9
Nanoparticles from polymers, functionalized polymers and polymersomes				
Polymer based Nanoparticles	Vivaspin®	30 kDa	Purification and concentration	22
Curdlan coated polymer nanoparticles	Vivaspin® 20	3 kDa	Washing	23
Docetaxel-carboxymethylcellulose Polymer Nanoparticles	Vivaspin®	10 kDa	Concentration	4
Functionalized Polymersomes	Vivaspin® 20	10 kDa	Concentration	3
Lipid Nanoparticles and Liposomes				
Liposomes and micelles	Vivaspin® 20	100 kDa	Concentration	25
Solid lipid Nanoparticles	Vivaflow® 50	100 kDa	Purification	26
Bacterial outer membrane vesicles	Vivaflow® 200	100 kDa	Buffer exchange and concentration	27
Bacterial outer membrane vesicles	Vivaspin® 500 and 20	100 kDa	Buffer exchange and concentration	28
Bacterial outer membrane vesicles	Vivaflow® 200	100 kDa	Buffer exchange and concentration	29
Bacterial outer membrane vesicles	Vivaspin®	100 kDa	Buffer exchange and concentration	30
Liposomes	Vivaspin®	100 kDa	External buffer exchange	2
Liposomes	Vivaflow® 50	100 kDa	Elimination of the free drug	31
Micelles				
Micelles	Vivaspin®	30 kDa	Separation of free substrate and concentration	4
Hydrophobic drug micelles based on polymers	Vivaflow®		Surfactant removal	14
Protein Nanoparticles				
Protein Nanoparticles	Vivaspin® 500	3 kDa	Separation of the free from the encapsulated drug (Drug binding quantification by subsequent UV-vis analysis)	32

References

- 1 Strebhardt, K. & Ullrich, A.: Paul Ehrlich's magic bullet concept: 100 years of progress. *8*, 473–480 (2008).
- 2 Jakoby, J., Beuschlein, F., Mentz, S., Hantel, C. & Süß, R.: Liposomal doxorubicin for active targeting: Surface modification of the nanocarrier evaluated in vitro and in vivo – challenges and prospects. *Oncotarget* *6*, (2015).
- 3 Klermund, L., Poschenrieder, S. T. & Castiglione, K.: Simple surface functionalization of polymersomes using non-antibacterial peptide anchors. *J. Nanobiotechnology* *14*, 48 (2016).
- 4 Mulder, W. J. M. et al.: Molecular imaging of macrophages in atherosclerotic plaques using bimodal PEG-micelles. *Magn. Reson. Med.* *58*, 1164–1170 (2007).
- 5 Murthy, S. K.: Nanoparticles in modern medicine: state of the art and future challenges. *Int. J. Nanomedicine* *2*, 129–41 (2007).
- 6 Voigt, R. & Fahr, A.: *Pharmazeutische Technologie für Studium und Beruf*. Deutscher Apotheker Verlag, 10th Edition (2010).
- 7 Unterweger, H. et al.: Development and characterization of magnetic iron oxide nanoparticles with a cisplatin-bearing polymer coating for targeted drug delivery. *Int J Nanomedicine* *9*, 3659–3676 (2014).
- 8 Park, W. et al.: Multi-modal transfection agent based on monodisperse magnetic nanoparticles for stem cell gene delivery and tracking. *Biomaterials* *35*, 7239–7247 (2014).
- 9 Chalmers, N. I. et al.: Use of quantum dot luminescent probes to achieve single-cell resolution of human oral bacteria in biofilms. *Appl. Environ. Microbiol.* *73*, 630–636 (2007).
- 10 Hoffman, L. W., Andersson, G. G., Sharma, A., Clarke, S. R. & Voelcker, N. H.: New insights into the structure of PAMAM dendrimer | gold nanoparticle nanocomposites. *Langmuir* *27*, 6759–6767 (2011).
- 11 Rademacher, T. & Williams, P.: Nanoparticle-peptide compositions. (2014).
- 12 Allard, E. & Larpent, C.: Core-shell type dually fluorescent polymer nanoparticles for ratiometric pH-sensing. *J. Polym. Sci. Part A Polym. Chem.* *46*, 6206–6213 (2008).
- 13 Prach, M., Stone, V. & Proudfoot, L.: Zinc oxide nanoparticles and monocytes: Impact of size, charge and solubility on activation status. *Toxicol. Appl. Pharmacol.* *266*, 19–26 (2013).
- 14 Zhang, Y. et al.: Therapeutic surfactant-stripped frozen micelles. *Nat Commun* *7*, 11649 (2016).
- 15 Klasson, A. et al.: Positive MRI contrast enhancement in THP-1 cells with Gd₂O₃ nanoparticles. *Contrast Media Mol. Imaging* *3*, 106–111 (2008).
- 16 Simonoska Crcarevska, M. et al.: Definition of formulation design space, in vitro bioactivity and in vivo biodistribution for hydrophilic drug loaded PLGA/PEO-PPO-PEO nanoparticles using OFAT experiments. *Eur. J. Pharm. Sci.* *49*, 65–80 (2013).
- 17 Boda, S. K. et al.: Cytotoxicity of Ultrasmall Gold Nanoparticles on Planktonic and Biofilm Encapsulated Gram-Positive Staphylococci. *Small* *11*, 3183–3193 (2015).
- 18 Schäffler, M. et al.: Blood protein coating of gold nanoparticles as potential tool for organ targeting. *Biomaterials* *35*, 3455–3466 (2014).
- 19 Arosio, D. et al.: Effective targeting of DC-sign by α-fucosylamide functionalized gold nanoparticles. *Bioconjug. Chem.* *25*, 2244–2251 (2014).
- 20 Miladi, I. et al.: Biodistribution of ultra small gadolinium-based nanoparticles as theranostic agent: application to brain tumors. *J. Biomater. Appl.* *28*, 385–94 (2013).
- 21 Faure, A. C. et al.: Control of the in vivo biodistribution of hybrid nanoparticles with different poly(ethylene glycol) coatings. *Small* *5*, 2565–2575 (2009).
- 22 Benita, S., Debotton, N. & Goldstein, D.: Nanoparticles for Targeted Delivery of Active Agent. (2008).
- 23 Tukulula, M. et al.: Curdlan-conjugated PLGA nanoparticles possess macrophage stimulant activity and drug delivery capabilities. *Pharm. Res.* *32*, 2713–2726 (2015).
- 24 Ernsting, M. J., Tang, W. L., MacCallum, N. W. & Li, S. D.: Preclinical pharmacokinetic, biodistribution, and anti-cancer efficacy studies of a docetaxel-carboxymethyl-cellulose nanoparticle in mouse models. *Biomaterials* *33*, 1445–1454 (2012).
- 25 Geelen, T., Paulis, L. E., Coolen, B. F., Nicolay, K. & Strijkers, G. J.: Passive targeting of lipid-based nanoparticles to mouse cardiac ischemia-reperfusion injury. *Contrast Media Mol. Imaging* *8*, 117–126 (2013).

- 26 Neves, A. R., Queiroz, J. F. & Reis, S.: Brain-targeted delivery of resveratrol using solid lipid nanoparticles functionalized with apolipoprotein E. *J. Nanobiotechnology* 14, 27 (2016).
- 27 Bogomolny, E. et al.: Analysis of bacteria-derived outer membrane vesicles using tunable resistive pulse sensing. *Prog. Biomed. Opt. Imaging – Proc. SPIE* 9338, 4–9 (2015).
- 28 Blenkiron, C. et al.: Uropathogenic *Escherichia coli* releases extracellular vesicles that are associated with RNA. *PLoS One* 11, 1–16 (2016).
- 29 Twu, O. et al.: *Trichomonas vaginalis* Exosomes Deliver Cargo to Host Cells and Mediate Host:Parasite Interactions. *PLoS Pathog.* 9, 22–24 (2013).
- 30 Tong, T. T., Mörgelin, M., Forsgren, A. & Riesbeck, K.: *Haemophilus influenzae* Survival during Complement-Mediated Attacks Is Promoted by *Moraxella catarrhalis* Outer Membrane Vesicles. *J. Infect. Dis.* 195, 1661–1670 (2007).
- 31 Prado, J. M. D., Antoranz, J. R. C., Barroeta, M. Á. E., Barroeta, B. E. & Diaz, M. C.: Liposomal formulations. (2009).
- 32 Achilli, E. et al.: Preparation of protein nanoparticle by dynamic aggregation and ionizing-induced crosslinking. *Colloids Surfaces A Physicochem. Eng. Asp.* 486, 161–171 (2015).

6. Vivaflow® and Vivaspin® Workflow in Protein Research Laboratories



Introduction

Efficiency and efficacy of a multiple cycle experimental procedure was performed using Vivaflow® tangential flow cassettes for initial concentration and diafiltration of a cell culture supernatant. This was followed by Vivapure® Ion Exchange spin columns for the protein purification step and finally Vivaspin® 20 ultrafiltration devices for the final sample concentration and desalting. An artificial mixture of proteins in a RPMI-1640 culture medium was created to mimic the type of product that many researchers culture using e.g. the UniVessel device. This procedure further reflects a method that can be adapted to a large number of protein purification protocols, selecting alternative MWCOs and device sizes where necessary.

Methods

Part 1 – Creating and concentrating the culture medium

2 bottles (4 g) of RPMI-1640 were dissolved into 1.8 L dd-H₂O and 4 g of sodium acetate was added.

The pH was adjusted to 7.2 using 1M HCl. 2 g of BSA and 1 g of Lysozyme was added as protein samples, meant to be separated by chromatography. The volume of the cell culture supernatant sample was brought up to 2 L using dd-H₂O. After every preparation, concentration and purification step, 1 mL sample was set aside for SDS gel analysis.

Ion exchange chromatography was the method of choice for purifying lysozyme from the cell culture supernatant, especially from the “contaminant” BSA. For this, the 2 L cell culture supernatant needed to be concentrated and then diafiltered to adjust the sample to the starting conditions needed for the ion exchange chromatography binding step.

For concentration and diafiltration, the Vivaflow® 200 was used with a 5 kDa MWCO PES membrane. Vivaflow® 200 is a ready-to-use laboratory cross-flow cassette in an acrylic housing, which allows caustic cleaning and 4–5 re-uses. Two cassettes can be run in parallel for the concentration of up to 5 L sample volumes. For the 2 L sample to be concentrated in this experiment, one cassette was sufficient. A Masterflex pump with an Easy-Load, size 16 pump head was used to run the Vivaflow® 200 cassette. Figure 1a. and 1b. show the Vivaflow® 200 set up with one or two cassettes.



The Vivaflow® 200 system was operated at 3 bar. Once 1.8 L of filtrate had been collected, the pump was stopped, the tubes removed from the cell culture medium concentrate and filtrate and the Vivaflow® system was purged with dd-H₂O. This solution now contained a 10-fold concentration of the constituent proteins from the original culture medium.

A BCA protein detection assay conveyed a 100% recovery of protein after this first concentration step. Table 1 indicates the time needed for the sample concentration.

Filtrate Volume (mL)	Time (hr:min:sec)
0	0:00:00
100	0:03:16
200	0:06:50
300	0:10:45
400	0:14:38
500	0:18:36
600	0:22:43
700	0:26:57
800	0:31:14
900	0:36:01
1,000	0:40:50
1,100	0:45:46
1,200	0:50:36
1,300	0:55:32
1,400	1:00:24
1,500	1:05:26
1,600	1:10:28
1,700	1:15:52
1,800	1:21:50




Figure 1a. Vivaflow® 200 – single cassette

Table 1: Vivaflow® 200, 5 kDa MWCO PES concentration speed




Figure 1b. Vivaflow® 50R – two cassettes

Part 2 – Buffer exchange of culture medium concentrate

The Vivaflow® 200 system was used for fast and easy diafiltration. To this end, the diafiltration reservoir, a Vivaflow® accessory, was filled with the 200 mL concentrated sample. Figure 2 shows the diafiltration set up. The Vivaflow® 200 system was set up as before, however attaching an additional tube to the diafiltration lid and placing this new inlet tube into a 25 mM sodium acetate (pH 5.5) buffer (needed to re-adjust the sample concentrate to the ionic starting conditions of the ion exchange chromatography step which was to follow). This enables concentration of the sample in the reservoir to the extent that the original buffer is removed and collected as waste (filtrate). Simultaneously, new buffer (25 mM sodium acetate) is drawn into the closed system, gradually leading to a buffer exchange while keeping the sample volume constant at 200 mL. The system was run at 3 bar. Once 1 L of buffer had been exchanged, diafiltration was stopped.

The 200 mL solution now contained the correct buffer to maintain the stability of the proteins of interest for the next part of the protocol and had the correct pH and salt concentration for the ion exchange binding step. BCA protein quantification again showed a 100% protein recovery.

Table 2 shows the time needed for diafiltration of 200 mL sample against 1,000 mL exchange buffer, again using Vivaflow® 200 with a 5 kDa MWCO PES membrane.



Figure 2: Diafiltration system set up for buffer exchange. Culture medium concentrate can be seen in the center of the image. 1 L 25 mM sodium acetate (exchange buffer) can be seen connected to the system on the left of the image.

Filtrate Volume (mL)	Time (hr:min:sec)
0	0:00:00
100	0:06:58
200	0:14:16
300	0:22:39
400	0:29:40
500	0:37:02
600	0:44:15
700	0:51:34
800	0:58:54
900	1:06:03
1,000	1:13:02

Table 2: Diafiltration of 200 mL concentrated cell culture supernatant containing the proteins lysozyme and BSA against 1,000 mL 25 mM sodium acetate.

Part 3 – Purification of Lysozyme, the protein of interest

The purification of Lysozyme was performed using a Vivapure® cation exchange membrane adsorber device (Vivapure® Maxi H S). The membrane adsorber matrix holds the active ligands and performs like a traditional cation exchanger. However, membrane adsorbers represent a special form of chromatography matrix. Unlike traditional resins, they make use of convective transport to bring proteins to the ion exchange surface; hence, binding, washing and elution is performed quickly and high binding capacities are achieved, even at high flow rates. This allows their use in fast and convenient centrifugal spin columns (Figure 3).




Figure 3: Electron micrograph of chromatography gel beads (upper right) in comparison to a Q ion exchange membrane adsorber (background), revealing 100-fold larger pore sizes of the membrane adsorber.

Two Vivapure® Maxi H S type devices (Figure 4) were each equilibrated with 10 mL of 25 mM sodium acetate (pH 5.5), by filling with 10 mL of this buffer and centrifuging for 5 min. in a swing bucket centrifuge at 500 g and discarding the flow through. Using the concentrated and buffer exchanged sample from Part 2, 10 mL samples were pipetted into each of the equilibrated Vivapure® devices and centrifuged again for 5 min. in a swing bucket centrifuge at 500 g. The Vivapure® devices were washed with 10 mL of 25 mM sodium acetate, discarding the flow through, followed by an elution step with 5 mL of 1 M NaCl in 25 mM sodium acetate. A BCA assay revealed 95% lysozyme recovery.



Figure 4: Vivapure® Maxi spin columns can be used in a centrifuge for fast and easy protein purification.

The eluate was then concentrated in a Vivaspin® 20 (5 kDa MWCO PES), Figure 5, and centrifuged at 5,000 g for 10 min. or until approximately 2 mL of concentrate had been collected. The device was then re-filled with 18 mL 50 mM potassium phosphate buffer (pH 7.2) to 20 mL for a final buffer exchange and desalting of the purified sample. The sample was again centrifuged until a final sample volume of 2 mL had been attained. A BCA assay revealed 97% lysozyme recovery.



Figure 5: Vivaspin® 20 ultrafiltration device, on the right with a pressure cap which allows pressurization of the device as an alternative to the regular centrifugal operation.

Part 4 – Analyzing the samples

The samples of the individual steps were analyzed by SDS gel, using reducing sample buffer (prepared by adding 50 µL 2-mercaptoethanol to 950 µL Laemmli sample buffer). For all steps, 5 µL of the 1 mL sample taken during the experiment were diluted with 95 µL reducing sample buffer, of which 20 µL were loaded onto a 12% Tris-HCl SDS gel (Figure 6).




Figure 6: Coomassie stained 12% Tris-HCl SDS gel loaded with 20 µL sample preparations. Lane 1: Marker (SDS Broad Range); Lane 2: Original sample; Lane 3: Original sample filtrate (Part 1); Lane 4: Marker; Lane 5: Buffer exchange concentrate (Part 2); Lane 6: Filtrate after binding (Part 3); Lane 7: Marker; Lane 8: Filtrate after elution (Part 3); Lane 9: Filtrate after concentrating and desalting (Part 3); Lane 10: Concentrate after concentrating and desalting.

Conclusion

The overall result shows that a standard and straightforward procedure can be followed to concentrate, purify, isolate and analyze a protein of interest from a cell culture, using Vivaflow[®] 200 tangential flow units for cell culture supernatant concentration and diafiltration, Vivapure[®] for ion exchange chromatography, followed by Vivaspin[®] 20 for final sample concentration and desalting.

In many cases dialysis, which is an overnight procedure would be performed instead of the much quicker alternative - ultrafiltration. Here, we show how time-saving and efficient ultrafiltration is for diafiltration and desalting applications, as well as for protein concentration.

The set up and completion of protein purification takes approx. 3.45 h using this method, starting from a culture supernatant, with high protein recoveries in each step (see Table 3). A total protein purification procedure can therefore be completed within 1 working day, including SDS gel analysis, utilizing this strategy, when adapted to individual needs.

Task	Time	Recovery
Vivaflow [®] 200 set up and concentration	1 hour 25 min.	100%
Vivaflow [®] 200 set up and diafiltration	1 hour 20 min.	100%
Vivapure [®] purification	45 min.	95%
Vivaspin [®] concentration and desalting	30 min.	97%
Total	3 hours 45 min.	92%

Products used in this experiment	Order No.
Vivaflow [®] 200, 5 kDa MWCO PES	VF20P1
500 mL Diafiltration Reservoir	VFA006
Vivapure [®] S H Maxi	VS-IX20SH08
Vivaspin [®] 20, 5 kDa MWCO PES	VS2011


Table 3: Processing times for a complete protein purification workflow.

Germany

Sartorius Lab Instruments GmbH & Co. KG
Otto-Brenner-Straße 20
37079 Göttingen
Phone +49 551 308 0

USA

Sartorius Corporation
565 Johnson Avenue
Bohemia, NY 11716
Phone +1 631 254 4249
Toll-free +1 800 635 2906

 For further information, visit
www.sartorius.com